



Professional in

- CLS (Import Range) Plasticware

Composite Lab Line Pvt Ltd

Sea Star® Vortexer®

Multi-purpose, maximum convenience

- Variable speed rotary dial provides finite speed control during vortex mixing
- Stable elastomeric device feet, provide vibration dampening and prevent movement during operation
- Innovative head can be removed for cleaning purposes



Specifications

- Speed Range: Variable 500-3000 rpm (shaking to vortexing)
- Orbit: 4.2 mm
- Control: 3 position switch-continuous, intermittent, manual, and off

Holding capacity (plus a center cup)

- Assay Plate x 1
- 50 mL conical tubes x 2
- 15 mL conical tubes x 2
- 1.5/2 mL microfuge tubes x 4
- 0.5 mL microfuge tubes x 6
- 0.2 mL microfuge tubes x 24



Cat No	Description	Price
CLS120214	Sea Star® Vortexer 220V AC, 50/60 Hz	28,000.00

Sea Star® Digital Orbital Shaker

Multipurpose compact universal shaker

- No-stop optical motor feedback, enables the unit to instantly adjust to weight load changes and unbalanced loads for truly reproducible results
- Suitable for use in incubators and cold rooms within the stated temperature range
- Soft-start feature reduces spill instances from full vessels
- A high weight loading of 4.5 kg, with an overload protection cut out



Specifications

- Speed range: 20 to 300 rpm 1 rpm increments
- Orbit diameter: 19.2 mm
- Timer range: 0 to 49 hours, 1 second increments or continuous
- Maximum load: 4.5 kg



Bungee cords allows you to create your own holding matrix for vessels of all shapes and sizes

Cat No	Description	Price
CLS4010A	Sea Star® Digital Orbital Shaker, 110–240 VAC, 50/60 Hz	95,000.00

Micro Tubes

Made from 100% virgin polypropylene, certified DNase, RNase and Pyrogen free and capable of withstanding high centrifugal forces. The natural colored tubes offer superior clarity, while the pure black version protects and prolongs the life of more light sensitive products.



Microtubes - 0.2, 0.5 and 1.5 mL

- Attached flat caps for easy one handed use, fully autoclavable
- 0.2 mL is suitable for Thermal Cyclers. The slender shape evenly distributes uniform heat
- 0.5 and 1.5 mL versions have a frosted writing area and a syringe pierceable cap

Cat No	Volume	Dimension	Noted Graduation	Pack	Price
CLS4423	0.2 mL	Max OD 6.9 mm	None	1000/pk	11,828.00
CLS4422	0.5 mL	Max OD 10.1 mm	0.1, 0.3, 0.5	1000/pk	7,941.00
CLS4323	1.5 mL	Max OD 12.9 mm	0.1, 0.5, 1.0	500/pk	3,717.00
CLS4323K	1.5 mL	Max OD 13.4 mm	0.1, 0.5, 1.0, 1.5	500/pk	4,292.00
CLS4325	1.5 mL	Max OD 12.9 mm	0.1, 0.5, 1.0	1000/pk	7,398.00

Centrifuge Tubes

Large Capacity Tubes - 15 and 50 mL

- Lot-to-lot traceability with easy tear packaging
- Writing area and flat polyethylene lids for easy labelling
- Imprinted graduations and a self-standing version (only in clear bulk pack)
- Reusable card used for the racked versions, making them more environmentally friendly




Cat No	Volume	Dimension	Noted Graduation	Pack	Price
CLS4426B	15 mL	OD 17 mm + 0.3 mm Length 120 mm + 2 mm	0.5 mL	500/pk	33,719.00
CLS4426R	15 mL	OD 17 mm + 0.3 mm Length 120 mm + 2 mm	0.5 mL	500/pk	36,172.00
CLS4428	15 mL	OD 17 mm + 0.3 mm Length 120 mm + 2 mm	0.5 mL	500/pk	49,000.00
CLS4427B	20 mL	OD 29 mm + 0.3 mm Length 114 mm + 2 mm	5.0 mL	500/pk	46,600.00
CLS4429	50 mL	OD 29 mm + 0.3 mm Length 114 mm + 2 mm	5.0 mL	500/pk	67,000.00
CLS4427S	50 mL	OD 29 mm + 0.3 mm Length 114 mm + 2 mm	5.0 mL	500/pk	47,850.00

Screw-Top Tubes with O-Ring Cap

Screw-on polypropylene lids with TPE O-Ring to create an excellent seal for secure storage.

- Clear polypropylene tubes provide excellent view of samples
- Skirted, self-standing base makes filling and retrieval easy
- Ideal for sample processing and storage application

Cat No	Volume	Dimension	Noted Graduation	Pack	Price
CLS10060	2.0 mL	Max OD 13.3 mm Closed height 46.7 mm	0.5, 1.0, 1.5	500/pk	19,726.00



Sprout® Mini-Centrifuge

Sprout® Package includes:

- Sprout® Mini-Centrifuge
- 6 place 1.5/2.0 mL tube rotor
- 16 place 0.2 mL tube rotor (singles, or 2 x 8 Strip)
- 6 tube adapters 0.2 mL
- 6 tube adapters 0.5 mL
- Storage case for rotors and adapters
- Tube rack
- Power cord

Specifications

- Maximum speed: 6,000 rpm
- Maximum RCF: 2000 x g
- Power: Regional specific versions
- Dimensions: 153 x 128 x 104 mm



Change out rotors with the push of a button



Cat No	Description	Price
CLS10022	Sprout® Mini-Centrifuge	20,000.00

Gusto® High-Speed Mini-Centrifuge

This compact centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement. Ideal for PCR®, microfilter cell separation, clinical & diagnostic applications, and HPLC protocols

- Compact footprint and brush-less motor for quiet operation and low vibration
- Select and set the digital speed display in xg or rpm
- Sample imbalance and tilt safety cut out
- Safety interlock lid prevents opening while rotor is turning
- Made from UV-resistant plastic, increasing its useful life under harsh light conditions
- Double insulated construction
- Push button lid release
- Manual override allows user to open lid in case of power failure



Specifications

- Speed: 500–12,500 rpm
- RCF: 15–9, 800 Xg
- Power: 110 to 240 VAC, 50/60 Hz with 4 interchangeable plug adapters

Gusto® Package includes:


- Gusto™ High-Speed Mini-Centrifuge
- 12 place 1.5/2.0 mL tube rotor
- 32 place 0.2 mL tube rotor (singles, or 4 x 8 Strip)
- 12 tube adapters 0.2 mL
- 12 tube adapters 0.5 mL
- Universal power adapter 110–240V 50/60 Hz, with 4 interchangeable plug type

Cat No	Description	Price
CLS10050	Gusto™, 110–240 VAC, 50/60 Hz	75,000.00

Universal Tube Rack

Quickly change tube sizes with this adaptable polypropylene rack. Open the well with a turn of the knob or directly insert a tube into the opening.

- Holds up to 12 similar size tubes
- Imprinted references for easy identification
- Generously-sized handles for easy transport
- Accepts tubes and vessels from 16–30 mm in size



Cat No	tubes	wells	Price
CLS24330A	15–50 mL	12	11,900.00

Grab Rack®

Rack features silicon grips to hold tubes in place

Turn the rack over and the tubes will not fall out! Great for discarding liquids while tubes stay in place.

- Tubes insert and extract with a gentle push or pull
- ABS rack with silicon gripping fingers. No assembly required
- Not recommended for cleaning in autoclave
- measures 259 x 127 x 49 mm



Cat No	Description	Tubes	Wells	Price
CLS24371A	Grab Rack	16 mL	50	5,517.00

ADAPT-A-RACK®

Capable of holding different size tubes in each available slot as the well literally adapts to your tube size to hold firmly in an upright position. The design enables you to hold 12 mm to 30 mm (5 mL to 50 mL) tubes simultaneously and still allow full viewing of the tube label

- Open sides allow for label viewing
- Flat base to hold tubes upright
- Linked racks should be supported underneath when carrying



available in solid or combination color

Cat No	Rack #1 Color Top/Bottom	Rack #2 Color Top/Bottom	Pack	Price
CLS120184	Blue/Green	Blue/Green	2 Rack/Box	2,512.00
CLS120185	Blue/Yellow	Blue/Yellow	2 Rack/Box	2,512.00
CLS120186	Pink/White	Pink/White	2 Rack/Box	2,512.00
CLS120187	Blue/Blue	Blue/Blue	2 Rack/Box	2,512.00
CLS120188	White/White	White/White	2 Rack/Box	2,512.00
CLS120189	Pink/Pink	Pink/Pink	2 Rack/Box	2,512.00
CLS120190	Yellow/Yellow	Yellow/Yellow	2 Rack/Box	2,512.00
CLS120191	Green/Green	Green/Green	2 Rack/Box	2,512.00

Microtube Racks 80-Well

Easily fill and transfer 80 tubes in compact arrangement (5 x 16 array).

- Polypropylene rack has molded-in coordinates to aid organization
- Measures 225 x 67 x 28 mm. Pack of 5



Cat No	Tubes	Well	Color	Price
CLS29025A	1.5/2.0 mL	80	Natural	4,731.00
CLS29025B	1.5/2.0 mL	80	Blue	4,731.00
CLS29025C	1.5/2.0 mL	80	Green	4,731.00
CLS29025D	1.5/2.0 mL	80	Pink	4,731.00
CLS29025E	1.5/2.0 mL	80	Yellow	4,731.00
CLS29025F	1.5/2.0 mL	80	Orange	4,731.00
CLS29025G	1.5/2.0 mL	80	Assorted Color	4,731.00
CLS29025J	1.5/2.0 mL	80	Red	4,731.00

3 Way Rack

Rack measures 164 x 123 x 62 mm.


- Made of polypropylene
- Finger slots are molded into the rack for easy pickup and moving
- Racks connect to one another on either end
- Autoclavable



Cat No	Tubes	Pack Size	Well	Color	Price
CLS4289B	0.2 mL	6	24	Assorted Color	9,124.00
	0.5 mL	6	14		
	1.5/2.0 mL	6	12		

PCR® Racks Reversible

- PCR® side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Opposite side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes
- Measures 206 x 131 x 54 mm.



Pack Size: 5

Cat No	Tubes	Well	Color	Price
	1.5/2.0 mL (both sides)	12		
CLS2344A	0.2 mL (one side)	168	Assorted Color	11,830.00
	0.5 mL (one side)	40		

Reversible Racks 96-Well

Flip polypropylene rack from 96 wells for 0.5 mL tubes on one side to 96 wells for 1.5/2.0 mL tubes on the other side.



Cat No	Tubes	Pack	Well	Color	Price
CLS2345A	1.5/2.0 mL 0.5 mL	5	96	Assorted Color	12,000.00
CLS2345B	1.5/2.0 mL 0.5 mL	5	96	Natural	12,000.00
CLS2345C	1.5/2.0 mL 0.5 mL	5	96	Blue	12,000.00
CLS2345D	1.5/2.0 mL 0.5 mL	5	96	Green	12,000.00
CLS2345E	1.5/2.0 mL 0.5 mL	5	96	Pink	12,000.00
CLS2345F	1.5/2.0 mL 0.5 mL	5	96	Yellow	12,000.00

Tube Rack 20-well

Interlock multiple racks to increase tube storage. Polypropylene racks lock together on each side.



Pack Size: 6

Cat No	Tubes	Well	Color	Price
CLS4284A	0.5/1.0/1.5/2.0 mL 0.5 mL	20	Blue	9,970.00

PCR® Racks 96-Well

Sturdy polypropylene rack accommodates individual tubes or strips of 8 or 12 tubes.



Pack Size: 5

Cat No	Tubes	Well	Color	Price
CLS23461A	0.2 mL	96	Assorted Color	8,450.00
CLS23461B	0.2 mL	96	Natural	8,450.00
CLS23461C	0.2 mL	96	Blue	8,450.00
CLS23461D	0.2 mL	96	Green	8,450.00

50-Well Microtube Storage Boxes

For smaller storage requirements of 1.5/2.0 mL microtubes or cryogenic vials.

Pack Size: 5



Cat No	Tubes	Pack	Well	Color	Price
CLS120032	1.5/2.0 mL	5	50	Natural	6,051.00
CLS120033	1.5/2.0 mL	5	50	Assorted Color	6,051.00

81-Well Microtube Storage Boxes

Store 1.5/2.0 mL microtubes or cryogenic vials in polypropylene storage boxes, with an easy open, friction fit lid.

- Black light tight version available for light sensitive samples
- Measures 130 x 130 x 47 mm and will fit into most standard freezer racking

Pack Size: 5



Cat No	Tubes	Well	Color	Price
CLS120034	1.5/2.0 mL	81	Natural	8,620.00
CLS120035	1.5/2.0 mL	81	Blue	8,620.00
CLS120036	1.5/2.0 mL	81	Green	8,620.00
CLS120037	1.5/2.0 mL	81	Pink	8,620.00
CLS120038	1.5/2.0 mL	81	Yellow	8,620.00
CLS120039	1.5/2.0 mL	81	Orange	8,620.00
CLS120040	1.5/2.0 mL	81	Assorted Color	8,620.00
CLS120041	1.5/2.0 mL	81	Black	8,620.00

100-Well Microtube Storage Boxes

For larger storage requirements of 1.5/2.0 mL microtubes or cryogenic vials.

- Black, light-tight version available for light sensitive samples

Pack Size: 5



Cat No	Tubes	Well	Color	Price
CLS120042	1.5/2.0 mL	100	Natural	9,901.00
CLS120043	1.5/2.0 mL	100	Black	9,901.00
CLS120044	1.5/2.0 mL	100	Assorted Color	9,901.00
CLS120201	1.5/2.0 mL	100	Blue	9,903.00
CLS120202	1.5/2.0 mL	100	Green	9,903.00
CLS120203	1.5/2.0 mL	100	Yellow	9,903.00
CLS120204	1.5/2.0 mL	100	Orange	9,903.00

Work2Store™ & Work2Store™ Micro Expanding Storage Boxes

Single rack, double function

Reduce sample handling time and the need for several racks when working with samples for cold temperature storage. Patent-pending opening and closing mechanism allows the rack to expand, creating needed room for sample manipulation. Then simply collapse the rack back to its compact form and it is ready for the freezer. Temperature range -80 to 121°C. Two versions now available the original for holding 1.5/2.0 mL cryo and micro tubes and now the new Micro reversible version for holding 0.5 mL and 0.2 mL tubes (see below)

- Instantly creates more working space
- Helps reduce potential cross contamination
- Stops opened microfuge lids from interfering with other wells
- Stackable, autoclavable and with alphanumeric well identifications
- Will fit into standard 2" (53 mm) freezer racking

Work2Store™



Work2Store™ Micro



Work2Store^a

- Clear lid locks the rack together for safe handling and transportation
- Wells have a center depression to stabilize conical tubes and a cryogenic vial lock for one handed cap removal
- Holds 64 x 1.5/2.0 mL cry and/or micro tubes

Work2Store^a Micro

- 0.2 mL tube side has been configured to work 8 strip, 12 strip or single tubes
- Reversible rack which will hold
- 64 x 0.5 mL micro tubes on one side and 96 x 0.2 mL micro tubes on the other side

Cat No	Description	Tubes	Well	Price
CLS120018	Work2Store™	1.5/2.0 mL	64	6,500.00
CLS120183	Work2Store™ Micro	0.5/0.2 mL	64/96	6,500.00

Cardboard Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Store in LN₂ or in mechanical freezers
- Bottom vent holes to allow quicker circulation and removal of cold air, vapor or LN₂
- All box footprints are 13.3 x 13.3 cm, select from two different heights
- Pack Size: 12 Pack



Cat No	Tubes	Color	Height	Price
CLS2860A	Vial Box and Lid	White	50 mm	8,189.00
CLS2860B	Vial Box and Lid	White	75 mm	8,871.00
CLS2860C	Vial Box and Lid	Assorted Color	50 mm	8,800.00
CLS2860CB	Vial Box and Lid	Blue	50 mm	8,800.00
CLS2860CR	Vial Box and Lid	Red	50 mm	8,800.00
CLS2860CY	Vial Box and Lid	Yellow	50 mm	8,800.00
CLS2860E	100-Cell Partition (10 x 10)	Max. tube dia 11.5 mm		6,824.00
CLS2860H	81-Cell Partition (9 x 9)	Max. tube dia 12 mm		6,824.00
CLS2860F	64-Cell Partition (8 x 8)	Max. tube dia 15 mm		6,500.00
CLS2860G	49-Cell Partition (7 x 7)	Max. tube dia 17.5 mm		6,500.00

Cardboard Micro Tube Boxes

- Store in mechanical freezers or vapor phase LN₂*
- Will fit into most standard freezer racking
- 0.2 mL tube (9.5 mm max) version has 169 wells in a 13 x 13 format
- 0.5 mL tube (12 mm max) version has 100 wells in a 10 x 10 format



Cat No	Description	Tubes	Well	Qty	Price
CLS120220	Cardboard Micro Tube Boxes	0.2 mL	169	PK/10	13,520.00
CLS120221	Cardboard Micro Tube Boxes	0.2 mL	100	PK/10	13,520.00

Flatpack Freezer Boxes

These easy self-assemble boxes come flat-packed, ensuring valuable lab space is not taken up storing the boxes before they are needed. Made from corrugated polypropylene, they are more durable than card boxes as they are not susceptible to moisture, so they can be frozen and thawed numerous times without a problem.



- Flat shipped boxes for convenient storage and quick assembly.
- Pre-assembled dividers included with the boxes.
- Store down to temperatures of -80°C
- Holds 81 x 12 mm (max. diameter) tubes in a 9 x 9 format
- Complete with well dividers in packs of 10

Cat No	Description	Tubes	Well	Qty	Price
CLS120262	Corrugated Polypropylene Box with Divider	12 mm	81	PK/10	11,600.00

Foam Tube Rack

Lightweight microtube racks are ideal for long-term studies or tube storage

- Holds 1.5 and 2.0 mL micro tubes. Pack of 10



Cat No	Well	Dimension	Price
CLS4282	100	336 x 95 x 73	13,180.00
CLS4283A	50	210 x 110 x 71	10,500.00

Cardboard Tube Racks and Partitions

- Partitions come in two sizes. Wells measure maximum of 17.5 mm and 29.5 mm dia.
- Economical cardboard racks are disposable
- Write-on surface for descriptions or self numbering systems
- Pack of 10



Cat No	Description	Price
CLS2860M	15 and 50 mL Box and Lid	12,400.00
CLS2860P	50 mL, 16 Cell Partion (4 Columns x 4 Row)	6,190.00
CLS2860T	15 mL, 36 Cell Partion (6 Columns x 6 Row)	7,740.00

Freezer Boxes

freezer boxes

Reusable, economical and moisture resistant plastic

- Made from polypropylene material
- Easy self-assemble boxes shipped lat-packed
- More durable than card boxes as they are not susceptible to moisture
- Significantly reduces the chance of mold
- Dividers included with the boxes
- Store down to temperatures of -80°C

Corrugated Polypropylene 0.2/0.5/1.5-2.0/5 mL

CLS120367



CLS120381



Cat. No.	Products	Place	Pack size	Price
CLS120379	Flatpack Freezer Boxes White, 0.2ml	144 place	10/pk	15,704.00
CLS120380	Flatpack Freezer Boxes Blue, 0.2ml	144 place	10/pk	15,704.00
CLS120381	Flatpack PP Freezer Boxes Red, 0.2ml	144 place	10/pk	15,704.00
CLS120382	Flatpack PP Freezer Boxes Purple, 0.2ml	144 place	10/pk	15,704.00
CLS120383	Flatpack PP Freezer Boxes White, 0.5ml	81 place	10/pk	15,704.00
CLS120384	Flatpack PP Freezer Boxes Blue, 0.5ml	81 place	10/pk	15,704.00
CLS120385	Flatpack PP Freezer Boxes Red, 0.5ml	81 place	10/pk	15,704.00
CLS120386	Flatpack PP Freezer Boxes Purple, 0.5ml	81 place	10/pk	15,704.00
CLS120262	Flatpack PP Freezer Boxes ,White 1.5/2ml	81 place	10/pk	9,980.00
CLS120376	Flatpack PP Freezer Boxes Blue, 1.5/2ml	81 place	10/pk	9,980.00
CLS120377	Flatpack PP Freezer Boxes Red, 1.5/2ml	81 place	10/pk	9,980.00
CLS120378	Flatpack PP Freezer Boxes Purple, 1.5/2ml	81 place	10/pk	9,980.00
CLS120387	Flatpack PP Freezer Boxes White, 5ml	25 place	10/pk	15,704.00
CLS120388	Flatpack PP Freezer Boxes Blue, 5ml	25 place	10/pk	15,704.00
CLS120389	Flatpack PP Freezer Boxes Red, 5ml	25 place	10/pk	15,704.00
CLS120390	Flatpack PP Freezer Boxes Purple, 5ml	25 place	10/pk	15,704.00

Cat. No.	Products	Place	Pack size	Price
CLS120391	Flatpack PP Freezer Boxes White, 15ml	36 place	10/pk	22,000.00
CLS120392	Flatpack PP Freezer Boxes Blue, 15ml	36 place	10/pk	22,000.00
CLS120393	Flatpack PP Freezer Boxes Red, 15ml	36 place	10/pk	22,000.00
CLS120394	Flatpack PP Freezer Boxes Purple, 15ml	36 place	10/pk	22,000.00
CLS120366	Flatpack PP Freezer Boxes White, 50ml	16 place	10/pk	21,100.00
CLS120367	Flatpack PP Freezer Boxes Blue, 50ml	16 place	10/pk	21,100.00
CLS120368	Flatpack PP Freezer Boxes Red, 50ml	16 place	10/pk	21,100.00
CLS120369	Flatpack PP Freezer Boxes Purple, 50ml	16 place	10/pk	21,100.00

Test Tube Rack

5 mL/15 mL/16mm Tube Rack*

This multi-purpose rack is designed with the flexibility to hold 20 - 5 mL, 20 - 15 mL or 16mm tubes OR to hold a40 - 5 mL, 40 - 15 mL or 16 mm tubes. The rack incorporates two moving hinged upper shelves that can be rotated to an up or down position to accommodate the size of the tube being used. In addition, the rack can also hold clinical tubes of a similar diameter (16 mm) and height.

- Rack holds 5 or 15 mL tubes (or tubes up to 16 mm in diameter) or a combination of 20 – 5 mL, 15 mL or 16 mm tubes.
- Unique tear drop shaped openings (top tray) for better retention of tubes
- Ergonomic handles for ease of use
- Imprinted grid on both sides of the upper shelves for quick identification and retrieval
- Shelves can be easily rotated to accommodate the different size tubes



Cat. No.	Product	Wells	Color	Pack Size	Price
CLS120167	Test Tube Rack, 5 mL/15mL	40	Purple Natural	1/pk	7,800.00
CLS120467	Test Tube Rack, 5 mL/15mL	40	Green Blue	1/pk	7,800.00

Trackrack

Take what you need and store the rest

- Segments box made of polypropylene
- Eight segments per box
- 1.5 -2.0 mL cryo and micro-centrifuge tubes (64 tubes)
- A cryogenic vial lock for one-handed cap removal



Cat. No.	Product	Well Dia.	Wells	Color	Pack Size	Price
CLS120356	TracRack, 1.5/2.0 mL	12.3 mm	40	Purple	1/pk	4,754.00

Transformer™ cube

- Polypropylene construction complete with two snap-in grids
- Imprinted grid, first-tube orientation — keyed lid
- Use without grids for larger sized product storage
- Fill with ice and use as a temporary low-temperature workstation
- Clear lid, translucent sides
- Autoclavable



Cat No	Tubes	Pack	Well	Color	Price
CLS23271A	15 mL	1	16	Purple	13,800.00
	50 mL	1	9	Purple	

Trackrack MICRO

Take what you need and store the rest

- Segments box made of polypropylene
- Eight reversible segments per box
- 64 x 0.5 mL tubes, 96 x 0.2 mL micro tubes, or a combination of both



TracRack® Outer Box

- Made from thin film polypropylene material
- Store down to temperatures of -80°C
- Moisture resistant, significantly reducing the chance of mold
- Removable segments allow outer shell to be used for storage of other small laboratory products
- 2 in. (5.3 cm) boxes fit in standard freezer racks

Cat. No.	Product	Well Dia.	Wells	Color	Pack Size	Price
CLS120357	TracRack, 0.2 mL, 0.5 mL	5.1 mm, 7.8 mm	40	Purple	1/pk	4,754.00

Microtube Tube Racks

Reinforced nylon tube racks offering superior chemical and heat resistance over most other plastic based racks.

- Handles provide easy and safe transportation and allow racks to be stacked
- Two formats available for smaller and larger requirements
- Temperature range -40°C to +140°C



Cat No	Tubes	Well	Dimension	Color	Pack	Price
CLS120241	1.5/2.0 mL	24	143 x 90 x 52 mm	Blue	8	13,800.00
CLS120242	1.5/2.0 mL	24	143 x 90 x 52 mm	White	8	15,574.00
CLS120243	1.5/2.0 mL	96	210 x 131 x 54 mm	Blue	4	13,520.00
CLS120244	1.5/2.0 mL	96	210 x 131 x 54 mm	White	4	13,520.00

Test Tube Racks

Tough and more highly resistant to acids, bases, solvents and heat

These durable racks are molded in a single, continuous piece so no assembly required and there are no detachable pieces. Delrin® fiberglass-reinforced polyoxymethylene (POM): Offers the highest chemical and heat resistance (-40 to +140°C)

- Autoclavable, available in half and full rack sizes for 5 tube diameters

Polypropylene:

Offers good chemical and heat resistance (0 to 110°C)

- Not recommended for autoclaving, available in the full size for 2 tube diameters

Both:

- Have molded coordinates for tube orientation
- Designed to be stable when used in water baths
- Have flat labeling area for easier identification
- Available in 5 different colors. Cases of 8



Cat No	Tubes	Well	Dimension	White	Blue	Green	Megenta	Yellow	Price (₹)
Full Size Darlin® Rack Autoclavable	13 mm	72	104 x 202 x 59	CLS27511A	CLS27511B	CLS27511C	CLS27511D	CLS120205	38880
	16 mm	72	127 x 250 x 70	CLS27512A	CLS27512B	CLS27512C	CLS27512D		
	20 mm/15 mL	40	100 x 252 x 83	CLS27513A	CLS27513B	CLS27513C	CLS27513D		38880
	25 mm	40	120 x 300 x 92	CLS27514A	CLS27514B	CLS27514C	CLS27514D		42260
	30/50 mL	24	110 x 282 x 85	CLS27515A	CLS27515B	CLS27515C	CLS27515D	CLS120207	42260
Full Size Darlin® Rack Autoclavable	13 mm	36	104 x 104 x 59	CLS27501A	CLS27501B	CLS27501C	CLS27501D		26130
	16 mm	36	127 x 127 x 70	CLS27502A	CLS27502B	CLS27502C	CLS27502D		27409
	20 mm/15 mL	20	100 x 127 x 83	CLS27503A	CLS27503B	CLS27503C	CLS27503D		29957
	25 mm	16	120 x 122 x 92	CLS27504A	CLS27504B	CLS27504C	CLS27504D		29957
	30 mm/50 mL	9	110 x 110 x 85	CLS27505A	CLS27505B	CLS27505C	CLS27505D		29957
Full Size Polypropylene Rack	13 mm	72	104 x 202 x 59	CLS27551A	CLS27551B	CLS27551C	CLS27551D		25251
	16mm	72	127 x 250 x 70	CLS27552A	CLS27552B	CLS27552C	CLS27552D		25251

Snap-Together Conical Tube Rack

- Shipped flat, easy-to-assemble ABS racks snap together securely
- Designed for stable use in water baths



One-Well 50 mL Connecting Tube Racks

- Molded circular base ring keeps conical or round-bottom tubes in place

Pack Size: 5



Cat No	Tubes	Pack	Well	Color	Price
CLS24306	15 mL		30	Orange	5,271.00
	50 mL		20		

Cat No	Description	Price
CLS23052	1-Well 50 mL Tube Racks, Assorted	8,955.00

One-Well 50 mL Single Tube Racks

Reinforced nylon tube racks offering superior chemical and heat resistance over most other plastic based racks.


- Handles provide easy and safe transportation and allow racks to be stacked
- Two formats available for smaller and larger requirements
- Temperature range -40°C to +140°C



Cat No	Description	Tubes	Well	Color	Price
CLS23051	Epoxy-Coated Wire Rack	50 mL	1	Yellow	3,250.00
CLS23053A	Polystyrene Friction-Fit Rack	50 mL	1	Blue	4,570.00

Fold and Snap Test Tube Rack


- Easy-to-assemble polypropylene racks fold and snap together
- Available in 4 different colors for 3 popular tube sizes



Cat No	Tube	Well	Color	Dimension	Price
CLS243073B	13 mm	90	Blue	245 x 105 x 63 mm	3,000.00
CLS243073W	13 mm	90	White	245 x 105 x 63 mm	3,000.00
CLS243073R	13 mm	90	Red	245 x 105 x 63 mm	3,000.00
CLS243073Y	13 mm	90	Yellow	245 x 105 x 63 mm	3,000.00
CLS243077B	17 mm	60	Blue	245 x 108 x 72 mm	3,000.00
CLS243077W	17 mm	60	White	245 x 108 x 72 mm	3,000.00
CLS243077R	17 mm	60	Red	245 x 108 x 72 mm	3,000.00
CLS243077Y	17 mm	60	Yellow	245 x 108 x 72 mm	3,000.00
CLS243071B	21 mm	40	Blue	245 x 108 x 72 mm	3,000.00
CLS243071W	21 mm	40	White	245 x 108 x 72 mm	3,000.00
CLS243071R	21 mm	40	Red	245 x 108 x 72 mm	3,000.00
CLS243071Y	21 mm	40	Yellow	245 x 108 x 72 mm	3,000.00

Mega Racks

Fully-assembled, polypropylene Mega Racks are great for high throughput labs requiring large sample storage.



Double Mega Rack

Cat No	Description	Size	Wells	Color	Price
CLS24322A	10–13 mm Tube Rack	single	216	Blue	9,200.00
CLS24322E	10–13 mm Tube Rack	single	216	White	9,200.00
CLS24323A	13–16 mm Tube Rack	single	120	Blue	9,200.00
CLS24323E	13–16 mm Tube Rack	single	120	White	9,200.00
CLS120113	10–13 mm Tube Rack	double	432	White	18,210.00
CLS120112	10–13 mm Tube Rack	double	432	Blue	18,210.00
CLS120115	13–16 mm Tube Rack	double	240	White	18,210.00
CLS120114	13–16 mm Tube Rack	double	240	Blue	18,210.00

Pop-Up™ Racks


- Designed to help organize 15 and 50 mL tubes
- Robust polypropylene construction. Not recommended for cleaning in autoclave
- 15 mL max tube dia. 17 mm
- 50 mL max tube dia. 30 mm
- Pack Size: 2



Cat No	Color	Tube	Well	Price
CLS24320B	Blue	15 mL 50 mL	21 12	7,340.00
CLS24320P	Pink	15 mL 50 mL	21 12	7,340.00
CLS24318	Purple	15 mL	45	7,340.00
CLS24319	Green	50 mL	18	7,340.00

Peg Racks


- Spacious peg design allows easy viewing of tube contents
- Use as a drying rack for glass plates or upside-down tubes
- Pegs remain straight even after repeated autoclaving
- Embossed indexing, ribbed grips, and drainage holes
- Pack of 2



Cat No	Tube	Well	Color	Dimension	Price
CLS24311A	14–17 mm	50	White	188 x 104 x 69 mm	5,070.00
CLS24311B	14–17 mm	50	Blue	188 x 104 x 69 mm	5,070.00
CLS24312A	10–13 mm	96	White	82 x 126 x 69 mm	5,600.00
CLS24312B	10–13 mm	96	Blue	82 x 126 x 69 mm	5,600.00


Rota-Rack® Duo Tube Rack

- Rotate and customize to suit your application needs. 4 combinations possible
- Blue section holds 48 x 0.5 mL and 48 x 1.5 mL microfuge tubes on the reverse side
- Purple section holds 96 x 0.2 mL tubes in a standard array pattern one and 48 x 1.5 mL tubes on the reverse side
- Tapered wells securely suspend tubes and makes removal easy
- Alphanumeric well identification



Cat No	Tube	Well	Color	Dimension	Price
CLS29041A	0.2 mL	96	Blue & Purple	296 x 118 x 34	8,600.00
	0.5 mL	48			
	1.5 mL	96			

Rota-Rack® Tube Racks



Cat No	Tube	Well	Color	Dimension	Price
CLS29040A	50 mL	9	Assorted Color	224 x 110 x 110 mm	10,500.00
	15 mL	8			
	20 mm	6			
	10-13 mm	10			
CLS29040B	15 mL	6	Assorted Color	210 x 100 x 100 mm	10,500.00
	1.5/2.0 mL	9			
	0.5 mL	12			
	0.2 mL	32			

Cube Rack



Interlocking multi-tube racks

- Heavy duty polypropylene rack.
- Pack of 5

Cat No	Tube	Well	Color	Price
CLS29050A	50 mL	4	Assorted Color	11,250.00
	15 mL	10		
	12 mL	12		
	1.5/2.0 mL	16		

4-Way Rack

Link together to cater for your application needs

- Pack of 5



Cat No	Color	Price
CLS29022A	Natural	9,630.00
CLS29022B	Blue	9,630.00
CLS29022C	Green	9,630.00
CLS29022D	Pink	9,630.00
CLS29022E	Yellow	9,630.00
CLS29022F	Orange	9,630.00
CLS29022G	Assorted Color	9,630.00

Puzzle Rack

Two-sided puzzle-shaped pieces hold an assortment of 1.5/2.0, 15 and 50 mL test tubes.

- Pack of 4



Cat No	Tube	Well	Color	Price
CLS29060A	15 mL	9	Assorted Color	13,880.00
	50 mL	4		
	1.5/2.0 mL	8 on top side		
	1.5/2.0 mL	4 on opposite side		

Transfer Pipette

These one-piece, unbreakable polyethylene plastic pipettes are well-suited for blood banking (2 mL size), urinalysis, hematology, wet chemistry, microbiology, and serology.

- Bulk packed versions come in handy dispenser box
- Sterile versions available single wrapped, or 20 wrapped at a time
- Lot-to-lot traceability on sterile versions. All in packs of 500

Cat No	Description	Sterile	Length	Graduations	Overall Volume	Price
CLS20635E	0.5 mL, Bulk	No	60 mm	None	1.1 mL	8,000.00
CLS206371A	1 mL, Bulk	No	150 mm	0.25 mL	5 mL	4,500.00
CLS206371B	1 mL, 20/wrap	Yes	150 mm	0.25 mL	5 mL	8,600.00
CLS206371C	1 mL, 1/wrap	Yes	150 mm	0.25 mL	5 mL	12,880.00
CLS206372A	2 mL, Bulk	No	155 mm	0.5 mL	5 mL	4,700.00
CLS206372B	2 mL, 20/wrap	Yes	155 mm	0.5 mL	5 mL	8,600.00
CLS206372C	2 mL, 1/wrap	Yes	155 mm	0.5 mL	5 mL	12,900.00
CLS206373A	3 mL, Bulk	No	155 mm	0.5 mL	7 mL	4,700.00
CLS206373B	3 mL, 20/wrap	Yes	155 mm	0.5 mL	7 mL	9,200.00
CLS206373C	3 mL, 1/wrap	Yes	155 mm	0.5 mL	7 mL	13,500.00



Wire Tube & Bottle Racks


- Chemically-resistant epoxy-coated steel wire racks
- Suitable for conical, round or flat bottom glass or plastic tubes and bottles

Bijou Bottle/Tube Rack




- Rack suitable for Bijou bottles and low profile tubes up to 25 mm dia
- Available in a 2 x 6 or a 5 x 10 format
- Mesh bottom support

Universal Bottle/Tube Rack



- Rack suitable for Universal, McCartney and Scintillation vials and low profile tubes up to 33 mm dia
- Available in a 2 x 6 or a 5 x 10 format
- Mesh bottom support

13 & 16 mm Tube Racks




- For 13 and 16 mm blood collection tubes and test tubes
- Two formats for each size 6 x 12 and 9 x 12, color coded
- Slotted bottom

15 mL/50 mL Tube Racks



- Support capped tubes up to 18 mm dia 3 x 5 format
- Support capped tubes up to 30 mm dia 2 x 4 format
- Slotted bottom

Urine Container Rack



- Ideal for urine containers up to 58 mm dia
- Slotted bottom
- Mesh bottom support

Centrifuge Bottle Rack



- Holds bottles and vessels up to 60 mm dia
- Slotted bottom

Cat No	Description	Pack	Dimension mm	Price
CLS120085	Tube rack, blue, 25 mm well 2 x 6, 12 place		164 x 64 x 32 mm	4,350.00
CLS120086	Tube rack, blue, 25 mm well 5 x 10, 50 place		271 x 145 x 32 mm	7,610.00
CLS120087	Tube rack, blue, 33 mm well 2 x 6, 12 place		212 x 80 x 60 mm	4,350.00
CLS120088	Tube rack, blue, 33 mm well 5 x 10, 50 place		355 x 189 x 60 mm	8,200.00
CLS120089	Tube rack blue, 13 mm well, 6 x 12, 72 place		205 x 105 x 65 mm	6,850.00
CLS120090	Tube rack red, 13 mm well, 9 x 12, 108 place		205 x 162 x 65 mm	6,850.00
CLS23072	Tube rack blue, 16 mm well, 6 x 12, 72 place		237 x 127 x 65 mm	6,850.00
CLS23108	Tube rack red, 16 mm well, 9 x 12, 108 place		241 x 188 x 65 mm	6,850.00
CLS23215	Tube rack blue, 18 mm well, 3 x 5, 15 place		127 x 76 x 76 mm	6,800.00
CLS232461	Tube rack blue, 30 mm well, 2 x 4, 8 place		147 x 80 x 89 mm	6,700.00
CLS120091	Tube/Bottle rack blue, 58 mm well, 6 x 4, 24 place		393 x 262 x 50 mm	13,041.00
CLS23224	Bottle rack blue, 60 mm well, 2 x 3, 6 place		205 x 135 x 104 mm	6,161.00

Floating Foam Tube Racks

- Closed-cell foam resists water absorption so racks can be washed and reused
- Racks include detachable carrying handles-easily lift samples out of baths
- Pack of 5



Cat No	Tubes	Well	Color	Dimension	Price
CLS2165A	50 mL	4	Green	165 x 137 x 65	4,157.00
CLS2165B	15 mL	8	Blue	172 x 57 x 65	3,720.00
CLS2165C	1.5-2.0 mL	24	Yellow	150 x 83 x 69	4,230.00
CLS2166	0.2/0.5/ 1.5-2.0 mL	18	Blue	82 mm dia	4,300.00

Pearl™ Pipettes Adjustable Volume

- Internal metal parts for increased accuracy and calibration stability
- Easy volume adjustment with locking tab to prevent accidental change
- Separate tip ejector and a large digital volume display



Cat No	No of Channels	Volume Range	Price
CLS120097	Single	0.1–2 µL	40954
CLS120098	Single	0.5–10 µL	40954
CLS120099	Single	2–20 µL	40954
CLS120100	Single	5–50 µL	40954
CLS120101	Single	10–100 µL	40954
CLS120102	Single	20–200 µL	40954
CLS120103	Single	30–300 µL	40954
CLS120104	Single	100–1000 µL	40954
CLS120105	Single	500–5000 µL	40954
CLS120106	Single	1000–10000 µL	40954
CLS120215	Eight	5-50 µL	132600
CLS120216	Eight	30-300 µL	132600

Flip & Grip™ Pipette Holder

- Will hold up to 12 single channel pipette, or a combo of 3 single and 3 multichannel
- Holds most common brands. Rubber grips prevent scratching pipette



- Holder rotates for easy pipette selection and storage
- Weighted base. Fixing screw provided for under shelf option
- Measures 332 mm H x 193 mm W (bench top style); 187 mm L x 193 mm W (under shelf style)



Cat No	Color	Price
CLS120026	Blue	20,955.00
CLS120056	Lime	20,955.00
S120057	Magenta	20,955.00

Pipette Stands

These classic 3, 4, and 6-place pipette stands are designed to hold most major pipette brands.

- The upper slots are 25 mm and the lower slots are 20 mm wide
- Acrylic stands are transparent/clear; ABS stands are solid opaque colors



Cat No	Type	Material	Color	Price
CLS206203	Three	Acrylic	Clear	10,750.00
CLS206204	Four	Acrylic	Clear	12,132.00
CLS20620	Six	Acrylic	Clear	15,720.00
CLS20613C	Six	Acrylic	Blue	16,221.00
CLS20613E	Six	ABS	Blue	11,997.00

RF 1000 Pipette Controller

- Tapered silicon pipette adapter allows a large range of pipette sizes and makes to be used. Replaceable when worn
- Accepts common, and readily available replacement filters unlike other pipette aids whose filters are custom. This will help reduce total cost of ownership over time
- Includes a magnetic/wall mountable holder and 0.45 µm filter
- Overall dimensions 180 x 45 x 120 mm.



SPECIFICATIONS:

- Power: Manual
- Filter: PTFE/25 mm dia / 0.45 µm filter in PP housing
- Pipettes: 1 to 50 mL, glass and plastic
- Pipette nozzle: polypropylene, autoclavable
- Dimensions: 180 x 45 x 120 mm

Cat No	Description	Price
CLS1000RFU	RF1000 with magnetic/mountable holder	20046

RF 3000 Pipette Controller

- Operates with 110–240 volt electrical supplies worldwide
- Long-life Li-Ion Battery, with quick charge—up to 20 hours continuous use
- UV resistant housing—longer lasting when exposed to UV sterilizing
- Universal grip and soft coated triggers with direction indicators, for right or left handed users
- Ergonomically balanced in the hand with a loaded pipette
- Bright LEDs indicate three operation modes, High, Low and Gravity, toggled by repeatedly pressing a large mode button that is in easy reach of your thumb when in use
- Five different colored nose cones included to help differentiate department, task or user, as well as a labeling area on the side of the unit
- Filters are readily available, so you are not tied into a pipette specific brand—helps reduced the total life costs of ownership verses other makes of pipette aids
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes
- Low battery indicator; can be recharged when in use



Operational mode indicator



Magnetic wall mount



Specifications:

- Power: Rechargeable Li-Ion battery (replaceable)
- Filter: PTFE/25 mm dia 0.45 µm filter
- Pipettes: 1 to 100 mL, glass and plastic
- Pipettes nozzle: Polypropylene, autoclavable
- Dimensions: 15 x 152 x 41 mm

Rf3000 Package includes:

- Charging stand and wall-mountable/magnetic holder
- Nozzle adapter for use with small pipettes
- Li-Ion rechargeable battery
- Three replacement 0.45 µm filters
- Five nose cones, Blue, Lime, Shamrock, Teal, Lilac
- AC adapter 100 to 240V, 50/60Hz with four interchangeable plugs



Cat No	Description	Price
CLS3000	RF3000 with Power Adapter	78931
Replacement Parts:		
CLS5000RFFIL	Filters, Pack of 5	6,947.00
CLS300501	Universal Power Adapter	15,626.00
CLS3000RFRWH	Wall Mount Holder	6,874.00

HS300502	Battery RF3000	7,946.00
CLS5000RFRP	Pipette Adapter, Silicone	4,779.00
CLS300504	Nose Cones, Pack of 5	15,012.00

Vacuum/Pressure Bulb

- Pack of 2



Cat No	Capacity	Price
CLS20631B	85 mL	2,354.00

Multi-Channel Reagent Reservoirs



- Tapered reservoir bottoms into a narrow trough to aid full sample recovery

Cat No	Volume	Material	Sterility	Qrt	Price
CLS20521A	55 mL	Clear PVC	Non-Sterile	100/pk	11,456.00
CLS20521B	55 mL	White Polystyrene	Non-Sterile	100/pk	11,456.00
CLS20521C	55 mL	White Polystyrene	Sterile	50/pk	21,121.00
CLS20721B	100 mL	White Polystyrene	Non-Sterile	100/pk	14,159.00
CLS20721C	100 mL	White Polystyrene	Sterile	50/pk	31,427.00

Dual Solution® Reservoir

- Non-sterile. Pack of 25 reservoirs and 25 lids



Cat No	Description	Price
CLS20821A	Dual Solution Reservoir	20,515.00

Pipette Racks

Angled four shelved compartments keep your 0.1 to 50 mL pipettes or cans within easy reach while preventing damage and clutter.

- Heavy duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- Magnetized ABS version has strong magnets on both sides and mounts securely to most metal surfaces*
- All racks measure 291 x 88 x 401 mm




Cat No	Description	Material	Color	Price
CLS20611	Pipette Rack	Acrylic	Clear	31,240.00
CLS20612	Pipette Rack	ABS	Blue	16,558.00
CLS20612M	Magnetized Pipette Rack	ABS	Blue	17,910.00
CLS20613	Pipette Rack	ABS	Red	16,560.00

Cargo™ Pipette Rack

Organize the Cargo™ Pipette Rack to meet your laboratory or personal needs. Click together horizontally or vertically to build the desired configuration, or use as individual holders, utilizing the three lid sections provided.

- Magnets on both sides of units connect securely and safely to most metal surfaces
- Angled shelves store 0.1 to 50 mL pipettes within arm's reach for use
- Protective lid encloses the unit to prevent dust build up and provides additional shelf space
- Link to additional sets for larger needs and requirements
- Assorted shades of blue: light, medium, dark. Pack of 3



Cat No	Color	Dimension	Price
CLS20615M	Blue	301 x 100 x 127.5	17,840.00

Pipette Tip Disposal Box

Convenient for one-handed disposal of used pipette tips. Eliminates any need to handle pipette tips used with hazardous fluids.

- The 3/8" (10 mm) thick acrylic provides protection from Beta radiation
- The large, hinged lid allows easy disposal of contents
- Suggest using in conjunction with small Biohazard Disposal Bags (page 71), if hazardous fluids are not acrylic friendly
- Measures 154 x 171 x 206 mm.




Cat No	Description	Price
CLS234634	Pipette Tip Disposal Box	28,662.00

Rubber Pipette Filler

This rubber pipette filler is designed for efficient one-handed operation. Control air evacuation, liquid uptake, and liquid dispensing with three valves, each containing a chemical-resistant stainless steel ball.

- Pinch release valve for smooth liquid flow
- Filler its standard size pipettes
- Order the red pipette filler for safe pipetting of caustics. It does not contain carbon black and will not smudge
- Filler measures 54 x 54 x 130 mm



Cat No	Description	Price
CLS20632A	Pipette Filler, Black	1,636.00
CLS20632B	Pipette Filler, Red	1,636.00
CLS20632C	Pipette Filler, Blue	1,636.00

Safety Bulb Pipette Filler

- Pack of 3



Cat No	Capacity	OD	Price
CLS20623	15 mL	41 mm	3,620.00

Latex (Natural Rubber) Dropper Bulbs

- Pack of 72



Cat No	Capacity	Length	OD	Hole Die	Price
CLS20622A	1 mL	43 mm	13 mm	5 mm	6,450.00
CLS20622B	2 mL	52 mm	13 mm	6 mm	6,450.00

Cardboard Slide Tray



Cat No	Capacity	Trim/imprinting	Dimension	Qty.	Price
CLS9902	20-Place	Black	208x341x9	72/cs	88322
CLS9913	20-Place	Blue	208x341x9	72/cs	88322
CLS9914	20-Place	Green	208x341x9	72/cs	88322
CLS120196	20-Place	Red	208x341x9	72/cs	88322
CLS120197	20-Place	Yellow	208x341x9	72/cs	88322
CLS9930	30-Place	Black	281x328x10	40/cs	85067

Cardboard Slide Mailers

1-Place



Cat No	Capacity	Dimension	Qty	Price
CLS9903	1-Place	95 x 42 x 5 mm	25/pk	5,050.00
CLS9904	2-Place	103 x 79 x 5 mm	36/pk	10,731.00

Plastic 5-Slide Mailers

Pack of 25



Cat No	Description	Color	Dimension	Price
CLS15982	Side-Opening Mailer	Natural	226 x 340 x 125	3,100.00
CLS15983G	End-Opening Mailer	*		3,100.00
CLS15986	End-Opening Mailer	Natural		3,041.00

Rubber Bulbs

• Pack of 1



Cat No	Capacity	Length	OD	Hole Die	Price (₹)
CLS20630B	15 mL	56 mm	34 mm	4 mm	534.00
CLS20630C	30 mL	66 mm	46 mm	6 mm	671.00
CLS20630D	60 mL	80 mm	54 mm	6 mm	790.00
CLS20630E	150 mL	103 mm	73 mm	6 mm	2,354.00

20-Place Slide Tray Racks

Search, identify, and easily locate the specimen slide you need!


- ABS rack holds five 20-place slide trays (not included see above)
- Eliminates the worry of breakage during retrieval-easy access to each of the trays
- Matching clasps provide area for identification labels and doubles as a handle



Cat No	Capacity	Color	Dimension	Price
CLS99040A	5 Slide Trays	Black	226 x 340 x 125	9,193.00

Plastic 2-Place Slide Mailer

- Measures 93 x 68 x 5 mm.
- Pack of 25



* Assorted Color

Cat No	Description	Color	Dimension	Price
CLS15983A	2-Place Slide Mailer	*	226 x 340 x 125	3,030.00
CLS15983B	2-Place Slide Mailer	Natural		3,030.00

Slide Staining Tray Set

- Pack of 4 base trays and 1 lid
- Stain, rinse, and dry slides on a single working tray
- Made from a polypropylene and polyethylene blend
- A clean and fresh working surface area every time
- Up to eight slides it comfortably onto the base and the deep well holds up to 38 mL
- Compact size and recessed handles for easy transport from work area to sink
- Two convenient pour spouts for quick disposal of unwanted liquid waste
- Dark lid protects slides for light sensitive applications
- Pack of 4 base trays and 1 lid
- Disposable




Cat No	Description	Dimension	Price (₹)
CLS15951A	8 slides	305 x 127 x 29.5 mm	9,062.00

Slide Holder

Hold slides during evaluation, cleaning, and drying

- Made from chemical-resistant polyethylene
- Solid, strong construction with rubber feet for bench top stability



Cat No	Description	Dimension	Price
CLS15981C	15	153 x 102 x 21 mm	9,327.00
CLS15981A	40	298 x 105 x 24 mm	10,037.00
CLS15981B	78	299 x 178 x 21 mm	12,506.00

Tungsten-Point Lab Pen


Cat No	Description	Price
CLS1846	Tungsten-Point Lab Pen	9,196.00



- Ideal for forensic, life science, agricultural, industrial and safety investigations to name a few
- Base measures 89 x 127 mm. Bendable neck and head measure 559 mm

Cat No	Description	Price (₹)
CLS15941B	Flexible Digital Microscope	76222

Handheld Digital Microscope



- Use with a computer or projector
- Uploads and charges via USB. Software compatible with Microsoft Windows® (98SE/2000/ME/XP/Vista/7), Mac OS® X 10.4. Computer with monitor required for use (not included)
- Up to 43x magnification using a 17" monitor (up to 53x using a 21" monitor)
- High-quality resolution with 640 x 480 pixels (1280 x 960 Mac). Direct x 9.0 or above required for Windows OS, VGA CMOS sensor
- 15 frames per second with auto white balance and auto exposure 2 head adapters for viewing different media (included)
- 4 LED lights illuminate dark settings or objects
- Removable pieces for easy cleaning. Simple plug and play functionality
- Measures 44 mm dia; 64 mm H

Cat No	Description	Price
CLS15941A	Handheld Digital Microscope	38106

Microscope Slide Dispenser



Cat No	Description	Price
CLS15984	Slide Dispenser	4,021.00

Lens Tissue Books

Non-abrasive cleaning for delicate optical surfaces Lint-free i ber paper is ideal for cleaning lens and optics.

- Pack consists of 20 books of 100 sheets. (4 inner packs of 5 books)



Cat No	Description	Price
CLS15920	Lens Tissue Books	36,361.00

Flexible Digital Microscope

Easy-to-use microscope connects to any PC and provides high-quality images from 10x to 200x magnification.

- Bendable gooseneck positions easily to view hard-to-see objects
- Built-in camera and microphone provide immediate assessment
- Use with a projector or PC. Charges and uploads via computer USB port
- Features true 1.3 MP resolution, 15 fps frame rate, focus wheel, 6 LED lights and measurement software (compatible with Windows XP®, Windows Vista®, Mac OS® 10.6 and above). Direct X 9.0 required
- Computer with monitor required for use (not included)
- Portable and ideal for investigating l at or 3D subject matter especially items often difficult to see round or behind



Tall Slide Box 100-Place

Sturdy ABS material for rugged use. Slide boxes are designed for storage and transport.

- Grooved slots hold individual slides in place
- Soft cork lining ensures that slides will not slip
- Identification card, lining and grooved slots are numbered
- Stackable with other Tall boxes and/or 100-place Premium slide boxes



Cat No	Capacity	Color	Dimension	Price
CLS15979A	100	Blue	208 x 175 x 60 mm	4,905.00

Standard Microscope Slide Boxes

Constructed of durable polypropylene, slide boxes hold either 25 or 100 standard slides (25 x 75 mm and 1 x 3").

- Stacking feet maximize stability
- Slide Identification card. No liner
- 25-place box measures 141 x 92 x 3 mm; 100-place box measures 210 x 169 x 37 mm



Cat No	Capacity	Color	Price
CLS15990A	25-Place	Blue	946.00
CLS15990B	25-Place	Red	946.00
CLS15990C	25-Place	White	946.00
CLS15991A	100-Place	Blue	1,284.00
CLS15991B	100-Place	Red	1,284.00
CLS15991C	100-Place	White	1,284.00




25-Place Slide Boxes

100-Place Slide Boxes

Color	Cat (cork-lined)	Price	Cat (cork-lined)	Price
Blue	CLS15994A	1,877.00	CLS15994M	2,215.00
Green	CLS15994B	1,877.00	-	-
Red	CLS15994C	1,877.00	CLS15994P	2,215.00
Yellow	CLS15994D	1,877.00	-	-
White	CLS15994E	1,877.00	CLS15994R	2,215.00
Grey	CLS15994F	1,877.00	-	-
Black	CLS15994G	1,877.00	-	-

True North® Microscope Slide Box



Cat No	Color	Slides	Price
CLS15988H	Purple	100	4,905.00

50-Place Slide Boxes

Color	Cat (cork-lined)	Price
Blue	CLS15996A	1,791.00
Red	CLS15996B	1,791.00
White	CLS15996C	1,791.00

25-Place Slide Boxes

Color	Cat (cork-lined)	Price (₹)	Cat (cork-lined)	Price
Blue	CLS15989A	1,623.00	CLS15989M	1,690.00
Green	CLS15989B	1,623.00	-	-
Red	CLS15989C	1,623.00	CLS15989P	1,690.00
White	CLS15989E	1,623.00	CLS15989R	1,690.00

Premium Microscope Slide Boxes

- Molded, rigid boxes securely hold slides in grooved slots
- Durable ABS plastic boxes are available with either cork or foam lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet
- [Pack size 1](#)

50-Place Slide Boxes



100-Place Slide Boxes

Storage Rack for Microscope Slide Boxes



Cat No	Bins	Dimension	Price
CLS159890B	15	711 x 171 x 228 mm	81,103.00

Sample Collection Cart

- Made of durable stainless steel and aluminum construction
- Handle for pulling the unit
- Painted surfaces utilize Microban® protection
- 5 wheel, weighted base for stability of unit
- Designed to hold the Droplet™ blood collection station
- Adjustable height – Range: 30 tp 55 inches (762 x 1397 mm)
- Designated tray area dimensions 16 x 9 x 6 inches (429 x 265 x 172 mm)
- Working surface dimensions 14.25 x 7.5 inches (362 x 191 mm)



Cat No	Dimension	Price
CLS120263	Sample Collection Cart 1	2,57,400.00
CLS120264	Deluxe Sample Collection Cart 2	3,58,280.00

Blood Collection Tray

Lightweight, roomy and easy-to-clean Droplet™ Blood Collection Tray. Keeping a clean and organized tray has never been easier. Wash the ABS tray as needed, and replace the insert periodically for a fresh start. Tray inserts come in two configurations and include an extra-deep compartment for a sharps container on the following page.



- Tray features a supply drawer, side pockets, and a ratcheting handle with soft rubber grip that rotates out of the way for easy access
- Each tray includes either a 13 mm or a 16 mm Heathrow OneRack® test tube rack and two disposable polystyrene inserts
- Additional inserts sold in packs of 4 in two different configurations

Cat No	Dimension	Price
CLS2200A	Droplet Blood Collection Tray with 2 Inserts (Style A) and a 36-Place, 13 mm Test Tube Rack	24,297.00
CLS2200B	Droplet Blood Collection Tray with 2 Inserts (Style B) and a 36-Place, 16 mm Test Tube Rack	24,297.00

Droplet™ Lid

Protect the contents of the Droplet Tray with this sturdy lid. Custom designed for a snug fit, the lid provides access to the tray handle in its vertical position.

- Can be secured with cable ties for extra security
- Made of polystyrene



Cat No	Color	Price
CLS2202	White	13,552.00

Disposable Droplet inserts



Insert A for 13 mm tube rack

Insert B for 16 mm tube rack

Cat No	Color	Price
CLS2201A	Insert A for 13 mm Tube Rack, Pack of 4	7,097.00
CLS2201B	Insert B for 16 mm Tube Rack, Pack of 4	7,097.00

Sharps Containers

Provides safe disposal of syringes, needles and other small sharps. Stackable to reduce storage space before use



- 1.0 L capacity (Quart) with a clear lid to view fill level
- Slide access lid for temporary and final closures for safety and security
- Built-in needle removal port, preventing direct contact with sharps
- Measures 115 W x 115 D x 192 H mm, packs of 18

Cat No	Description	Color	Price
CLS120177	1.0 L sharps container	Yellow	18,210.00
CLS120178	1.0 L sharps container	Red	18,210.00

X-tourn™ Tourniquets

High-quality disposable tourniquets that will not rope or curl when being used, so you can rely on the performance of every one. Latex free, no rubber proteins, resins or fatty acids that can produce hypersensitive reactions.



- Made from supple and smooth polyisoprene
- Less abrasive than other brands, uniform thickness prevents shredding and breaking
- Each measures 2.54 x 45.72 cm

Cat No	Capacity	Color	Pack	Price
CLS120108	Tourniquets	Green	100/pk	8,271.00
CLS120109	Tourniquets	Green	25/roll	3,860.00
CLS120110	Tourniquets	Orange	250/pk	19,300.00
CLS120111	Tourniquets	Orange	25/roll	3,860.00

ToteMax™ Blood Collection Trays

Unique phlebotomy tray design and patented hook feature allows the tray to be on or over most bed rail, footboard and handrails outside a typical patient room. Made from food grade resins for a durable non-breakable design



Cat No	Color	Price
CLS120021	Blue/Red/Pink/Purple/Green	23,160.00

Blood Collection Tray

- Customize the compartments to suit your personal needs with dividers provided.
- Lightweight easy-to-clean ABS construction
- Comfort grip handle tested to hold up to 10 Kg!




Cat No	Description	Tube Rack	Price
CLS120261	Blood Collection Tray, White	17 mm rack	36,956.00

Droplet® Sample Storage Tray

Hard wearing storage trays that can be used for specimen containers with the insert provided, or for general lab products without the insert.


- Lightweight and sturdy to hold up to 24 full specimen containers (up to 58 mm in diameter). Specimen containers not included.
- Removable insert to allow outer shell to be used for storage or transportation of laboratory products



Cat No	Color	Price
CLS120172	Droplet® Sample Storage Tray	15,532.00

Urine Container Rack

- Chemically-resistant epoxy-coated steel wire racks
- Suitable for 1 at bottom containers up to 58 mm dia
- Slotted bottom
- Mesh bottom support



Cat No	Description	Dimension	Price
CLS120091	Tube/Bottle rack blue, 58 mm well, 6 x 4, 24-place	393 x 262 x 50 mm	13,041.00

Bouffant Head Covers

Latex-free bouffant head covers, designed for roominess and comfort. Gentle elastic provides secure coverage without binding. Made of comfortable SPP (spunbond polypropylene).



Cat No	Description	Color	Pack	Price
CLS120058	Bouffant Head Cover 24"	Blue	Pk/100	8,247.00

Shoe Cover

Slip resistant, Latex Free
Made of strong, low-linting polypropylene, slip-resistant tread



Cat No	Description	Color	Pack	Price
CLS120059	Shoe Cover, Slip Resistant, PP	Blue	Pk/100	17,087.00

Clamp-Ette® Pipette Holders

- Chemical-resistant polypropylene and its most major pipette brands.
- Available in a single or triple format



Cat No	Description	Color	Price
CLS2062A	Clamp-Ette	Blue	3,839.00
CLS20162A	Tri Clamp-Ette	Blue	8,450.00

Benchtop Counters

Mechanical, benchtop counters are perfect for blood cells (label included), bacterial colonies or other repetitive counting.

- Each key records up to 999
- Total window keeps track of total number of strokes made on all other keys; bell rings when total reaches 100



Blood cell key labels


- CLS5360A: Baso., Eosin., Mono., Lymph., Poly., and Total, measures 250 x 76 x 56 mm
- CLS5360B: Baso, Eosin., Myelo., Juven., Stab., Seg., Lymph., Mono., and Total, measures 320 x 76 x 56 mm

Cat No	Description	Dimension	Price
CLS5360A	6-Key Benchtop Counter	250 x 76 x 56	81,696.00
CLS5360B	9-Key Benchtop Counter	320 x 76 x 56	118,335.00

Arctic Squares® Freezer Rack

Inovative horizontal drawer system for upright freezers

Constructed of corrosion-resistant stainless steel with four drawers to hold 16 boxes.



- Drawers feature a unique stop, that when fully extended pivot down allowing for easy box access
- Polycarbonate drawer handles provide a user-friendly grip. Includes a slot and label for drawer/sample identification

Cat No	Orientation	Capacity	Shelves	Box Height	Dimension	Price
CLS2864	Horizontal	16 Boxes	4	53 mm	139 x 571 x 238	78,709.00

Assay Plate Freezer Rack

For Assay Plate and PCR® rack storage in upright freezers



Conig	Cat No	Open (Dimension)	Overall (Dimension)	Price
4 x 4	CLS28624	135 x 111 x 56 mm	543 x 112 x 230 mm	41,886.00

Freezer Racks

Upright Freezer Rack Configurations



Upright Freezer Rack

Chest Freezer Rack



Chest Freezer Racks

Conig	Cat No	Box Size	Dimension	Price
7 x 1	CLS2862AB	2" (53 mm)	140 x 141 x 390	22,618.00
10 x 1	CLS2862A	2" (53 mm)	140 x 141 x 557	21,654.00
11 x 1	CLS2862AC	2" (53 mm)	140 x 141 x 613	27,975.00
12 x 1	CLS2862AD	2" (53 mm)	140 x 141 x 669	29,761.00
13 x 1	CLS2862AE	2" (53 mm)	140 x 141 x 724	30,952.00
5 x 1	CLS2862DA	3" (75 mm)	140 x 141 x 406	20,833.00
8 x 1	CLS2862D	3" (75 mm)	140 x 141 x 649	20,341.00
9 x 1	CLS2862DB	3" (75 mm)	140 x 141 x 730	27,975.00

Upright Freezer Racks 2" boxes

Conig	Cat No	Box Size	Dimension	Price
4 x 3	CLS2862HA	2" (53 mm)	418 x 140 x 225	35,559.00
4 x 4	CLS2862H	2" (53 mm)	556 x 140 x 225	37,848.00
4 x 5	CLS2862HD	2" (53 mm)	695 x 140 x 225	51,784.00
5 x 3	CLS2862HE	2" (53 mm)	418 x 140 x 280	37,499.00
5 x 4	CLS2862HB	2" (53 mm)	556 x 140 x 280	51,784.00
7 x 4	CLS2862HC	2" (53 mm)	556 x 140 x 392	60,713.00

Upright Freezer Racks 3" boxes

Conig	Cat No	Box Size	Dimension	Price (₹)
3 x 3	CLS2862KA	3" (75 mm)	418 x 140 x 245	34,523.00
3 x 4	CLS2862K	3" (75 mm)	556 x 140 x 245	36,746.00
4 x 4	CLS2862KB	3" (75 mm)	556 x 140 x 326	49,999.00
5 x 4	CLS2862KC	3" (75 mm)	556 x 140 x 407	57,141.00

Cube Freezer Rack - Upright

For 15 and 50 mL tube box storage
Constructed of corrosion-resistant stainless steel with 6 (2 x 3) rack spaces.



Conig	Cat No	Box Size	Dimension	Price
2 x 3	CLS28622	152 x 157 x 134	459 x 159 x 271	40,839.00

Cube Freezer Rack - Chest

For 15 and 50 mL tube box storage
Constructed of corrosion-resistant stainless steel with 4 (4 x 1) rack spaces.



Conig	Cat No	Box Size	Dimension	Price
4 x 1	CLS28621	152 x 156 x 133	154 x 157 x 539	19,372.00

Slide Box Freezer Rack

For slide boxes, Reversible PCR® and 50-Well EPS racks in upright freezers



Conig	Cat No	Box Size	Dimension	Price
3 x 3	CLS28623	180 x 212 x 73	541 x 217 x 221	56,022.00

Cryogenic Vial Holders

- Well tabs allow one-handed opening and closing of vials that have a slot feature
- Embossed index numbers and squared, upper left corner provide easy orientation
- Ribbed grips make handling with gloves easier



Pack of 4

Cat No	Wells	Color	Price
CLS23229B	50	Blue	21,121.00

Arctic Squares® Polycarbonate Cryostorage Boxes




Safely store vials from -196 to +121°C in mechanical freezers as well as liquid nitrogen! Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease.

- Hinged lid with built-in stop allows easy one-handed access to samples (not on 25 place)
- Forward-sloped base and high contrast indexing on the transparent lid ensure quick visual orientation
- Magnetic tipped forceps supplied with styles B, D, E and F. Help pick vials and remove boxes from freezer racking
- All boxes have vent and drainage holes, are stackable and autoclavable

Cat No	Description	Base Color	Vial Size	l x m x h ml	Qty	Price
CLS23325A	5 x 5 Array Boxes	Ruby Red	up to 2.0	76 x 76 x 53	8/pk	14,484.00
CLS23381B	9 x 9 Array Standard Boxes	Assorted Color	up to 2.0	133 x 133 x 53	4/pk	11,845.00
CLS23381C	9 x 9 Array Tall Boxes	Amethyst Purple	up to 5.0	133 x 133 x 96	5/pk	18,443.00
CLS23300D	10 x 10 Array Boxes	Arctic Blue	up to 5.0	133 x 133 x 53	4/pk	11,845.00

PVC Cryogenic Cane Sleeves

- Clear plastic sleeves protect vials in aluminum canes during handling and storage
- Safe for cryogenic freezing
- More durable than cardboard sleeves
- Transparency allows sample identification without removing sleeve
- Measure 14.5 mm dia x 275 mm long. Pack of 100



Cat No	Description	Price
CLS23210	Clear Cane Sleeves	17,572.00

Cryogenic Vial Canes



Aluminum canes hold vials for freezer storage or immersion in L₂N dewars.


- For use in liquid and vapor phase; use approved vials only
- Canes bend to fit all major cryogenic vial brands. Pack of 12

Cat No	Length	No of Vials	Vials Size	Price
CLS23205	290 mm	5	up to 2.0 mL	5,576.00
CLS23206	300 mm	6	up to 5.0 mL	5,576.00

Mini Coolers

Ideal for keeping PCR® reagents, enzymes, and other biochemical reagents cool while working at the bench. Filled with insulating gel, the Mini Coolers can also be used for controlled sample freezing and for longer temperature retention samples.

- Polycarbonate construction with rubber non-slip feet
- Gridded and keyed one-position lids for easy sample identification
- 0.5 mL tube adapters supplied as standard with the 0.5–2.0 mL tube size versions
- Temperature-retention graph supplied with each of the 7 styles
- Versions available for 0.5 and 1.5/2.0 mL microcentrifuge tubes and cryogenic vials and a 96-place version for 0.2 mL tubes and assay plates



Cat No	Temp	Tube Array	No. Tubes	Tube Sizes	Color	Lid	Dimensions (mm)	Temp Retention	Price
CLS120070	0°C	3 x 4	12	0.5–2 mL	Green	Clear	145 x 105 x 100	1°C for 4 hr**	35,493.00
CLS120071	0°C	4 x 8	32	0.5–2 mL	Green	Clear	245 x 105 x 100	1°C for 5 hr**	60,065.00
CLS120072	0°C	4 x 8	32	0.5–2 mL	Green	Gel Filled	245 x 105 x 100	1°C for 7.5 hr**	65,526.00
CLS120073	-20°C	3 x 4	12	0.5–2 mL	Blue	Clear	145 x 105 x 100	-15°C for 2.5 hr**	35,493.00
CLS120074	-20°C	4 x 8	32	0.5–2 mL	Blue	Clear	245 x 105 x 100	-15°C for 3.5 hr**	60,065.00
CLS120075	-20°C	4 x 8	32	0.5–2 mL	Blue	Gel Filled	245 x 105 x 100	-15°C for 5 hr**	65,526.00
CLS120076	-20°C	8 x 12	96	0.2 mL	Blue	Clear	225 x 165 x 110	-15°C for 3.5 hr	65,526.00

Cool Containers

Conveniently nest lid beneath bucket or pan

- The close-it lid keeps evaporation to a minimum and features a comfortable recessed knob
- The lid and bottom of the bucket/pan are keyed to it together
- For low-temperature use to -196°C at 1 bar
- Taper design for easy stacking
- Precise pour spout and generous handles
- In four, color-coded, user-friendly sizes



Cat No	Shape	Size	Thickness	Qty	Price
CLS28714K	Bucket	4 L	181 (H) x 234 (D)	Black Onyx	24,410.00
CLS28721G	Pan	1 L	267 x 195.5 x 114	Emerald	17,780.00
CLS28724B	Pan	4 L	381 x 254 x 140	Arctic Blue	26,976.00
CLS28729R	Pan	9 L	521 x 292 x 163	Ruby	36,509.00

Freezer Thermometer

Protect frozen and chilled samples from potential spoilage. Monitor temperatures inside your freezer or refrigerator with the added convenience of sound and flashing LED lights to alert out-of-range temperatures.

- Suction cup attaches the probe inside freezer/refrigerator while the thermometer secures to the outside with hook-and loop tape
- Large digital display shows freezer/refrigerator temperature plus ambient temperature
- Set freezer range: 14.9° to -9.4°F (-9.5° to -23°C); set refrigerator range: 35.6° to 46.4°F (2° to 8°C)
- Switches toggle for freezer or refrigerator monitoring and °F or °C
- Includes one AAA battery
- Unit measures 136 x 19 x 28 mm; probe cable measures 0.9 m; probe measures 5 mm OD x 17 mm L



Cat No	Color	Price
CLS28630	Freezer Thermometer	8565

smartSPATULAS®

Ideal for keeping PCR® reagents, enzymes, and other biochemical reagents cool while working at the bench. Filled with insulating gel, the Mini Coolers can also be used for controlled sample freezing and for longer temperature retention samples.



- Polycarbonate construction with rubber non-slip feet
- Gridded and keyed one-position lids for easy sample identification
- 0.5 mL tube adapters supplied as standard with the 0.5–2.0 mL tube size versions
- Temperature-retention graph supplied with each of the 7 styles
- Versions available for 0.5 and 1.5/2.0 mL microcentrifuge tubes and cryogenic vials and a 96-place version for 0.2 mL tubes and assay plates



smartSPATULAS® autoclavable

Cat No	Description	use with tubes	Qty.	Price
CLS120005	Micro Antistatic	0.2 mL	300/pk	7,171.00
CLS120008	Eco	15 and 1.5 mL	300/pk	10,748.00
CLS120003	Standard Blue	15 and 1.5 MI	300/pk	10,124.00
CLS120004	Standard Opaque	15 and 1.5 mL	300/pk	10,124.00
CLS120007	Std. Opaque Sterile	15 and 1.5 MI	300/pk	40,954.00
CLS120006	Macro Green	50 mL	150/pk	11,226.00

smartSPATULAS® autoclavable

Cat No	Description, mm	l, mm	Finish	Price
CLS15906	Flat end, 30 L; spoon, 12 L x 4 W	150	Glazed	1,014.00
CLS15907	Round lat end, 50 L x 8 W Tapered Flat end, 45 L x 6.5 W	195	Mirror	1,014.00
CLS15909	Flat end, 50 L; spoon, 38 L x 15 W	225	Mirror	1,081.00

Forceps

Disposable forceps feature corrugated handles and serrated tips. ABS. Pack of 6.



Cat No	Length	Tip	Color	Price
CLS15998C	111 mm	Straight	Dark Blue	4021


Weighing Dishes

Pack of 100



Color	Size	Dimension mm	Regular Material	Price
Aluminum	Small	43 dia x 12 H	CLS14521A	1,654.00
	Medium	57 dia x 14 H	CLS14521B	1,860.00
	Large	70 dia x 16 H	CLS14521C	2,549.00
G Aluminum	Large	101 dia x 8 H	CLS14522	3,238.00

Weighing Boats



- A. Standard Weighing Boats
- B. Diamond-Shaped Weighing Boats
- C. Hexagonal Weighing Boats
- D. Pour-Boat Weighing Dishes
- E. Square Weighing Boats

Pack of 500

Color	Size	Dimension mm	regular material	Price (₹)	antistatic material	Price
White Blue Black	Small	46 x 46 x 8	CLS1420A	3281.00	CLS1420AA	5,750.00
	Medium	85 x 85 x 24	CLS1420B	5256.00	CLS1420BB	9,000.00
	Large	140 x 140 x 22	CLS1420C	8988.00	CLS1420CC	14,434.00
B White Black	5 mL	55 x 35 x 6	CLS1424AA	2724.00	CLS1425AA	4,375.00
	25 mL	71 x 46 x 14	CLS1424A	4505.00	CLS1425A	7,751.00
	30 mL	78 x 56 x 14	CLS1424B	4505.00	CLS1425B	7,751.00
	100 mL	119 x 90 x 19	CLS1424C	9417.00	CLS1425C	15,812.00
C White	Small	38 x 25 x 9	N/A		CLS14251A	3,686.00
	Medium	64 x 51 x 19	N/A		CLS14251B	6,545.00
	Large	102 x 76 x 22	N/A		CLS14251C	18,086.00
	X-Large	127 x 89 x 31	N/A		CLS14251D	23,081.00
D White	Small	43 x 58 x 13	CLS1419A	8,612.00	CLS120225	9,157.00
	Medium	83 x 132 x 26	CLS1419B	17,879.00	CLS120226	18,528.00
	Large	108 x 183 x 26	CLS1419C	27,525.00	CLS120227	28,647.00
E White	Medium	90 x 90 x 7	CLS1422	9,129.00	N/A	

Wash Bottles Safety Labeled Self Venting

- Wide bottle mouth makes filling easier
- Non-metallic vent that prevents the risk of solvent leakage
- With the exception of the bottle designed for Sodium Hypochlorite, which is a white LDPE bottle for light protection and has a HDPE closure, the bottles are made from LDPE, with a PP closure
- Six named chemicals, a self labeling** and assorted pack of bottles to choose from
- 500 mL bottle size, supplied in packs of 6



Cat No	Description	Color	Price
CLS120250	Acetone	Red	10,910.00
CLS120251	Distilled Water	Natural	10,910.00
CLS120252	Isopropanol	Yellow	10,910.00
CLS120253	Ethanol	White	10,910.00
CLS120254	Methanol	Green	10,910.00
CLS120255	Sodium Hypochlorite	Blue	10,910.00
CLS120256	Self labeling**	Natural	10,910.00
CLS120257	Assorted*		10,910.00

Wash Bottles Color-Coded

- Wide bottle mouth makes filling easier
- Non-metallic vent that prevents the risk of solvent leakage
- With the exception of the bottle designed for Sodium Hypochlorite, which is a white LDPE bottle for light protection and has a HDPE closure, the bottles are made from LDPE, with a PP closure
- Six named chemicals, a self labeling** and assorted pack of bottles to choose from
- 500 mL bottle size, supplied in packs of 6



Cat No	Color	Price
CLS120245	White	6,670.00
CLS120246	Red	6,670.00
CLS120247	Yellow	6,670.00
CLS120248	Green	6,670.00
CLS120249	Blue	6,670.00

Weighing Papers

- Nitrogen-free for a very smooth surface and good sample recovery
- Strong enough to support dense material
- Available in 3 sizes



Cat No	Dimensions	Qty.	Price
CLS120116	3 x 3" (76.2 x 76.2 mm)	500/pk	7,688.00
CLS120117	4 x 4" (101.6 x 101.6 mm)	500/pk	11,348.00
CLS120118	6 x 6" (152.4 x 152.4 mm)	500/pk	22,729.00

Eyewash Stations

- Frame is equipped with top and bottom holes for easy mounting to any lat surface. Screws provided
- Each station includes attached mirror, instruction cards, and one or two 1 L (32 oz) dousing bottles



Cat No	Description	Bottles	h x w x d, mm	Price
CLS1019B	Single Bottle Eyewash Station	One 32 oz (1 L)	425 x 167 x 105	6,725.00
CLS1020B	Double Bottle Eyewash Station	Two 32 oz (1 L)	425 x 355 x 105	10,780.00
CLS1021B	Replacement Bottle	32 oz (1 L)	86 dia x 275 H	2,197.00

Transport Box

- Polycarbonate construction and available in three colors for user/owner identification
- Silicone water-tight seal and three lid clasps provide secure closure
- Carry handle folds neatly into body space allowing the units to be stacked
- Lid opens a full 180°, allowing total access to contents and making it easier to clean
- Use the separators provided to create your own compartment sizes, or use empty to transport larger products and equipment
- Autoclavable and measures 380 x 196 x 160 mm



Cat No	Color	Price
CLS120052	Clear with Blue handles	23,160.00
CLS120077	Red with Red handles	23,160.00
CLS120078	Yellow with Yellow handle	23,160.00

Grabbit™ Temp Mitts

- Oversized surfaces with nubs for secure grip
- Mitts withstand temperatures from -70 to 500°F (-57 to 260°C)



Cat No	Description	Dimension, mm	Price
CLS1034A	Full-Hand Mitt	177 x 100 x 32	4,596.00

High-Temperature Gloves

- Maintain flexibility and comfort in most applications
- Two sizes available; wrist length and forearm length
- Machine washable



Cat No	Gauntlet	Overall Length	Price
CLS1031B	250 mm	485 mm	9,631.00

Bottle Carrier

- Made from 304 grade Stainless Steel, with a wire mesh bottom for increased stability
- Handles fold down for easier storage
- Autoclavable and measures 370 x 180 x 380 mm with handles extended



Cat No	Dimension	Price
CLS120092	CLS 120092 370 x 180 x 380 (with handles extended)	13,584.00

Laboratory Safety Razor Blade

- Cleanly slice agarose gels* right on the light box without scratching the surface
- Ideal for cutting soft tissue and plant matter
- Adjust the size required by simply cutting with scissors!
- Supplied in packs of 100



Cat No	Dimension	Price
CLS120259	19 mm H x 38 mm W x 0.22 mm thick	6,516.00
CLS120260	25 mm H x 50 mm W x 0.22 mm thick	7,337.00

Eco-smartFunnel™ Paper Funnel

- Cleanly slice agarose gels* right on the light box without scratching the surface
- Ideal for cutting soft tissue and plant matter
- Adjust the size required by simply cutting with scissors!
- Supplied in packs of 100



Cat No	Dimension	Price
CLS120219	Eco-smartFunnel™ 186 mL cap, Top Dia. 6 cm2, Bottom Dia. 0.64 cm2, Height 12.8 cm	10,192.00

Lead Rings

- Maintain flexibility and comfort in most applications
- Two sizes available; wrist length and forearm length
- Machine washable



Cat No	Type	lask sizes	id, mm	od, mm	weight, g	Price
CLS8882A	O-Shape	125 to 500 mL	48	77	277	6,252.00
CLS8882B	O-Shape	250 to 1000 mL	51	90	508	6,928.00
CLS8882C	O-Shape	500 to 2000 mL	57	98	635	8,279.00
CLS8882D	O-Shape	1000 to 4000 mL	70	115	970	9,124.00
CLS8881A	C-Shape	125 to 500 mL	42	75	227	7,265.00
CLS8881B	C-Shape	250 to 1000 mL	51	90	472	8,279.00
CLS8881C	C-Shape	500 to 2000 mL	54	98	599	9,124.00
CLS8881D	C-Shape	1000 to 4000 mL	66	113	907	9,800.00

Clearly Safe® Safety Glasses Dispensers

Ergonomically angled for easy dispensing Translucent acrylic dispensers provide easy access to safety accessories.

- Hinged lid protects glasses from dust and debris
- Set on counter or wall mount (screws included)
- Large Dispenser holds approx. 20 glasses, Small Dispenser approx. 10
- Great for holding other personal items, too!



Cat No	Description	l x w x d, mm	Color	Price
CLS1041	Small Safety Glass Dispenser	227 x 162 x 156	Blue	19,695.00
CLS1040A	Large Safety Glass Dispenser	203 x 203 x 406	Blue	35,655.00

Biohazard Disposal Bags

- Epoxy coated steel wire holder available for bags HS1003D & HS10320
- Autoclavable to 250°F (121°C)



Cat No	Description	Thickness	Dimension	Price
CLS10320	Bag, Case of 500	1.57	203 x 305	22776.00
CLS10321	Bag, Case of 200	1.57	356 x 483	20576.00
CLS10322	Bag, Case of 200	1.57	483 x 584	31788.00
CLS10323	Bag, Case of 200	1.57	635 x 838	76066.00
CLS10044H	Bag Holder			5372.00

Autoclavable Bags

Clear bags imprinted with blue biohazard labeling designed for safe disposal.

- 4 inner packs of 50 bags, sold as a case of 200




Cat No	w x l, mm	Price
CLS1002A	310 x 660	22,071.00
CLS1002B	415 x 600	24,523.00
CLS1002C	610 x 810	45,980.00

Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or lat tubing from 0.01 to 0.8 mm thick. Sliding blade removes excess material 19 mm from the seal.

- Timing knob ensures safe sealing. Simply turn knob to set time (1-8 seconds), press handle down and remove when light indicates cycle is complete
- Heating element measures 310 mm long and is only hot when sealing



Cat No	Description	Plug Type	Price
CLS28271	220/230 VAC	Cont. Europe/Israel	53,560.00

e Count Colony Counter

eCount™ combines a multi-function electronic counter with a Sharpie® pen for marking to prevent missing or double-counting colonies.

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED lash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different
- Lock counter of to allow use of the marker without registering more counts



Cat No	Description	Price
CLS120000	Electronic Colony Counter Package	82,929.00

Sharpie® Marker Pens

- Dual-Tip Colored Sharpie® Pens
- Bright selection of colors offers instant visual recognition
- Combines a fine tip at one end and an ultrafine tip at the other end, perfect for a larger range of labeling requirements



Standard Black Sharpie® Pens

- Strong fine tip will keep its shape much longer than other brands of pens
- Long lasting, supplied in packs of 12

Cat No	Description	Pack	Price
CLS15094	Sharpie® Dual-Tip Pen Set (Assorted)	8/pk	7,970.00
CLS15092ZP	Sharpie® Standard Fine Pens (Black)	12/pk	8785.00

Mini -Light Box & Mini-Magnifier

Mini Light Box saves bench space and can be stored in a drawer.

- Viewing area (10 x 127 mm) is perfect for quick screening of 96-well plates, Petri dishes, electrophoresis gels, slides, and transparencies
- Powered by four AA batteries (included), Measures 161 x 135 x 44 mm
- Mini Magnifier enlarges objects for easier viewing
- Swiveling, molded lens. 1.75X magnification
- Lens measures 185 x 100 mm; stand measures 207 mm H




Cat No	Description	Price
CLS15093	Mini Light Box, With Batteries	35,708.00
CLS15092D	Mini Magnifier	25,189.00

smartPicks® Inoculation Picks

Available in three lengths and two needle sizes to suit a wide range of application needs. The Mini its neatly into a 1.5 mL microfuge tube allowing the lid to be closed for full sample transfer/recovery.

- Naturally biodegradable, made from 100% USA grown White Birch wood
- Manufactured in an FDA-registered factory and comply with FDA standards for food contact
- Autoclavable and available in semi or needle point styles for delicate operations



Cat No	Description	Qty	Price
CLS120009	Mini, Semi Point, 1.28", Fits 1.5 mL tube (closed)	2000/pk	21,577.00
CLS120010	Mini, Needle Point, Fits 1.5 mL tube (closed)	2000/pk	21,577.00
CLS120011	Standard, Semi Point, 2.5" Single Tip	1500/pk	23,638.00
CLS120012	Standard, Needle Point, 2.5" Single Tip	1500/pk	23,638.00
CLS120013	Large, Semi Point, 6" Single Tip	500/pk	15,333.00
CLS120014	Large, Needle Point, 6" Single Tip	500/pk	15,333.00

Inoculating Needle/Loops

Transfer bacterial cultures without risk of contamination.

- Use loop end to uniformly streak a gel surface
- Use needle end for removal of colonies
- Choose between flexible polypropylene or the more rigid ABS versions
- Disposable, Sterile, 20 per zip bag with lot-to-lot traceability
- Measures 200 mm L x 4 mm dia. Pack of 1000



Cat No	Description	Material	Price
CLS81121A	1 µL Loop	Flexible PP	15,333.00
CLS81121B	1 µL Loop	Rigid ABS	15,333.00
CLS81121C	10 µL Loop	Flexible PP	15,333.00
CLS81121D	10 µL Loop	Rigid ABS	15,333.00

Cell Spreaders

Delta™



L-Shape



Stainless Steel



Delta™: Flexible rippled handle and domed spreading bar to help reduce agar tears. Gamma sterilized

L-Shape: Lot-to-lot traceability, 10 wrap comes in a zip bag

Stainless Steel: Reusable. Flame and Autoclave sterilizable

Cat No	Description	Sterile	bar x dia, mm	Qty	Price
CLS8151	Delta™	Yes, single wrap	172 x 37 x 4	500/pk	20,242.00
CLS8171A	L-Shape	Yes, single wrap	146 x 38 x 4	500/pk	26,364.00
CLS8171B	L-Shape	Yes, 10 Wrap	146 x 38 x 4	500/pk	20,232.00
CLS86655	Stainless Steel	No	203 x 33 x 4	each	2,826.00

Petri Dish Dispenser

- Sturdy acrylic holds up to 30 x 100 mm dishes in two easy-access partitions
- Bench top or wall mountable (screws provided)
- Measures 222 x 107 x 236 mm



Cat No	Description	Price
CLS23452	Acrylic Petri Dish Dispenser	17,914.00

Wire Petri Dish Holder


- Holds up to 6 x 100 mm Petri Dishes
- Epoxy-coated steel wire
- Measures 127 x 86 x 107 mm



Cat No	Description	Price
CLS23455	Wire Petri Dish Holder	4,562.00

Droplet® Sample Storage Tray

- Large, usable, easy to clean surface area enables products of all shapes and sizes to be contained
- Stackable with other trays support as many as four high
- Integrated label holder and cards to enable easily changed identification
- Lightweight and sturdy insert to hold up to 24 full specimen containers (up to 58 mm in diameter). Specimen containers not included



Cat No	Description	Price
CLS120172	Droplet® Sample Storage Tray	15,533.00

Drawer Organizer

Polystyrene organizer is perfect for storing and organizing small laboratory supplies.

- Protect small and fragile items
- Fit in most standard drawers or store on benches and shelves
- Dimensions: 460 x 307 x 63 mm



Cat No	Description	Price
CLS2538A	3-Compartment Tray: 2 compartments measure 430 x 80 mm 1 compartment measures 430 x 110 mm	10,306.00

MagLab™ Series



Create lab storage space-of your bench top

Easily organize common of your bench top with color and features four most metal surfaces, or create your own holding point with the metal plate provided.

the MagLab™ storage system. Each rack is made of robust ABS plastic in a bright translucent color and features four powerful magnets for a secure hold onto most metal surfaces, or create your own holding point with the metal plate provided.

- MagRack™ holds four 50 mL tubes and nine 15 mL tubes
- MagWipe™ holds a standard (110 x 210 mm) Kimwipes® box
- MagPette™ holds two pipettors (its most popular styles)
- MagBox™ holds miscellaneous accessories
- Use the MagLab™ self-adhesive mounting plates for use on non-metallic surfaces
- The MagLab™ Series is available individually or as a pack of four

Cat No	Description	Color	l x w x h, mm	Price
CLS23500	MagLab Series, 4/pk	Teal/Green/Purple/Orange	522 x 393 x 102	15,476.00
CLS23501	MagRack	Teal	131 x 112 x 64	4,464.00
CLS23502	MagWipe	Green	130 x 91 x 98	4,464.00
CLS23503	MagPette	Purple	131 x 112 x 76	4,464.00
CLS23504	MagBox	Orange	130 x 78 x 102	4,464.00

Clearly Safe® Workstation Storage Bin

The essential combination accessory holder
Organize your most commonly used laboratory disposables in one central location.

- Acrylic unit has room to store gloves, Kimwipes®, and Parafilm®
- Use bins to store tubes, tips, etc.
- Measures 265 x 137 x 440 mm




Cat No	Description	Price
CLS234523	Workstation Storage Bin	35,486.00

Benchtop Dispensing Bins

Work great with gloved hands, bright colors!

Ergonomically designed with large, angled openings and oversized lip.

- Lid remains open when tilted back
- Nonskid rubber feet and pre-drilled holes for wall mounting
- Green bin includes a removable divider
- Mounting screws included



Cat No	Color	Compartments	w x h x d, mm	Price
CLS23411	Neon Green	Dual	180 x 155 x 165	17,231.00
CLS23412	Neon Blue	Single	155 x 155 x 170	17,742.00

Pipette Workstation

Keep pipette and disposable tips together in one convenient stand.

- Acrylic rack holds up to four pipettors on top shelf
- Dual bins house two different sizes of tips. Self-closing lids provide dust-free storage
- Oversized lip on lid allows for easy access even when wearing gloves
- Measures 244 x 165 x 246 mm



Cat No	Description	Price
CLS234631	Pipette Workstation	30,414.00

Workstation Dispensing Bins

See inventory levels of pipette tips, tubes, fittings, finger cots, and other bulk small items at a glance.

- Flip-top covers on generous openings make filling and dispensing easy and keep products dust-free
- 4 sizes available to cater to a wide variety of storage needs
- Includes non-skid rubber feet



Cat No	No Bins	l x w x d, mm	Price
CLS23401	One	180 x 175 x 235	16,380.00
CLS23402	Two	139 x 160 x 368	19,108.00
CLS23403	Three	210 x 155 x 212	21,837.00
CLS23404	Four	278 x 158 x 217	23,202.00

Modular Glove Box Holders

The Anti-Microbial version (blue) is ideal for use in clean rooms, hospitals, and microbiological work areas. The Microban additive prevents growth of molds and bacteria.

- High-impact, durable resin is resistant to harsh cleaning solutions
- Modular design allows single boxes or horizontal stacking for multiple box holder
- Universal fit for most major brands of glove boxes
- Mounting screws included
- Each unit measures 254 x 156 x 97 mm. Pack of 3



Cat No	Description	Color	Price
CLS120258	Modular Glove Box Holder	White	18,210.00
CLS23491A	Modular Glove Box Holder, Anti-Microbial	Blue	23,818.00

Glove Box Holders

Free up benchtop space by mounting holders on wall for easy dispensing of gloves.

- Clear acrylic holders include a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves
- Sturdy, epoxy-coated steel wire racks are an economical solution
- All holders fit most major brands of boxes
- Mounting screws included



Cat No	Description	Material	w x h x d, mm	Price
CLS23471A	1-Box Holder	Acrylic	141 x 260 x 110	11,089.00
CLS23471B	2-Box Holder	Acrylic	295 x 260 x 110	16,890.00
CLS23471C	3-Box Holder	Acrylic	295 x 385 x 110	19,961.00
CLS23457	1-Box Holder	Wire	140 x 209 x 108	5,576.00
CLS23458	3-Box Holder	Wire	250 x 454 x 108	11,152.00

Wall Holder 4 Position Adjustable

Not just for glove boxes! Ideal for gloves, face masks, shoe covers... etc. The adjustable holding points provide true flexibility to hold up to four different boxed protective products and other lab essentials (size dependent)

- Ideal for holding different boxed sizes of gloves, mop caps, ear plugs, face masks, shoe covers etc.
- Adjustable universal fit, easy to clean, with a low profile
- Wall- or door-mountable, mounting hardware included
- Measures 152 x 762 x 76 mm



Cat No	Description	Color	Price
CLS120107	Adjustable Wall Holder	152 x 762 x 76 mm	37,497.00

Lab Wipe Holders

Clear acrylic holders in boxes of small or large lab wipes.

- Mount on a shelf or wall
- The smaller version attaches with double-stick tape and the larger version with screws
- Mounting materials included




Cat No	Description	l x w x d, mm	Price
CLS234521	Small Holder	118 x 122 x 84	4,095.00
CLS234522	Large Holder	320 x 127 x 97	8,530.00

Storage Racks for 80-Well Racks

These racks help you keep large numbers of microtubes organized on your benchtop.

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0 mL microtubes
- Plenty of room to slide the racks in and out



Cat No	Description	l x w x d, mm	Price
CLS23263A	Vertical Rack	157 x 236 x 234	44,144.00
CLS23263B	Horizontal Rack	310 x 236 x 122	44,144.00

Floating Microtube Racks

- Detachable legs provide two options: Attach legs so rack stands upright on bench top, or detach legs so samples can be quickly popped out by pushing down onto any surface.
- Not recommended for cleaning in autoclave.
- Measure 60 mm high including legs.
- Pack of 4.



Cat No	Shape	Tubes	Well	Color	Price
CLS2132A	Square (102 mm)	1.5/2.0 mL	16	White	6,694.00
CLS2134A	Round (98 mm dia)	1.5/2.0 mL	20	White	6,694.00

Soft Cover Dispenser

Translucent acrylic dispensers provide easy access to safety accessories.

- Hinged lid protects contents from dust and debris
- Set on counter or wall mount (screws included)
- Great for holding other personal items, too!



Cat No	Description	w x h x d, mm	Color	Price
CLS1042C	Soft Cover Dispenser	227 x 162 x 156	Natural	48,898.00

Earplug Dispenser

Translucent acrylic dispensers provide easy access to safety accessories.

- Holds approximately 200 pairs of foam earplugs
- Set on counter or wall mount (screws included)



Cat No	Description	w x h x d, mm	Color	Price
CLS1040B	Earplug dispenser	227 x 162 x 156	Green	35,655.00

Three-in-One Dispenser

Three separate lift-up lids make loading and dispensing supplies easy.

- Middle compartment has reach-in hole for dispensing
- Durable acrylic is see-through so product levels are always visible



Cat No	Description	w x h x d, mm	Color	Price
CLS1042	Three-in-One dispenser	406 x 406 x 203	Grey	48,898.00

Tubby® Container

Sturdy storage tub with lid and divider

- The lid snaps in place to secure the contents and for safe transport



Cat No	Description	Qty	w x h x d, mm	Price
CLS23453	Tubby	Pk/5	330 x 200 x 115	8,955.00

Wire Drying Rack

Saves glassware from chips and cracks
Drain and dry test tubes, beakers, flasks, and other labware.

- Holds up to 90 pieces of labware on 18 x 127 mm pins and 72 x 76 mm pins
- Constructed of spot-welded, vinyl-coated steel wire for durability
- Easily accessible and removable bottom drainage pan



Cat No	Description	Price
CLS23243A	Wire Drying Rack	28,886.00

Foil Rolls


Designed for easy use and minimal waste

- Ideal for autoclaving, freezing, incubating, protecting labware, lining, sealing, and general storage



Cat No	Description	Roll Size	Price
CLS23534A	102 mm (4-inch) length Foil Roll	152.4 m	3,372.00
CLS23536A	153 mm (6-inch) length Foil Roll	152.4 m	4,904.00

Test Tube Baskets



Ideal for washing and storing test tubes
Strong baskets for use with mechanical washers and autoclaves.

- Two versions available, epoxy-coated steel, or the lighter weight aluminum style
- Will withstand repeated autoclaving and use in the lab

	Cat No	Material	l x w x h, mm	Price
A	CLS20341A	Aluminum	152 x 152 x 152	6,590.00
B	CLS20341B	Aluminum	256 x 152 x 152	7,401.00
C	CLS20341C	Aluminum	332 x 228 x 178	9,631.00
D	CLS120082	Coated steel	125 x 125 x 125	4,893.00
E	CIS120083	Coated steel	150 x 150 x 150	5,977.00
F	CLS120084	Coated steel	332 x 228 x 178	9,237.00

Label Station

Organize your labels for easy accessibility

- Acrylic station includes a writing surface, a cutting edge for easy dispensing, and non-skid rubber feet
- Counter-sunk holes provided for wall mounting



Cat No	Description	Color	l x w x h, mm	Price
CLS23460	Label Station	Blue	280 x 120 x 89	21,307.00

ABS Dispenser for Parafilm® Sealing Film

Unique design for maximum safety
Store, dispense, and cut Parafilm®, tape, labels, and Tough-Tags®.

- ABS plastic dispenser designed with a safety razor that slides and cuts with no danger to fingers. Supplied with two easily exchanged cutting blades



Cat No	Description	Color	w x h x d, mm	Price
CLS234525B	ABS Dispenser	Blue	120 x 156 x 171	12,165.00
CLS234525C	ABS Dispenser	Green	120 x 156 x 171	12,165.00
CLS234525D	ABS Dispenser	Red	120 x 156 x 171	12,165.00
CLS234525R	Replacement blades for ABS dispenser, pack of 5			1,585.00

Parafilm® Sealing Film

Parafilm® M is a semi-transparent, flexible, thermoplastic film with two-way stretch.

- Will stretch 200% and adhere to irregular shapes and surfaces to provide a secure seal



Cat No	Film Width	Roll Length	w x h x d, mm	Price
CLS234526A	50 mm (2-inch)	75 m	135 x 58 x 132	9,188.00
CLS234526B	100 mm (4-inch)	38 m	132 x 135 x 112	9,498.00
CLS234526C	100 mm (4-inch)	75 m	135 x 107 x 132	16,506.00

DuraSEAL™

Durable heat-resistant, solvent-proof Laboratory stretch ilm

- Heat-resistant and microwave compatible for temperatures up to 100°C (boiling water)
- Highly resistant to solvents and caustic agents
- Film can stretch 500–800% before breaking or tearing
- Forms an airtight barrier and is twice as long as some other brands!



Stretches 500-800%



Cat No	Size	Ideal Use	Price
CLS-1-150	25 mm x 46 m/1.0" x 150'	Tubes/Vials	8,944.00
CLS-2-500	50 mm x 152 m/2.0" x 500'	Beakers	14,144.00
CLS-4-250	100 mm x 76 m/4.0" x 250'	Jars/Bottles	19,968.00
CLS-4-500	100 mm x 152 m/4.0" x 500'	Jars/Bottles	26,000.00
CLS-5-150	127 mm x 46 m/5.0" x 46'	Large Containers	17,264.00

PetriSEAL™ and ContainerSEAL™

Sealing tape for Petri Dishes, Microplates, lids and caps

- Seal remains flexible even with intermittent exposure to 100°C
- Tape stretches, conforms and seals to irregular perimeters
- Waterproof and chemically resistant to most solvent and caustic agents



Cat No	Size	Description	Price
CLS-SEAL-108B	13 mm x 33 m/0.5" x 108'	PetriSEAL Blue	15,600.00
CLS-SEAL-108C	13 mm x 33 m/0.5" x 108'	PetriSEAL Clear	15,600.00
CLS-SEAL-108R	13 mm x 33 m/0.5" x 108'	PetriSEAL Red	15,600.00
CLS-SEAL-108Y	13 mm x 33 m/0.5" x 108'	PetriSEAL Yellow	15,600.00
CLS-SEAL-108W	13 mm x 33 m/0.5" x 108'	PetriSEAL White	15,600.00
CLS-SEAL-75R	8.5 mm x 33 m/0.75" x 108'	ContainerSEAL Red	20,800.00
CLS-SEAL-75W	8.5 mm x 33 m/0.75" x 108'	ContainerSEAL White	20,800.00
CLS-SEAL-75Y	8.5 mm x 33 m/0.75" x 108'	ContainerSEAL Yellow	20,800.00
CLS-SEAL-XLR	25 mm x 33 m/1.0" x 108'	ContainerSEAL Red XL	20,800.00
CLS-SEAL-XLW	25 mm x 33 m/1.0" x 108'	ContainerSEAL White XL	20,800.00
CLS-SEAL-XLY	25 mm x 33 m/1.0" x 108'	ContainerSEAL Yellow XL	20,800.00

Mini Bin™ Foil Dispenser

Safely cut foil to any size

- Compact to fit in limited work space areas and easy to transport
- Safe to use—load the foil roll, press the cover to activate the safety guard
- Foam strip prevents foil retracting back into unit
- For use with 102 mm or 153 mm (4 or 6 inch) length foil rolls (sold separately)



Cat No	w x h x d, mm	Price
CLS23520A	190 x 100 x 112	17,857.00

Laboratory Notebooks

For recording research data, experimental procedures, and design ideas, use these notebooks to keep information safe and well organized.

- Numbered, acid-free, pages include spaces for title, project number, book number, date, and witness signature
- Lined notebooks feature college-ruled pages for easy report writing
- Covers for all versions are waterproof and chemical-resistant
- All notebooks contain periodic and conversion tables in the back
- Measures 218 x 281 mm



Cat No	Description	Page Format	Color	Qty	Price
CLS8610A	100-Page Notebook	Grid	Black	12/cs	74,090.00
CLS8610B	200-Page Notebook	Grid	Blue	12/cs	92,591.00
CLS8610D	200-Page Notebook	Lined	Red	12/cs	90,136.00

Time Tracker® Mini

Based on the Original Time Tracker, the Mini is an economical option intended for easy use.

- Total alarm time ranges from 5 to 120 minutes in 5-minute increments
- Visual indicator features lighted dome. Green light during count down; yellow for the designated "warning" time; red when time has expired
- Volume control and audible warning alarm
- Requires 3 AAA batteries, not included. Measures 83 x 83 x 121 mm



Cat No	Description	Price
CLS24653	Time Tracker Mini	12,870.00

Time Tracker®

Makes it easy to monitor your process from a distance! Designed like a traffic signal so the meaning of each light is intuitive. Green light means less than half done; yellow means more than half done, and red means finished.

- Programmable time
- Six different audible alerts. With LCD display
- LCD with 8 mm high characters, 10 characters across. Operates on 110 VAC adapter (included)
- Measures 112 x 118 x 230 mm



Cat No	Description	Price
CLS24651	Time Tracker	21,311.00

Heat Resistant Gloves

Terry cloth gloves provide durability and long wear with heat and abrasion resistance

- Applications include: general use for handling utensils, handling items when autoclaving, light/medium fabrication, metal and glass handling, etc.
- Gloves are heat resistant when dry (do not use when wet, near steam, with sharp objects, or with LN2)
- Made of 100% knitted cotton terry cloth that is washable
- Glove is designed with a straight thumb allowing more comfort and flexibility with heat barrier
- Gauntlets provide extra safety for wrist and forearm protection



Specifications

- Overall Length: Approximately 480 mm (18.9 in)
- Width (Palm): Approximately 135 mm (5.3 in)

Cat No	Description	Size	Color	Qty	Price
CLS120483	Heat Resistant Gloves	Universal	Brown	2 gloves per pack	9,630.00

smartBOATS

Made from FSC® Certified (FSC® Mix) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

A truly green solution to sample weighing.

- Made of grease proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Control contents like a polystyrene weigh boat
- An alternative to regular weighing papers
- Available in a small and medium in both brown for light colored samples and white for dark samples



Cat. No.	Size	L x W x H	Color	Pack Size	Price
CLS120312	Small	35 x 35 x 15.8 mm	Brown	500/pk	10,977.00
CLS120313	Medium	47.6 x 47.6 x 14.2 mm	Brown	500/pk	12,693.00
CLS120314	Small	35 x 35 x 15.8 mm	White	500/pk	10,977.00
CLS120315	Medium	47.6 x 47.6 x 14.2 mm	White	500/pk	12,693.00



- A. Pocket Timer/Stopwatch/Clock**
 - Four independent timing channels with visual and audible alarms
 - Use as a stopwatch, a count up/down timer, or a clock
 - Includes spring clip and magnetic back
- B. Triple Timer/Clock**
 - Three channels to count up/down with
 - Simultaneous display, time-out, and 60 second audible alarm
 - Channels may be started independently or simultaneously
 - Timer automatically counts up after counting down to zero
 - Memory recall
- C. Digital Sand Timer/Clock**
 - Watch the digital sand gradually disappear over time just like mechanical sand timers
 - Audible and visual alarms
 - Flip one way to display count up/down timer, flip the other way to display the clock in 12- or 24-hour mode
 - Memory recall
- D. Big-Digit Dual Timer/Clock**
 - Display clock and timer simultaneously, or two independent timers
 - Use one timer to count up, or use two timers to count down independently
 - Audible and visual alarm
- E. Multi-Function Stopwatch/Clock**
 - Quick start, quick stop timer features time, split, hourly chime, calendar, alarm
 - 100 cm neck cord
- F. Four Channel Timer/Clock**
 - Independently or simultaneously set time on four channels
 - Audible and visual alarms signify elapsed time
 - Easily toggle between channels and clock with top buttons
 - Includes spring clip with magnetic back

Cat No	No of Channels	Timing Capacity	Resolution	Count Direction	Display Digits, h	Clock	Alarm	Battery	l x h x d, mm	Price	
A	CLS24490D	Four	100 Hr.	1/100 Sec.	Up/down	6 mm	Yes	Yes	One LR44 1.5V	61 x 71 x 25	5,742.00
B	CLS24490M	Three	20 Hr.	1/100 Sec.	Down	10.5 mm	Yes	Yes	One AAA	86 x 66 x 22	7,570.00
C	CLS24671	One	99 Min.	1 Sec.	Up/down	15 mm	Yes	Yes	Three LR44	58 x 114 x 35	5,820.00
D	CLS24780	Two	24 Hr.	1/100 Sec.	Up/down	25.5 mm	Yes	Yes	One AAA	142 x 97 x 20	9,630.00
E	CLS24645	One	60 Min.	1/100 Sec.	Up	8.5 mm	Yes	Yes	One LR44 1.5V	87 x 78 x 21	7,941.00
F	CLS24670	Four	99 Hr.	1/100 Sec.	Down	14 mm	Yes	Yes	Two LR44	73 x 59 x 13	5,819.00



Code	Product Description	Label	Pk Size	Price	Code	Product Description	Label	Pk Size	Price
CLS-122 BLACK	Solvent Resistant Pen, black	Laboratory Marker	5	3790.00	CLS-CLTT-3000	Clear Laser Tough-Tags 1.05 x 0.50" 1,920/pk	Labels	5	50050.00
CLS-130 RED	Solvent Resistant Pen, red	Laboratory Marker	5	3790.00	CLS-CMBO-1000	Laser Cryo-Babies/Tough-Spot Set	Labels	5	46200.00
CLS-GD-110760	Gel Drying Frames 14 x 14cm	Gel Related	5	81713.50	CLS-COLB-1000	Color Coded Tough-Tags, Blue 500/roll	Labels	5	41085.00
CLS-GD-110761	Gel Drying Frames 24 x 24cm	Gel Related	5	160688.00	CLS-COLG-1000	Color Coded Tough-Tags, Green 500/roll	Labels	5	41085.00
CLS-C-110762	Cellophane 14 x 14cm-50/pk	Gel Related	5	10955.00	CLS-COLR-1000	Color Coded Tough-Tags, Red 500/roll	Labels	5	41085.00
CLS-C-110763	Cellophane 24 x 24cm-50/pk	Gel Related	5	15520.00	CLS-COLW-1000	Color Coded Tough-Tags, White 500/roll	Labels	5	41085.00
CLS-TP-1000	Autoclave Indicator Tape 0.75" x 2160"	Tape	5	11165.00	CLS-COLY-1000	Color Coded Tough-Tags, Yellow 500/roll	Labels	5	41085.00
CLS-TP-2000	Autoclave Indicator Tape 1.00" x 2160"	Tape	5	13475.00	CLS-CRTH-1081	TT Cryo-Babies 0.81 x 0.28" 2,000/roll	Labels	5	62254.50
CLS-TT-1000	Chameleon Tough-Tags® 0.94" x 0.50" 1,000/Roll	Labels	5	34650.00	CLS-CRTH-1094	TT Cryo-Babies 0.94" x 0.50" 2,000/roll	Labels	5	62254.50
CLS-TT-2000	Chameleon Tough-Tags® 1.28" x 0.50" 1,000/Roll	Labels	5	42350.00	CLS-CRTH-1128	TT Cryo-Babies 1.28 x 0.50 2,000/roll	Labels	5	77000.00
CLS-AGM-200	Agarmate, 200g	Gel Related	5	81255.00	CLS-CRTH-1000	TT Cryo-Babies 1.05 x 0.50" 2,000/roll	Labels	5	56980.00
CLS-ALUM-1000	Alumina Seal, 100/pk	Sealing Film	5	62085.00	CLS-CRTH-2100	TT Cryo-Babies 1.375 x 0.625" 2,000/roll	Labels	5	65450.00
CLS-ARS-150	Autoradiography Ruler Set	Gel Related	5	26400.00	CLS-CRTH-2300	TT Cryo-Babies 1.05 x 0.50" 4,000/roll	Labels	5	1659200.00
CLS-BT-250	Bake Tags 1.50 x 0.75" 250/roll	Labels	5	60280.00	CLS-CRTH-3000	TT Cryo-Babies 1.50 x 0.50" 2,000/roll	Labels	5	84045.00
CLS-BE-1	Breathe-Easy, 100/pk	Sealing Film	5	72945.00	CLS-CRTH-3081	TT Cryo-Babies 0.81" x 0.28" 4,000/roll	Labels	5	174720.00
CLS-BE-1000	Breathe-Easy Strips	PetriDish Accessories	5	54780.00	CLS-CRTH-3094	TT Cryo-Babies 0.94" x 0.50" 4,000/roll	Labels	5	174720.00
CLS-BE-2000	Breathe-EASIER, 100/pk	Sealing Film	5	54780.00	CLS-CRTH-3100	TT Cryo-Babies 1.28" x 0.50" 4,000/roll	Labels	5	165760.00
CLS-TS-1001	7/16" Bullseye Tough-Spots, Dark Blue	Labels	5	28413.00	CLS-CRTH-3128	TT Cryo-Babies 1.28" x 0.50" 4,000/roll	Labels	5	212800.00
CLS-TS-1002	7/16" Bullseye Tough-Spots, Dark Green	Labels	5	28413.00	CLS-CRTH-3150	TT Cryo-Babies 1.28" x 0.50" 4,000/roll	Labels	5	308000.00
CLS-TS-1003	7/16" Bullseye Tough-Spots, Copper	Labels	5	28413.00	CLS-CRTH-3155	TT Cryo-Babies 1.28" x 0.50" 4,000/roll	Labels	5	233137.00
CLS-TS-1004	7/16" Bullseye Tough-Spots, Plum	Labels	5	28413.00	CLS-CRTH-3225	TT Cryo-Tags 2.00 x 0.25" 4,000/roll	Labels	5	179200.00
CLS-TS-1000	3/8" BioHazard Tough-Spots 1,000/roll	Labels	5	19250.00	CLS-CRTH-4000	TT Cryo-Tags 1.50 x 0.75" 2,000/roll	Labels	5	103565.00
CLS-BM-1000	Breathe-EASY Membranes, 100/pk		5	32865.00	CLS-CRTH-5000	TT Cryo-Tags 2.00 x 0.25" 2,000/roll	Labels	5	59675.00
CLS-BM-2000	Breathe-EASIER Membranes, 100/pk		5	32865.00	CLS-CRY-0500	Cryo-Babies 0.94 x 0.50" 1,000/roll WHITE	Labels	5	23100.00
CLS-BOAT-1000	Weigh-Better Boats 150/PK	Weighing Boats	5	24194.50	CLS-CRY-1000	Cryo-Babies 1.28 x 0.50" 1,000/roll WHITE	Labels	5	23100.00
CLS-BOAT-2000	Weigh-Better Boats 24/case	Weighing Boats	5	710480.00	CLS-CRY-1000-B	Cryo-Babies 1.28 x 0.50" 1,000/roll BLUE	Labels	5	26895.00
CLS-BOC-1000	Benchtop Bottle Cooler, 0.5 liter	Benchtop Coolers	5	8580.00	CLS-CRY-1000-G	Cryo-Babies 1.28 x 0.50" 1,000/roll GREEN	Labels	5	26895.00
CLS-BOC-2000	Benchtop Bottle Cooler, 1.0 liter	Benchtop Coolers	5	21175.00	CLS-CRY-1000-O	Cryo-Babies 1.28 x 0.50" 1,000/roll ORANGE	Labels	5	26895.00
CLS-WP-1000	Bio-Pure Alcohol Wipes		5	10043.00	CLS-CRY-1000-R	Cryo-Babies 1.28 x 0.50" 1,000/roll RED	Labels	5	26895.00
CLS-WP-2000	Bio-Pure Surfactant Alcohol Wipes		5	10043.00	CLS-CRY-1000-Y	Cryo-Babies 1.28 x 0.50" 1,000/roll YELLOW	Labels	5	26895.00
CLS-SRB-4000	4" Soft Rubber Brayer	Sealing Film	5	17908.00	CLS-CRY-2000	Cryo-Tags 1.50 x 0.75" 1,000/roll WHITE	Labels	5	23100.00
CLS-CB-1501	1°C Cool Brick	Cooler Accessories	5	26400.00	CLS-CRY-2000-B	Cryo-Tags 1.50 x 0.75" 1,000/roll BLUE	Labels	5	28240.00
CLS-CB-1520	-10°C Cool Brick	Cooler Accessories	5	30800.00	CLS-CRY-2000-G	Cryo-Tags 1.50 x 0.75" 1,000/roll GREEN	Labels	5	28240.00
CLS-CB-2501	1°C Cool Brick	Cooler Accessories	5	26400.00	CLS-CRY-2000-O	Cryo-Tags 1.50 x 0.75" 1,000/roll ORANGE	Labels	5	28240.00
CLS-BAG-1000	Cool Safe Bag for easy transport	Cooler Accessories	5	17600.00	CLS-CRY-2000-R	Cryo-Tags 1.50 x 0.75" 1,000/roll RED	Labels	5	28240.00
CLS-NC-A960	Novel Carriers for 2ml Tubes, Set of 2	Cooling Chambers	1	332800.00	CLS-CRY-2000-Y	Cryo-Tags 1.50 x 0.75" 1,000/roll YELLO	Labels	5	28240.00
CLS-CHAM-A100	96-Well Cooling Chamber for Automation	Cooling Chambers	1	57600.00	CLS-CSAF-1000	CoolSafe System for 0.2ml tubes	Benchtop Coolers	1	79360.00
CLS-CHAM-1000	Chamber for 0.2ml PCR plates or tubes	Cooling Chambers	1	53504.00	CLS-CSAF-2000	CoolSafe System for 0.5ml tubes	Benchtop Coolers	1	79360.00
CLS-CHAM-1020	Double Chamber for 0.2ml tubes	Cooling Chambers	1	94464.00	CLS-CSAF-3000	CoolSafe System for 1.5ml tubes	Benchtop Coolers	1	67072.00
CLS-CHAM-2000	Chamber for 0.5ml tubes	Cooling Chambers	1	49920.00	CLS-CSAF-3050	CoolSafe System for 1.5 Cryogenic tubes	Benchtop Coolers	1	67072.00
CLS-CHAM-3000	Chamber for 1.5ml tubes	Cooling Chambers	1	43520.00	CLS-CSAF-3100	CryoVial CoolSafe System for 1.2 Cryogenic tubes	Benchtop Coolers	1	102400.00
CLS-CHAM-3006	Mini Chamber 1.5-2.0ml Tubes, 6 Holes	Cooling Chambers	5	101200.00	CLS-CSAF-3200	CryoVial CoolSafe System for 2.0 Cryogenic tubes	Benchtop Coolers	1	102400.00
CLS-CHAM-3015	Mini Chamber 1.5-2.0ml Tubes, 15 Holes	Cooling Chambers	5	172800.00	CLS-CSAF-3840	CoolSafe System for 384-well plates	Benchtop Coolers	1	85760.00
CLS-CHAM-3030	Mini Chamber 1.5 -2.0ml Tubes, 30 Holes	Cooling Chambers	5	192000.00	CLS-CSAF-4000	CoolSafe System for 5ml tubes	Benchtop Coolers	1	69880.00
CLS-CHAM-3050	Chamber for 1.5ml Cryogenic tubes	Cooling Chambers	1	45824.00	CLS-CSAF-5000	CoolSafe System for multi-tube format	Benchtop Coolers	1	67072.00
CLS-CHAM-3100	CryoVial Cooling Chamber for 1.2ml Cryogenic tubes	Cooling Chambers	1	83200.00	CLS-CSAF-9600	CoolSafe System for 96-well Microplates	Benchtop Coolers	1	59260.00
CLS-CHAM-3200	CryoVial Cooling Chamber for 2.0ml Cryogenic tubes	Cooling Chambers	1	83200.00	CLS-CSEAL-75R	ContainerSEAL, red	Tape	5	35545.00
CLS-CHAM-3300	CryoVial Cooling Chamber for 3.0-4.0ml Cryogenic tubes	Cooling Chambers	1	96000.00	CLS-CSEAL-75W	ContainerSEAL, white	Tape	5	35545.00
CLS-CHAM-3840	Chamber for 384-well microplates	Cooling Chambers	1	74240.00	CLS-CSEAL-75Y	ContainerSEAL, yellow	Tape	5	35545.00
CLS-CHAM-4000	Chamber for 5ml tubes	Cooling Chambers	1	48640.00	CLS-CSEAL-XLR	ContainerSEAL XL Red 1" x 108'	Tape	5	34650.00
CLS-CHAM-5000	Chamber for multi-tube format	Cooling Chambers	1	48640.00	CLS-CSEAL-XLW	ContainerSEAL XL White 1" x 108'	Tape	5	34650.00
CLS-CHAM-6000	Chamber for 15ml tubes, 6 holes	Cooling Chambers	1	48640.00	CLS-CSEAL-XLY	ContainerSEAL XL Yellow 1" x 108'	Tape	5	34650.00
CLS-CHAM-6050	Chamber for 15ml tubes, 12 holes	Cooling Chambers	1	48640.00	CLS-CSF-BOX	CoolSafe Box	Cooler Accessories	5	25960.00
CLS-CHAM-6400	64-Place Step Chamber	Cooling Chambers	1	80640.00	CLS-DS1-150	Test Tube Seal 1" x 150'	Sealing Film	5	18920.00
CLS-CHAM-7000	Chamber for 50ml tubes, 2 holes	Cooling Chambers	1	40448.00	CLS-DS2-500	Duraseal, 2" x 500'	Sealing Film	5	29920.00
CLS-CHAM-7050	Chamber for 50ml tubes, 4 holes	Cooling Chambers	1	45937.80	CLS-DS4-250	Duraseal, 4" x 250'	Sealing Film	5	42240.00
CLS-CHAM-8000	80-well Chamber for 1.5ml tubes	Cooling Chambers	1	13200.00	CLS-DS4-500	Duraseal, 4" x 500'	Sealing Film	5	55000.00
CLS-CHAM-9600	Microplate Chamber	Cooling Chambers	1	12760.00	CLS-DS5-150	Duraseal, 5" x 150'	Sealing Film	5	36520.00
CLS-CHAM-LIFT	6 1/8"L x 3 1/8" W x 1/2" H Dry Bath Lifting Block	Cooling Chambers	1	13200.00	CLS-DSC1-150	Test Tube Seal 1" x 150'	Sealing Film	5	23760.00
CLS-CHDN-1000	Chill Down Agarose Cooling chamber	Gel Related	1	12760.00	CLS-DSC2-500	Duraseal, 2" x 500'	Sealing Film	5	42240.00
CLS-CLEAR-170	Cryo-Clear Laser Labels 1.28 x 0.50" 1,700/pk,	Labels	5	48510.00	CLS-DSC4-500	Duraseal, 4" x 500'	Sealing Film	5	73040.00
CLS-CLEAR-238	Cryo-Clear Laser Labels 0.94 x 0.50" 2,380/pk	Labels	5	48510.00	CLS-DSC5-150	Duraseal, 5" x 150'	Sealing Film	5	66000.00
CLS-CLTT-1000	Clear Laser Tough-Tags 1.50 x 0.75" 1,200/pk	Labels	5	50050.00	CLS-DSH-100	DuraSharp Shears	Sealing Film	5	9680.00
CLS-CLTT-2000	Clear Laser Tough-Tags 1.50 x 0.50" 1,600/pk	Labels	5	50050.00	CLS-DTCR-1000	DT Cryo-Tags 1.05 x 0.50" 1,000/roll WHITE	Labels	5	33110.00
					CLS-DTCR-1000-BDT	Cryo-Tags 1.05 x 0.50" 1,000/roll BLUE	Labels	5	33110.00
					CLS-DTCR-1000-GDT	Cryo-Tags 1.05 x 0.50" 1,000/roll GREEN	Labels	5	33110.00
					CLS-DTCR-1000-ODT	Cryo-Tags 1.05 x 0.50" 1,000/roll ORANGE	Labels	5	33110.00



Code	Product Description	Label	Pk Size	Price	Code	Product Description	Label	Pk Size	Price
CLSDTCR-2000	DT Cryo-Tags 1.50 x 0.50" 750/roll WHITE	Labels	5	30800.00	CLSLBTP-2200	Labeling Tape, Green 3/4" x 2,160"/ Roll	Tape	5	24640.00
CLSDTCR-200-B	DT Cryo-Tags 1.50 x 0.50" 750/roll BLUE	Labels	5	32340.00	CLSLBTP-2300	Labeling Tape, Blue 3/4" x 2,160" / Roll	Tape	5	24640.00
CLSDTCR-2000-G	DT Cryo-Tags 1.50 x 0.50" 750/roll GREEN	Labels	5	32340.00	CLSLBTP-2400	Labeling Tape, Red 3/4" x 2,160" / Roll	Tape	5	24640.00
CLSDTCR-2000-Y	DT Cryo-Tags 1.50 x 0.50" 750/roll YELLOW	Labels	5	32340.00	CLSLCRY-1000	Laser Cryo-Tag Sheet, 20/pk	Labels	5	38500.00
CLSDTCR-3000	DT Cryo-Tags 1.50 x 0.75" 750/roll WHITE	Labels	5	33110.00	CLSLCRY-1100	Laser Cryo-Tags 1.69 x 0.75" 1,040/pk	Labels	5	46200.00
CLSDTCR-3000-B	DT Cryo-Tags 1.50 x 0.75" 750/roll BLUE	Labels	5	34820.00	CLSLCRY-1100-B	Laser Cryo-Tags 1.69 x 0.75" Blue	Labels	5	46200.00
CLSDTCR-3000-G	DT Cryo-Tags 1.50 x 0.75" 750/roll GREEN	Labels	5	34820.00	CLSLCRY-1100-G	Laser Cryo-Tags 1.69 x 0.75" Green	Labels	5	46200.00
CLSDTCR-3000-Y	DT Cryo-Tags 1.50 x 0.75" 750/roll YELLOW	Labels	5	34820.00	CLSLCRY-1100-O	Laser Cryo-Tags 1.69 x 0.75" Orange	Labels	5	46200.00
CLSDTCR-4000	DT Cryo-Spots. 3/8" diam., 1,000/roll WHITE	Labels	5	50165.00	CLSLCRY-1100-R	Laser Cryo-Tags 1.69 x 0.75" Red	Labels	5	46200.00
CLSDTCR-4001	DT Cryo-Spots. 3/8" diam., 1,000/roll BLUE	Labels	5	50165.00	CLSLCRY-1100-Y	Laser Cryo-Tags 1.69 x 0.75" Yellow	Labels	5	46200.00
CLSDTCR-4002	DT Cryo-Spots. 3/8" diam., 1,000/roll GREEN	Labels	5	50165.00	CLSLCRY-1200	Laser Cryo-Tags 1.50 x 0.75" 1,200/pk	Labels	5	46200.00
CLSDTCR-4003	DT Cryo-Spots. 3/8" diam., 1,000/roll YELLOW	Labels	5	50165.00	CLSLCRY-1200-B	Laser Cryo-Tags 1.50 x 0.75" Blue	Labels	5	46200.00
CLSDTCR-5000	DT Cryo-Tags 2.00" x 0.25" 1,000/roll WHITE	Labels	5	64083.50	CLSLCRY-1200-G	Laser Cryo-Tags 1.50 x 0.75" Green	Labels	5	46200.00
CLSDTCR-6000	DT Cryo-Tag 1.50" x .50" and 3/8" Cryo Spot Combo 500 sets/roll	Labels	5	44660.00	CLSLCRY-1200-O	Laser Cryo-Tags 1.50 x 0.75" Orange	Labels	5	46200.00
CLSDTCR-6000-B	DT Cryo-Tag 1.50" x .50" and 3/8" Cryo Spot Combo 500 sets/roll	Labels	5	44660.00	CLSLCRY-1200-R	Laser Cryo-Tags 1.50 x 0.75" Red	Labels	5	46200.00
CLSDTCR-6000-G	DT Cryo-Tag 1.50" x .50" and 3/8" Cryo Spot Combo 500 sets/roll	Labels	5	44660.00	CLSLCRY-1200-Y	Laser Cryo-Tags 1.50 x 0.75" Yellow	Labels	5	46200.00
CLSDTCR-6000-Y	DT Cryo-Tag 1.50" x .50" and 3/8" Cryo Spot Combo 500 sets/roll	Labels	5	44660.00	CLSLCRY-1258	Laser Cryo-Tags 2.625 x 1.00" 600/pk	Labels	5	46200.00
CLSDYMO-0450	Dymo Label Writer 450	Printer	5	324000.00	CLSLCRY-1258-B	Laser Cryo-Tags 2.625 x 1.00" 600/pk Blue	Labels	5	46200.00
CLSLFLAG-1000	Laser Flag Tough-Tags 2.50 x 0.50" 1,500/pk	Labels	5	52580.00	CLSLCRY-1258-G	Laser Cryo-Tags 2.625 x 1.00" 600/pk Green	Labels	5	46200.00
CLSGM-100	Flow Mesh, 20 sheets	Gel Related	5	39600.00	CLSLCRY-1258-O	Laser Cryo-Tags 2.625 x 1.00" 600/pk Orange	Labels	5	46200.00
CLSGFRK-100	Fluorescent Ruler Kit	Gel Related	5	23285.00	CLSLCRY-1258-R	Laser Cryo-Tags 2.625 x 1.00" 600/pk Red	Labels	5	46200.00
CLSGH-100	Gel Handler, 5 mil thick - 10 shts	Gel Related	5	31315.00	CLSLCRY-1258-Y	Laser Cryo-Tags 2.625 x 1.00" 600/pk Yellow	Labels	5	46200.00
CLSGH-200	Gel Handler, 10 mil thick - 5 shts	Gel Related	5	31315.00	CLSLCRY-1300	Laser Cryo-Tags 1 1/8 x 1 3/4" 840/pk	Labels	5	46200.00
CLSGLWR-1000	Glow-Writer Pen	Gel Related	5	48180.00	CLSLCRY-1400	Laser Cryo-Babies 1.0" x 1.0" 1,600/pk	Labels	5	46200.00
CLSGNIB-1000	Glow Writer Nibs, 4/pack	Gel Related	5	8030.00	CLSLCRY-1700	Laser Cryo-Babies 1.28 x 0.50" 1,700/pk	Labels	5	46200.00
CLSGOO-100	Goo-Off, lemon scented - 1 oz.	Laboratory Cleaning Solutions	5	6390.00	CLSLCRY-1700-B	Laser Cryo-Babies 1.28 x 0.50" Blue	Labels	5	46200.00
CLSGOO-200	Goo-Off, lime scented - 1 oz.	Laboratory Cleaning Solutions	5	6390.00	CLSLCRY-1700-G	Laser Cryo-Babies 1.28 x 0.50" Green	Labels	5	46200.00
CLSGOO-300	Goo-Off, cherry scented - 1 oz.	Laboratory Cleaning Solutions	5	6390.00	CLSLCRY-1700-O	Laser Cryo-Babies 1.28 x 0.50" Orange	Labels	5	46200.00
CLSGRDT-1000	DT WetGrip 1.05 x 0.50" 1,000/roll	Labels	5	36190.00	CLSLCRY-1700-R	Laser Cryo-Babies 1.28 x 0.50" Red	Labels	5	46200.00
CLSGRDT-2000	DT WetGrip 1.50 x 0.50" 750/roll	Labels	5	38500.00	CLSLCRY-1700-Y	Laser Cryo-Babies 1.28 x 0.50" Yellow	Labels	5	46200.00
CLSGRDT-3000	DT WetGrip 1.50 x 0.75" 750/roll	Labels	5	34650.00	CLSLCRY-2380	Laser Cryo-Babies 0.94 x 0.50" 2,380/pk	Labels	5	46200.00
CLSGRDT-4000	DT WetGrip 3/8" diam., 1,000/roll	Labels	5	52475.00	CLSLCRY-2380-B	Laser Cryo-Babies 0.94 x 0.50" Blue	Labels	5	46200.00
CLSGRDT-5000	DT WetGrip 2.00" x 0.25" 1,000/roll	Labels	5	60934.50	CLSLCRY-2380-G	Laser Cryo-Babies 0.94 x 0.50" Green	Labels	5	46200.00
CLSGRDT-6000	DT WetGrip 1.50" x .50" and 3/8" Tough Spot Combo 500 sets/roll	Labels	5	46200.00	CLSLCRY-2380-O	Laser Cryo-Babies 0.94 x 0.50" Orange	Labels	5	46200.00
CLSGRID-1000	Slide Grids, 20mm x 20mm	Microscopy	5	25564.00	CLSLCRY-2380-R	Laser Cryo-Babies 0.94 x 0.50" Red	Labels	5	46200.00
CLSGRLA-1000	Microscope Slide WetGrip 7/8" x 7/8" 2,400/pk	Labels	5	52475.00	CLSLCRY-2380-Y	Laser Cryo-Babies 0.94 x 0.50" Yellow	Labels	5	46200.00
CLSGRLA-1050	Laser WetGrip 1.50 x 0.75" 1,500/pk	Labels	5	52475.00	CLSLCRY-2500	Laser Cryo-Tags 2.50" x 0.50" 1,200/pk	Labels	5	46200.00
CLSGRLA-2016	Laser WetGrip 1.28 x 0.50" 1,700/pk	Labels	5	52475.00	CLSLCRY-3020	Laser Cryo-Tags 3.0 x 2.0" 240/pk	Labels	5	46200.00
CLSGRLA-2380	Laser WetGrip 0.94 x 0.50" 2,380/pk	Labels	5	52475.00	CLSLTBK-1000	Light-Bloc Labels 4.75" x 2.25" 100/roll	Labels	5	30800.00
CLSGRRL-1000	WetGrip 1.28 x 0.50" 1,000/roll	Labels	5	26950.00	CLSLTBK-2000	Light-Bloc Labels 4.75" x 3.75" 100/roll	Labels	5	34650.00
CLSGRRL-2000	WetGrip 1.50 x 0.75" 1,000/roll	Labels	5	26950.00	CLSLTTI-1000	Laser Tough-Tags 1.00 x 1.00" 1,600/pk	Labels	5	50160.00
CLSGUIDE-300	Pipetter Guide, single place	Microwell Plate Accessories	5	202860.00	CLSLTTI-2000	Laser Tough-Tags 1 5/8 x 5/8" 1,600/pk	Labels	5	50160.00
CLSGUIDE-400	Pipetter Guide, double place	Microwell Plate Accessories	5	203933.50	CLSLTTM-1000	Laser Teeny Tough-Tags 0.91 x 0.32" 3,850/pk	Labels	5	50160.00
CLSHALF-1000	Clear over White Laser Tough-Tags 520/pk	Labels	5	52550.00	CLSMINI-1700	Mini Laser Cryo-Babies 1.28" x 0.50", 1,020/pk	Labels	5	36190.00
CLSHOTS-1000	HotSeal, 100/pk	Sealing Film	5	54780.00	CLSMISL-1000	Microscope Slide Tough-Tags 7/8" x 7/8" 2,400/pk	Labels	5	50160.00
CLSIBUC-1000	2.5 Liter Triple Density Ice Bucket Icy Blue	Ice Pan Ice Bucket	5	26950.00	CLSMISL-1000-B	Microscope Slide Tough-Tags, Blue	Labels	5	50160.00
CLSIBUC-1100	2.5 Liter Triple Density Ice Bucket Icy Green	Ice Pan Ice Bucket	5	26950.00	CLSMISL-1000-G	Microscope Slide Tough-Tags, Green	Labels	5	50160.00
CLSIBUC-2000	5 Liter Triple Density Ice Bucket Icy Blue	Ice Pan Ice Bucket	5	34650.00	CLSMISL-1000-O	Microscope Slide Tough-Tags, Orange	Labels	5	50160.00
CLSIBUC-2100	5 Lite Triple Density Ice Bucket, Icy Green	Ice Pan Ice Bucket	5	34650.00	CLSMISL-1000-R	Microscope Slide Tough-Tags, Red	Labels	5	50160.00
CLSID-100	Identi-Kit, 174 symbols	Gel Related	5	31955.00	CLSMISL-1000-Y	Microscope Slide Tough-Tags, Yellow	Labels	5	50160.00
CLSIPAN-1000	1 Liter Triple Density Ice Pan Icy Blue	Ice Pan Ice Bucket	5	11550.00	CLSMIWE-1000	Slide Wells, 0.25mm thick	Microscopy	5	38720.00
CLSIPAN-1100	1 Liter Triple Density Ice Pan Icy Green	Ice Pan Ice Bucket	5	11550.00	CLSMIWE-2000	Slide Wells, 0.50mm thick	Microscopy	5	47080.00
CLSIPAN-2000	4 Liter Triple Density Ice Pan Icy Blue	Ice Pan Ice Bucket	5	26950.00	CLSNEW-100	New Ice, 45/pack	Reusable Ice Packs	5	8580.00
CLSIPAN-2100	4 Liter Triple Density Ice Pan Icy Green	Ice Pan Ice Bucket	5	26950.00	CLSNOC-1000	NovaChamber, red	Rack	5	56100.00
CLSIPAN-3000	9 Liter Triple Density Ice Pan Icy Blue	Ice Pan Ice Bucket	5	37730.00	CLSNOC-2000	NovaChamber, teal	Rack	5	56100.00
CLSIPAN-3100	9 Liter Triple Density Ice Pan Icy Green	Ice Pan Ice Bucket	5	37730.00	CLSNOC-3000	NovaChamber, yellow	Rack	5	56100.00
CLSLBTP-1000	Labeling Tape, White 3/4" x 500"/Roll	Tape	5	8085.00	CLSNOC-4000	NovaChamber, pink	Rack	5	56100.00
CLSLBTP-1100	Labeling Tape, Yellow 3/4" x 500"/Roll	Tape	5	8085.00	CLSNOC-5000	NovaChamber, purple	Rack	5	56100.00
CLSLBTP-1200	Labeling Tape, Green 3/4" x 500"/ Roll	Tape	5	8085.00	CLSNORC-1000	NovaRack, red	Rack	5	39600.00
CLSLBTP-1300	Labeling Tape, Blue 3/4" x 500"/ Roll	Tape	5	8085.00	CLSNORC-2000	NovaRack, teal	Rack	5	39600.00
CLSLBTP-1400	Labeling Tape, Red 3/4" x 500" / Roll	Tape	5	8085.00	CLSNORC-3000	NovaRack, yellow	Rack	5	39600.00
CLSLBTP-2000	Labeling Tape, White 3/4" x 2,160"/Roll	Tape	5	24640.00	CLSNORC-4000	NovaRack, pink	Rack	5	39600.00
CLSLBTP-2100	Labeling Tape, Yellow 3/4" x 2,160"/Roll	Tape	5	24640.00	CLSNORC-5000	NovaRack, purple	Rack	5	39600.00
					CLSNOST-1000	NovaRack stacking adapter	Rack	5	8580.00
					CLSPDHL-1000	Petri Dish Support N Seal	Petri Dish Accessories	5	43305.00
					CLSPDOR-1032	PetriOrienter, 32-square grid	--	5	12045.00
					CLSPDOR-1050	PetriOrienter, 50-square grid	--	5	12045.00



Code	Product Description	Lable	Pk Size	Price	Code	Product Description	Lable	Pk Size	Price
CLSPDOR-1070	PetriOrienter, 70-square grid	--	5	12045.00	CLSRNBW-2400	Rainbow Pack for SIDE-1000	Labels	5	46200.00
CLSPDOR-1100	PetriOrienter, 100-square grid	--	5	12045.00	CLSRNBW-2500	Rainbow Pack for LCRY-1258	Labels	5	46200.00
CLSPDOR-2006	PetriOrienter, 6-sector pie	--	5	12045.00	CLSRNBW-3000	Rainbow 1.28 x 0.50" Laser Tough-Tags	Labels	5	50160.00
CLSPDOR-2008	PetriOrienter, 8-sector pie	--	5	12045.00	CLSRNBW-3100	Rainbow 0.94 x 0.50" Laser Tough-Tags	Labels	5	50160.00
CLSPDOR-2012	PetriOrienter, 12-sector pie	--	5	12045.00	CLSRNBW-3200	Rainbow 1.50 x 0.75" Laser Tough-Tags	Labels	5	50160.00
CLSPDOR-2016	PetriOrienter, 16-sector pie	--	5	12045.00	CLSRNBW-3300	Rainbow Microscope Slide Tough-Tags	Labels	5	50160.00
CLSPDST-1000	PetriDish Stand 88mm diameter	--	5	45650.00	CLSRock-1000	Rocker Tape on a roll, 1.50" x 20'	Tapes	5	26400.00
CLSPDST-2000	PetriDish Stand 100mm diameter	--	5	45650.00	CLSRock-2000	Rocker Tape Sheets, 12" x 12", 3/pk	Tapes	5	31680.00
CLSPeel-1000	Breathable Peel-n-Seal, 1.50" dia. 250ml Flask 50/pack	Sealing Film	5	27390.00	CLSRsPO-1000	Printed Tough-Spots, 1 to 10	Labels	5	19250.00
CLSPeel-1100	Breathable Peel-n-Seal, 1.63" dia. 500ml Flask 50/pack	Sealing Film	5	27390.00	CLSRsPO-2000	Printed Tough-Spots, 1 to 100	Labels	5	19250.00
CLSPeel-1101	Breathable Peel-n-Seal, 2.13" dia. 500ml Flask 50/pack	Sealing Film	5	27390.00	CLSRsPO-3000	Printed Tough-Spots for microplates	Labels	5	19250.00
CLSPeel-2000	Non-breathable Peel-n-Seal, 1.50" dia. 50/pack	Sealing Film	5	27390.00	CLSRWTS-1000	3/8" Radioactive Warning Tough-Spots 1,000/roll	Labels	5	19250.00
CLSPeel-2100	Non-breathable Peel-n-Seal, 1.63" dia. 50/pack	Sealing Film	5	27390.00	CLSSCS-1	3/4" Binder Clips, 10/pk	Gel Related	5	21780.00
CLSPeel-2200	Non-breathable Peel-n-Seal, 2.13" dia. 50/pack	Sealing Film	5	27390.00	CLSSCS-3	2" Binder Clips, 8/pk	Gel Related	5	32670.00
CLSPICK-1000	Toothpick Dispenser w/polypropylene case	Cell Culture Control System Accessories	5	45650.00	CLSSett-1000	SlenderEdge Laser Tough-Tags, 2,475/pk	Labels	5	50160.00
CLSPICK-2000	Kokeshi Style Toothpicks, 350/pk.	Cell Culture Control System Accessories	5	4950.00	CLSSFBK-0250	250ml Sure-Footed Thermo-Beaker, 50/cs	--	5	63910.00
CLSPSeal-108B	PetriSEAL, blue	PetriDish Accessories	5	27390.00	CLSSIDE-1000	Laser Sidewall Cryo-Tags 1.50 x 0.25" 3,120/pk	Labels	5	46200.00
CLSPSeal-108C	PetriSEAL, clear	PetriDish Accessories	5	27390.00	CLSSIDE-1100	Laser Sidewall Cryo-Tags 1.50 x 0.25" Blue	Labels	5	46200.00
CLSPSeal-108R	PetriSEAL, red	PetriDish Accessories	5	27390.00	CLSSIDE-1200	Laser Sidewall Cryo-Tags 1.50 x 0.25" Green	Labels	5	46200.00
CLSPSeal-108W	PetriSEAL, white	PetriDish Accessories	5	27390.00	CLSSIDE-1300	Laser Sidewall Cryo-Tags 1.50 x 0.25" Orange	Labels	5	46200.00
CLSPSeal-108Y	PetriSEAL, yellow	PetriDish Accessories	5	27390.00	CLSSIDE-1400	Laser Sidewall Cryo-Tags 1.50 x 0.25" Red	Labels	5	46200.00
CLSPSTK-1000	PetriStickers, 32-square grid 36/pk	PetriDish Accessories	5	15975.00	CLSSIDE-1500	Laser Sidewall Cryo-Tags 1.50 x 0.25" Yellow	Labels	5	46200.00
CLSPSTK-1050	PetriStickers, 50-square grid 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-1000	Lab-Sorb 16 oz jar	Waste Disposal	5	28105.00
CLSPSTK-1070	PetriStickers, 70-square grid 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-2000	Kit (Includes Lab-Sorb 16 oz jar plus small bags (25)	Waste Disposal	5	40150.00
CLSPSTK-1100	PetriStickers, 100-square grid 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-2100	Kit (Includes Lab-Sorb 16 oz jar plus medium bags (25)	Waste Disposal	5	40150.00
CLSPSTK-2000	PetriStickers, 6-sector pie 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-2200	Kit Includes Lab-Sorb 16 oz jar plus Large bags	Waste Disposal	5	40150.00
CLSPSTK-2010	PetriStickers, 8-sector pie 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-2300	Variety kit Lab-sorb plus variety size bags (15sm, 10md, 5lg)	Waste Disposal	5	40150.00
CLSPSTK-2020	PetriStickers, 12-sector pie 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-3000	50 small bags	Waste Disposal	5	40150.00
CLSPSTK-2030	PetriStickers, 16-sector pie 36/pk	PetriDish Accessories	5	15975.00	CLSSORB-3100	50 medium bags	Waste Disposal	5	24090.00
CLSPSTK-3000	PetriStickers, 3 Color Grid 18/pk	PetriDish Accessories	5	14300.00	CLSSORB-3200	50 large bags	Waste Disposal	5	24090.00
CLSRAD-10	Radtape, 4" x 6" - 10/pk	Gel Related	5	51040.00	CLSSPOT-1000	3/8" Laser Tough-Spots, White 3,840/pk	Labels	5	24090.00
CLSRCS-50	Rapid Coomassie Stain, 100ml	Gel Staining	5	83600.00	CLSSPOT-1100	3/8" Laser Tough Spots, Blue	Labels	5	46200.00
CLSRBP-1000	Bio-Pure 100 ml Biodegradable Reservoir-Sterile	Reagent Reservoir	5	73040.00	CLSSPOT-1200	3/8" Laser Tough-Spots, Green	Labels	5	46200.00
CLSRBP-2000	Bio-Pure 25ml w/divider Biodegradable Reservoir-Sterile	Reagent Reservoir	5	59345.00	CLSSPOT-1300	3/8" Laser Tough-Spots, Orange	Labels	5	46200.00
CLSRBP-3000	Bio-Pure 25ml Biodegradable Reservoir- Sterile	Reagent Reservoir	5	59345.00	CLSSPOT-1400	3/8" Laser Tough-Spots, Red	Labels	5	46200.00
CLSRBP-5000	Bio-Pure 50ml Biodegradable Reservoir - Sterile	Reagent Reservoir	5	73040.00	CLSSPOT-1500	3/8" Laser Tough-Spots, Yellow	Labels	5	46200.00
CLSRSE-1000	100ml Reservoir, 100/case, sterile	Reagent Reservoir	5	70300.00	CLSSPOT-2000	1/2" Laser Tough-Spots 3,840/pk	Labels	5	46200.00
CLSRSE-1050	100ml Reservoir, 50 individually wrapped, sterile	Reagent Reservoir	5	62540.00	CLSSPOT-2100	1/2" Laser Tough-Spots, Blue	Labels	5	46200.00
CLSRSE-1100	100ml Reservoir, 100/case, Non-sterile	Reagent Reservoir	5	65120.00	CLSSPOT-2200	1/2" Laser Tough-Spots, Green	Labels	5	46200.00
CLSRSE-2000	25ml Reservoir with divider, 100/case, sterile	Reagent Reservoir	5	60255.00	CLSSPOT-2300	1/2" Laser Tough-Spots, Orange	Labels	5	46200.00
CLSRSE-2050	25ml Reservoir with divider, 50 individually wrapped, sterile	Reagent Reservoir	5	52955.00	CLSSPOT-2400	1/2" Laser Tough-Spots, Red	Labels	5	46200.00
CLSRSE-2100	25ml Reservoir with divider, 100/case, Non-sterile	Reagent Reservoir	5	51040.00	CLSSPOT-2500	1/2" Laser Tough-Spots, Yellow	Labels	5	46200.00
CLSRSE-3000	25ml Reservoir, 100/case, sterile	Reagent Reservoir	5	60258.00	CLSSPOT-3000	3/4" Laser Tough-Spots 2,160/pk	Labels	5	46200.00
CLSRSE-3050	25ml Reservoir, 50 individually wrapped, sterile	Reagent Reservoir	5	54560.00	CLSSPOT-4000	1.00" Laser Tough-Spots 1,260/pk	Labels	5	46200.00
CLSRSE-3100	25ml Reservoir 100/case, Non-sterile	Reagent Reservoir	5	51040.00	CLSSPOT-5000	7/16" Laser Tough-Spots 4,680/pk	Labels	5	46200.00
CLSRSE-4000	100ml Sidecar Reservoir, 100/case, sterile	Reagent Reservoir	5	74865.00	CLSSPRL-1000	3/8" Tough-Spots , 5-up 5,000/roll	Labels	5	76615.00
CLSRSE-4050	100ml Sidecar Reservoir, 50 individually wrapped, sterile	Reagent Reservoir	5	67560.00	CLSSPRL-1100	3/8" Tough-Spots, 5-up Assorted Colors 5,000/roll	Labels	5	76615.00
CLSRSE-4100	100ml Sidecar Reservoir, 100/case, Non-sterile	Reagent Reservoir	5	65120.00	CLSSPRL-2000	1/2" Tough-Spots, 5-up 5,000/roll	Labels	5	76615.00
CLSRSE-5000	50ml Reservoir, 100/case, sterile	Reagent Reservoir	5	70300.00	CLSSPRL-2100	1/2" Tough-Spots, 5-up Assorted Colors 5,000/roll	Labels	5	76615.00
CLSRSE-5050	50ml Reservoir, 50 Individually wrapped, sterile	Reagent Reservoir	5	62084.00	CLSSTK-57R	Sticky Stand, regular 5" x 7"	Test tube accessories	5	31315.00
CLSRSE-5100	50ml Reservoir, 100/case, Non-sterile	Reagent Reservoir	5	60255.00	CLSSTK-57SS	Sticky Stand, super 5" x 7"	Test tube accessories	5	31315.00
CLSRSE-C025	Cooler for 25ml Reservoir	Reagent Reservoir	5	27390.00	CLSSUE-100	Sweet Sue, balsam fir, 10/pk	Odor Control Products	5	12848.00
CLSRSE-C050	Cooler for 50ml Reservoir	Reagent Reservoir	5	27390.00	CLSSUE-200	Sweet Sue, wild cherry, 10/pk	Odor Control Products	5	12848.00
CLSRSE-C100	Cooler for 100ml Reservoir	Reagent Reservoir	5	27390.00	CLSSUE-300	Sweet Sue, sage & citrus, 10/pk	Odor Control Products	5	12848.00
CLSRest-1000	Angled Bottle Rest	Microwell Plate Accessories	5	59935.00	CLSSUE-400	Sweet Sue for animal facilities, 10/pk	Odor Control Products	5	12848.00
CLSRNBW-1000	Rainbow Pack for SPOT-1000	Labels	5	46200.00	CLSSUEP-1000	Sweet Sue Pouch, 6/pk	Odor Control Products	5	3575.00
CLSRNBW-1100	Rainbow Pack for SPOT-2000	Labels	5	46200.00	CLSSYN-100	Synergel, 100g	Gel Related	5	87120.00
CLSRNBW-2000	Rainbow Pack for LCRY-1100	Labels	5	46200.00	CLSTH-20	Thirsty Stix, 6.0" x 0.20"	Absorbent Materials	5	22825.00
CLSRNBW-2100	Rainbow Pack for LCRY-1200	Labels	5	46200.00	CLSTHER-1094	TT Tough-Tags 0.94 x 0.50" 2,000/roll	Labels	5	56980.00
CLSRNBW-2200	Rainbow Pack for LCRY-1700	Labels	5	46200.00	CLSTHER-1128	TT Tough-Tags 1.28 x 0.50" 2,000/roll	Labels	5	63910.00
CLSRNBW-2300	Rainbow Pack for LCRY-2380	--	5	46200.00	CLSTHER-1137	TT Tough-Tags 1.375 x 0.625" 2,000/roll	Labels	5	63910.00
					CLSTHER-1150	TT Tough-Tags 1.50 x 0.75" 2,000/roll	Labels	5	165005.00
					CLSTHER-1225	TT Tough-Tags 2.0 x 0.25" 2,000/roll	Labels	5	63910.00
					CLSTHER-3094	TT Tough-Tags 0.94 x 0.50" 4,000/roll	Labels	5	173935.00



Code	Product Description	Label	Pk Size	Price	Code	Product Description	Label	Pk Size	Price
CLSTHER-3128	TT Tough-Tags 1.28 x 0.50" 4,000/roll	Labels	5	191080.00	CLSTTLW-2016-Y	Laser Tough-Tags 1.28 x 0.50" 2,125/pk, Yellow	Labels	5	50160.00
CLSTHER-3096	TT Tough-Tags 1.375 x 0.625" 4,000/roll	Labels	5	191080.00	CLSTTLW-KK-20	KissCut Tough-Tags 1.28 x 0.50" 2,520/pk	Labels	5	50160.00
CLSTHER-3097	TT Tough-Tags 1.50 x 0.75" 4,000/roll	Labels	5	105460.00	CLSTTLY-1000	Tough-Tags 1.28 x 0.50" Yellow	Labels	5	26755.00
CLSTHER-3225	TT Tough-Tags 2.0 x 0.25" 4,000/roll	Labels	5	192973.00	CLST-TOPS-100	Titer Tops, 100/pk	Sealing Film	5	54780.00
CLSTHPR-4000	Zebra GK420T Thermal Transfer Printer	Printer	1	253440.00	CLST-TOPS-50	Titer Tops, 50/pk	Sealing Film	5	31955.00
CLSTHRI-5000	Resin Ribbon 2.50" x 244' for TLP 2844	Printer Supplies	5	20460.00	CLSTTSA-9000	Tough-Tag Station	Label Accessories	5	51040.00
CLSTHUM-1000	0.5ml Thumbs-Up, Clear 1000/pack	Microtubes	5	52800.00	CLSTTSB-1000	Tough-Tags 0.94 x 0.50" Blue	Labels	5	23100.00
CLSTHUM-1100	0.5ml Thumbs-Up, Blue 1000/pack	Microtubes	5	52800.00	CLSTTSC-1000	Tough-Tags 0.94 x 0.50" Clear	Labels	5	23100.00
CLSTHUM-1200	0.5ml Thumbs-Up, Green 1000/pack	Microtubes	5	52800.00	CLSTTSG-1000	Tough-Tags 0.94 x 0.50" Green	Labels	5	23100.00
CLSTHUM-1300	0.5ml Thumbs-Up, Pink 1000/pack	Microtubes	5	52800.00	CLSTTSH-1000	Laser Tough-Tag sheet, 25/pk	Labels	5	48510.00
CLSTHUM-1400	0.5ml Thumbs-Up, Purple 1000/pack	Microtubes	5	52800.00	CLSTTSO-1000	Tough-Tags 0.94 x 0.50" Orange	Labels	5	23100.00
CLSTHUM-1500	0.5ml Thumbs-Up, Yellow 1000/pack	Microtubes	5	52800.00	CLSTTSR-1000	Tough-Tags 0.94 x 0.50" Red	Labels	5	23100.00
CLSTHUM-2000	1.5ml Thumbs-Up, Clear 500/pack	Microtubes	5	27390.00	CLSTTSW-1000	Tough-Tags 0.94 x 0.50" 1,000/roll	Labels	5	23100.00
CLSTHUM-2100	1.5ml Thumbs-Up, Blue 500/pack	Microtubes	5	27390.00	CLSTTSW-2240	Laser Tough-Tags 0.94 x 0.50" 2,975/pk	Labels	5	50380.00
CLSTHUM-2200	1.5ml Thumbs-Up, Green 500/pack	Microtubes	5	27390.00	CLSTTSW-2240-B	Laser Tough-Tags 0.94 x 0.50" 2,975/pk, Blue	Labels	5	50380.00
CLSTHUM-2300	1.5ml Thumbs-Up, Pink 500/pack	Microtubes	5	27390.00	CLSTTSW-2240-G	Laser Tough-Tags 0.94 x 0.50" 2,975/pk, Green	Labels	5	50380.00
CLSTHUM-2400	1.5ml Thumbs-Up, Purple 500/pack	Microtubes	5	27390.00	CLSTTSW-2240-O	Laser Tough-Tags 0.94 x 0.50" 2,975/pk, Orange	Labels	5	50380.00
CLSTHUM-2500	1.5ml Thumbs-Up, Yellow 500/pack	Microtubes	5	27390.00	CLSTTSW-2240-R	Laser Tough-Tags 0.94 x 0.50" 2,975/pk, Red	Labels	5	50380.00
CLSTNY-TWIN	Twin-Tip Teeny Writer, black	Laboratory Marker	5	4993.00	CLSTTSW-2240-Y	Laser Tough-Tags 0.94 x 0.50" 2,975/pk, Yellow	Labels	5	50380.00
CLSTNY-WTR	Teeny Writer, black	Laboratory Marker	5	4993.00	CLSTT-SWALL	Sidewall Tough-Tags 1.50 x 0.25" 1,000/roll	Labels	5	23100.00
CLST-SPOTS	3/8" Tough-Spots 1,000/roll	Labels	5	19250.00	CLSTTSG-KK-20	KissCut Tough-Tags 0.94 x 0.50" 3,200/pk	Labels	5	46200.00
CLST-SPOTS-100	1.00" Tough-Spots 500/roll	Labels	5	21829.00	CLSTTSY-1000	Tough-Tags 0.94 x 0.50" Yellow	Labels	5	21560.00
CLST-SPOTS-44	7/16" Tough-Spots 1,000/roll	Labels	5	19250.00	CLSTT-TNY	Teeny Tough-Tags 0.81 x 0.28" 1,500/roll	Labels	5	21560.00
CLST-SPOTS-50	1/2" Tough-Spots 1,000/roll	Labels	5	19250.00	CLSVIBE-1000	7" x 8.5", 2/pk	Vibration Dampening	5	21560.00
CLST-SPOTS-50B	1/2" Tough-Spots, Blue	Labels	5	19250.00	CLSVIBE-2000	8.5" x 11" 2/pk	Vibration Dampening	5	31955.00
CLST-SPOTS-50C	1/2" Tough-Spots, Coffee	Labels	5	19250.00	CLSVIBE-3000	Combo Pack, 2/pk	Vibration Dampening	5	31955.00
CLST-SPOTS-50D	1/2" Tough-Spots, Gold	Labels	5	19250.00	CLSWASH-1000	Wash-Off Labels 2.625" x 1.00" 750/pk	Labels	5	91520.00
CLST-SPOTS-50G	1/2" Tough-Spots, Green	Labels	5	19250.00	CLSWASH-2000	Wash-Off Labels 4.00" x 1.00" 500/pk	Labels	5	91520.00
CLST-SPOTS-50L	1/2" Tough-Spots, Lavender	Labels	5	19250.00	CLSWASH-3100	1"x1" /25mm x 25mm 250/roll Warning: Flammable	Label	5	30800.00
CLST-SPOTS-50M	1/2" Tough-Spots, Maroon	Labels	5	19250.00	CLSWASH-3200	1"x1" /25mm x 25mm 250/roll Warning: Poison	Label	5	30800.00
CLST-SPOTS-50O	1/2" Tough-Spots, Orange	Labels	5	19250.00	CLSWASH-3300	1"x1" /25mm x 25mm 250/roll Warning: Biohazard	Label	5	30800.00
CLST-SPOTS-50P	1/2" Tough-Spots, Pink	Labels	5	19250.00	CLSWASH-3400	1"x1" /25mm x 25mm 250/roll Warning: Corrosive	Label	5	30800.00
CLST-SPOTS-50R	1/2" Tough-Spots, Red	Labels	5	19250.00	CLSWASH-3500	1"x1" /25mm x 25mm 250/roll Warning: Radioactive	Label	5	30800.00
CLST-SPOTS-50S	1/2" Tough-Spots, Silver	Labels	5	19250.00	CLSWASH-3600	1"x1" /25mm x 25mm 250/roll Warning: Ready for Use	Label	5	30800.00
CLST-SPOTS-50Y	1/2" Tough-Spots, Yellow	Labels	5	19250.00	CLSWASH-3700	1"x1" /25mm x 25mm 250/roll Warning: Do Not Use	Label	5	30800.00
CLST-SPOTS-75	3/4" Tough-Spots 500/roll	Labels	5	21829.00	CLSWASH-4100	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Flammable	Label	5	38500.00
CLST-SPOTS-B	3/8" Tough-Spots, Blue	Labels	5	19250.00	CLSWASH-4200	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Poison	Label	5	38500.00
CLST-SPOTS-C	3/8" Tough-Spots, Coffee	Labels	5	19250.00	CLSWASH-4300	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Biohazard	Label	5	38500.00
CLST-SPOTS-D	3/8" Tough-Spots, Gold	Labels	5	19250.00	CLSWASH-4400	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Corrosive	Label	5	38500.00
CLST-SPOTS-G	3/8" Tough-Spots, Green	Labels	5	19250.00	CLSWASH-4500	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Radioactive	Label	5	38500.00
CLST-SPOTS-L	3/8" Tough-Spots, Lavender	Labels	5	19250.00	CLSWASH-4600	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Ready for Use	Label	5	38500.00
CLST-SPOTS-M	3/8" Tough-Spots, Maroon	Labels	5	19250.00	CLSWASH-4700	3.00"x2.00" / 76mm x 51mm 250/roll Warning: Do Not Use	Label	5	38500.00
CLST-SPOTS-O	3/8" Tough-Spots, Orange	Labels	5	19250.00	CLSWELL-1000	96 Well Orienter, 4/pk	Microwell Plate Accessories	5	20075.00
CLST-SPOTS-P	3/8" Tough-Spots, Pink	Labels	5	19250.00	CLSWELL-2000	48 Well Orienter, 4/pk	Microwell Plate Accessories	5	20075.00
CLST-SPOTS-R	3/8" Tough-Spots, Red	Labels	5	19250.00	CLSWELL-3000	24 Well Orienter, 4/pk	Microwell Plate Accessories	5	20075.00
CLST-SPOTS-S	3/8" Tough-Spots, Silver	Labels	5	19250.00	CLSWELL-384H	Multi-colored 384 Well Orienter, 2/pk	Microwell Plate Accessories	5	24090.00
CLST-SPOTS-Y	3/8" Tough-Spots, Yellow	Labels	5	19250.00	CLSWELL-384V	Multi-colored 384 Well Orienter, 2/pk	Microwell Plate Accessories	5	24090.00
CLSTS-TNY	3/16" Tough-Spots, White 1,500/roll	Labels	5	32340.00	CLSWELL-96H	Multi-colored 96 Well Orienter, 2/pk	Microwell Plate Accessories	5	24090.00
CLSTTGP-1000	Tough-Tags 1.50 x 0.75" 1,000/roll	Labels	5	36190.00	CLSWELL-96V	Multi-colored 96 Well Orienter, 2/pk	Microwell Plate Accessories	5	24090.00
CLSTTGP-1050	Laser Tough-Tags 1.50 x 0.75" 1,500/pk	Labels	5	50160.00	CLSWEST-1000	96-Well Orienter Stickers, 50/pk	Microwell Plate Accessories	5	23285.00
CLSTTGP-1050-B	Laser Tough-Tags 1.50 x 0.75" 1,500/pk, Blue	Labels	5	50160.00	CLSWIPE-1000	Tough Wipes, 50/Bottle	Laboratory Cleaning Solutions	5	18260.00
CLSTTGP-1050-G	Laser Tough-Tags 1.50 x 0.75" 1,500/pk, Green	Labels	5	50160.00	CLSWPAJ-1000	Adjustable Well Plate Stand	Microwell Plate Accessories	5	65120.00
CLSTTGP-1050-O	Laser Tough-Tags 1.50 x 0.75" 1,500/pk, Orange	Labels	5	50160.00	CLSWPNP-1000	96 Well Pad, 50 sheets	Microwell Plate Accessories	5	10955.00
CLSTTGP-1050-R	Laser Tough-Tags 1.50 x 0.75" 1,500/pk, Red	Labels	5	50160.00	CLSWPNP-2000	48 Well Pad, 50 sheets	Microwell Plate Accessories	5	10955.00
CLSTTGP-1050-Y	Laser Tough-Tags 1.50 x 0.75" 1,500/pk, Yellow	Labels	5	50160.00	CLSWPNP-3000	24 Well Pad, 50 sheets	Microwell Plate Accessories	5	10955.00
CLSTTLB-1000	Tough-Tags 1.28 x 0.50" Blue	Labels	5	26757.50	CLSWPST-1000	Well Plate Stand	Microwell Plate Accessories	5	31955.00
CLSTTLC-1000	Tough-Tags 1.28 x 0.50" Clear	Labels	5	26757.50	CLSWPST-2000	Dual Microplate Stand, angled	Microwell Plate Accessories	5	57200.00
CLSTTLG-1000	Tough-Tags 1.28 x 0.50" Green	Labels	5	26757.50	CLSWPST-3000	Dual Microplate Stand, lat	Microwell Plate Accessories	5	45650.00
CLSTTLO-1000	Tough-Tags 1.28 x 0.50" Orange	Labels	5	26757.50	CLSWRIT-1000	Write Protect, 1.25 x 0.50"	Tape	5	22825.00
CLSTTLR-1000	Tough-Tags 1.28 x 0.50" Red	Labels	5	26757.50	CLSWRIT-2000	Write Protect, 1.60 x 0.63"	Tape	5	22825.00
CLSTTLW-1000	Tough-Tags 1.28 x 0.50" 1,000/roll	Labels	5	19250.00	CLSWSDR-1000	Wash-N-Dry Cover Slip Rack	Microscopy	5	24090.00
CLSTTLW-2016	Laser Tough-Tags 1.28 x 0.50" 2,125/pk	Labels	5	50160.00	CLSWSDR-2000	Wash-N-Dry Cover Slide Rack	Microscopy	5	35200.00
CLSTTLW-2016-B	Laser Tough-Tags 1.28 x 0.50" 2,125/pk, Blue	Labels	5	50160.00	CLSWSDR-BKER	200ml Tall Form Beacker	Microscopy	5	14080.00
CLSTTLW-2016-G	Laser Tough-Tags 1.28 x 0.50" 2,125/pk, Green	Labels	5	50160.00	CLSYLN-1000	.875" x .875" Xylene Resistant Labels	Labels	5	56980.00
CLSTTLW-2016-O	Laser Tough-Tags 1.28 x 0.50" 2,125/pk, Orange	Labels	5	50160.00	CLSYLN-2000	1.69" x .625" Xylene Resistant Labels	Labels	5	56980.00
CLSTTLW-2016-R	Laser Tough-Tags 1.28 x 0.50" 2,125/pk, Red	Labels	5	50160.00					

Microscope slides

- Produced from the finest quality optical white glass
- Better optical properties than most microscope slides
- Packed in high quality plastic boxes
- Dimensions: 25 x 75mm (+/- 0.5mm)
- Thickness: 1mm (+/- 0.05mm)



Cat No	Product	Pack Size	Price
CLS1380-10	Microscope slides, ground edges, 90 corners, plain	1440	28200.00
CLS1380-20	Microscope slides, ground edges, 90 corners, frosted	1440	28200.00
CLS1380-30	Microscope slides, ground edges, 90 corners, frosted 1 end, 2 sides	1440	28200.00

CLS1380-50	Microscope slides, ground edges, 90 corners, frosted aqua	1440	28200.00
------------	---	------	----------

Available Color: Aqua, Blue, Gold, Green, Lilac, Orange, Pink, Tan, White, Yellow

CLS1384-50W	Microscope slides, 45 beveled edges, clipped corners	1440	28200.00
-------------	--	------	----------

CLS1365-10W	Microscope slides for use on automated instruments, ground edges, clipped corners	1440	28200.00
-------------	---	------	----------

Available Color: White

CLS1358A	Charged microscope slides, ground edges, 90 corners	1440	63582.00
----------	---	------	----------

Available Color: Aqua, Blue, Green, Lilac, Orange, Pink, Tan, White, Yellow

CLS1354P	Charged microscope slides, 45 beveled edges, frosted white	1440	41064.00
----------	--	------	----------

Available Color: Frosted Pink, Frosted White, Frosted Yellow

Plain Microscope slides

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Packaging: 72/box, 20 boxes/case (10 gross)



Cat No	Product	Pack Size	Price
CLS1301	Plain, ground edges, 90 corners	1440	16194.00
CLS1321	Plain, ground edges, safety corners	1440	16194.00
CLS1331	Plain, beveled edges with clipped corners	1440	21450.00

Frosted Microscope Slides

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Packaging: 72/box, 20 boxes/case (10 gross)



Cat No	Product	Pack Size	Price
CLS1304	Frosted 1 end, 1 side, ground edges, 90 corners	1440	18681.00
CLS1308	Frosted 1 end, 2 sides, ground edges, 90 corners	1440	19841.00
CLS1324	Frosted 1 end, 1 side, ground edges with safety corners	1440	18682.00
CLS1328	Frosted 1 end, 2 sides, ground edges with safety corners	1440	19842.00
CLS1334	Frosted 1 end, 1 side, beveled edges with clipped corners	1440	23986.00
CLS1338	Frosted 1 end, 2 sides, beveled edges with clipped corners	1440	25147.00

Frosted Color

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Frosted one end, one side Packaging: 72/box, 20 boxes/case (10 gross)



Ground Edge with 90 Corners

Cat No	Product	Pack Size	Price
CLS1304W	Frosted, ground edge, 90 corners	1440	18622.00

Available Color: White, Blue Green, Orange, Pink, Tan, Yellow

Ground Edge with Safety Corners

Cat No	Product	Pack Size	Price
1324W	Frosted, ground edge with safety corners	1440	18622.00

Available Color: White, Aqua, Blue, Gold, Green, Lilac, Orange, Pink, Yellow

Beveled Edge with Clipped Corners

Cat No	Product	Pack Size	Price
CLS1334W	Frosted, beveled edge with clipped corners	1440	23290.00

Available Color: White, Blue Green, Pink, Yellow

Charged Microscope Slides

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Frosted one end, one side Packaging: 72/box, 20 boxes/case (10 gross)



Ground Edge

Cat No	Product	Pack Size	Price
CLS1358	Charged microscope slide, Aqua	1440	63583.00

Available Color: Aqua, Blue, Green, Lilac, Orange, Pink, Tan, White, Yellow

Beveled Edge with Clipped Corners

CLS1354P	Charged microscope slide, beveled	1440	41064.00
----------	-----------------------------------	------	----------

Available Color: Pink, White, Yellow

Microscope Slides with Cavity

Our 25 x 75 x 1mm cavity microscope slides are pre-cleaned and ready to use. Two styles are available: single cavity or double cavity.



Cat No	Product	Pack Size	Price
CLS1341	Microscope slide with single cavity 72/box, 20 boxes/case	1440	71839.00
CLS1341-72	Microscope slide with single cavity	72	5569.00
CLS1344	Microscope slide with double cavity 72/box, 20 boxes/case	1440	140000.00
CLS1344-72	Microscope slide with double cavity	72	11502.00

Cover Glass

Our cover glass is packaged: 1 oz per box, 10 boxes per carton (10 oz). Each lip top box contains a desiccant and is vacuum sealed in aluminum to prevent sticking.




Cat No	Product	Unit	Pack Size	Price
CLS1401-10	18 x 18mm, No. 1	10 oz	248	4664.00
CLS1401-20	18 x 18mm, No. 2	10 oz	168	4664.00
CLS1402-10	20 x 20mm, No. 1	10 oz	192	4664.00
CLS1402-20	20 x 20mm, No. 2	10 oz	154	4664.00
CLS1404-10	22 x 22mm, No. 1	10 oz	168	4664.00
CLS1404-15	22 x 22mm, No. 1.5	10 oz	148	4664.00
CLS1404-20	22 x 22mm, No. 2	10 oz	136	4664.00
CLS1407-10	22 x 32mm, No. 1	10 oz	114	4664.00
CLS1407-20	22 x 32mm, No. 2	10 oz	84	4664.00
CLS1411-10	22 x 40mm, No. 1	10 oz	86	4664.00
CLS1411-20	22 x 40mm, No. 2	10 oz	75	4664.00
CLS1414-10	22 x 50mm, No. 1	10 oz	70	4664.00
CLS1414-20	22 x 50mm, No. 2	10 oz	63	4664.00
CLS1418-10	22 x 60mm, No. 1	10 oz	60	4664.00
CLS1418-20	22 x 60mm, No. 2	10 oz	50	4664.00
CLS1405-10	24 x 24mm, No. 1	10 oz	141	4664.00
CLS1405-20	24 x 24mm, No. 2	10 oz	102	4664.00
CLS1408-10	24 x 32mm, No. 1	10 oz	104	4664.00
CLS1408-20	24 x 32mm, No. 2	10 oz	78	4664.00
CLS1412-10	24 x 40mm, No. 1	10 oz	82	4664.00
CLS1412-20	24 x 40mm, No. 2	10 oz	71	4664.00
CLS1415-10	24 x 50mm, No. 1	10 oz	65	4664.00
CLS1415-15	24 x 50mm, No. 1.5	10 oz	56	4664.00
CLS1415-20	24 x 50mm, No. 2	10 oz	56	4664.00
CLS1419-10	24 x 60mm, No. 1	10 oz	55	4664.00
CLS1419-20	24 x 60mm, No. 2	10 oz	46	4664.00

Multi-Well Plates

96-Well Multi-Well Plates, Non-Treated

- Individually wrapped to protect from air contaminants.
- Available in "U", "V" and Flat bottoms



Round Bottom (U)
Each well is 300uL

Conical Bottom (V)
Each well is 250uL

Flat Bottom (-)
Each Well is 375uL

Cat No	Product	Pack Size	Price
CLS120030	U-bottom, individually wrapped, non-Sterile	50	12066.00
CLS120038	U-bottom, individually wrapped, Sterile	50	13108.00
120130	V-bottom, individually wrapped, non-Steri	50	12066.00
120138	V-bottom, individually wrapped, Sterile	50	13872.00
120330	Flat-bottom, individually wrapped, non-Sterile	50	12835.00
120338	Flat-bottom, individually wrapped, Sterile	50	13872.00
129930	Lid, individually wrapped, non-Sterile	150	21680.00
129938	Lid, individually wrapped, Sterile	150	25129.00

Sealing Film

Designed for the sealing of multi-well plates, our adhesive sealing films minimize evaporation, prevent spillage and eliminate well-to-well contamination.

The secure sealing of all wells eliminates "edge effects" in sensitive ELISA assays. Additionally, our films can be used to seal tissue culture plates for short term storage, incubation, transport, and the containment of biohazardous



Cat No	Product	Pack Size	Price
CLS129940	EZ-Seal sealing ilm, non-sterile	100	POR
CLS129948	EZ-Seal sealing ilm, STERILE	100	POR

Pipette Tips

Standard



Cat No	Product	Pack Size	Price
0.1 - 10uL			
CLS151154	For Pipetman P2-P10, natural, 31mm, 1000/ pack	1000	4231.00
CLS151154R-96	For Pipetman P2-P10, natural, 31mm, 96/rack, 10 racks per case	960	8872.00
CLS151154RS-96	For Pipetman P2-P10, natural, 31mm, 96/rack, 10 racks per case, STERILE	960	11011.00
0.5 - 20uL			
CLS151145	Ultra-micro tip, natural, 45mm, 1000/pack	1000	4322.00
CLS151145R-96	Ultra-micro tip, natural, 45mm, 96/rack, 10 racks per case	960	9213.00
CLS151145RS-96	Ultra-micro tip, natural, 45mm, 96/rack, 10 racks per case, STERILE	960	13377.00
1-200uL			
CLS152143	Universal, yellow, 1000/pack	1000	2730.00
CLS152143R	Universal, yellow, 200/rack, 5 racks/stack	1000	5510.00
CLS151140	For use with MLA & Ovation, natural, 1000/pack	1000	2984.00
CLS151140R	For use with MLA & Ovation, natural, 200/rack, 5 racks/stack	1000	5514.00
CLS151143	For use with Eppendorf, yellow, 1000/pack	1000	2984.00
CLS151241	For use with Oxford 8000 series, natural, 1000/pack	1000	3767.00
CLS151241R	For use with Oxford 8000 series, natural, 200/ rack, 5 racks/stack	1000	12285.00
CLS151142	For use with Oxford SlimLine, natural, 1000/pack	1000	2691.00
CSL10-500uL			
CLS151351	For use with Beckman & Oxford Lancer, natural, 70mm, 1000/pack	1000	6433.00
200-500uL			
CLS151242	For Oxford 8000 series, yellow, 1000/pack	1000	8340.00

Cat No	Product	Pack Size	Price
100-1000uL			
C151146	Universal, blue, 1000/pack	1000	3380.00
C151146R	Universal, blue, 200/rack, 5 racks/stack	1000	6137.00
C151341	For use with Selectapette, natural, 1000/pack	1000	10801.00
200-1000uL			
C151147	For use with MLA pipettors, natural, 1000/pack	1000	5369.00
C151147R	For use with MLA pipettors, natural, 200/rack, 5 racks/stack (1000/rack)	1000	9750.00
C151243	For use with Oxford 8000 series, green, 1000/pack	1000	9100.00
C151144	For use with Oxford Benchmate, Oxford Slimline, Finnpiquette, Labsystems, blue, 1000/pack	1000	3735.00
1000-5000uL (1mL - 5mL)			
C151244	For use with MLA, Oxford, & Pipetman, natural, 116mm, 250/pack	250	3735.00
C151246	For Finnpiquette, Brand, Labsystems, SMI & EDP2, 250/pack	250	9330.00
C151247	For Biohit Proline, Oxford & Eppendorf Research, natural, 250/pack	250	7725.00
1000-10,000uL (1mL - 10mL)			
C151249	For use with Pipetman P10, natural, 250/pack	250	5992.00
C151250	For use with Finnpiquette, Brand, Gilson & Socorex, 100/bag	100	7816.00
C151247	For Biohit Proline, Oxford & Eppendorf Research, natural, 250/pack	250	7725.00

Cat No	Product	Pack Size	Price
C151151RS-96	Universal, graduated, natural, 58mm, 96/rack, 10 racks/case, STERILE	1000	12239.00
(Graduated at 10uL, 20uL, 50uL, 100uL and 200uL)			
100 - 1000uL			
C151152	Universal, graduated, natural, 70mm, 1000/pack	1000	5778.00
C151152R	Universal, graduated, natural, 70mm, 200/rack, 5 racks/stack	1000	8872.00
C151152R-100	Universal, graduated, natural, 70mm, 100/rack, 10 racks/case	1000	11306.00
C151152RS-100	Universal, graduated, natural, 70mm, 100/rack, 10 racks/case, STERILE	1000	16516.00
C151152B	Universal, graduated, blue, 70mm, 1000/pack	1000	5778.00
C151152BR	Universal, graduated, blue, 70mm, 200/rack, 5 racks/stack	1000	8872.00
C151152BR-100	Universal, graduated, blue, 70mm, 100/rack, 10 racks/case	1000	11306.00
C151152BRS-100	Universal, graduated, blue, 70mm, 100/rack, 10 racks/case, STERILE	1000	16516.00
(Graduated at 100uL, 250uL, 500uL, 750 and 1000uL)			
100 - 1250uL			
C151153	Universal, graduated, natural, 84mm (extra long), 1000/pack	1000	5687.00
C151153R-100	Universal, graduated, natural, 84mm (extra long), 100/rack, 10 racks/case	1000	13832.00
(Graduated at 100uL, 250uL, 500uL, 750 and 1000uL)			

Reference (Graduated)

- Low retention
- Universal fit - designed for use on all popular pipettors
- Pipette tips are fully autoclavable (121°C/15 minutes)



Cat No	Product	Pack Size	Price
1 - 200uL			
C151150	Universal, graduated, natural, 50mm, 1000/pack	1000	4982.00
C151150R	Universal, graduated, natural, 50mm, 200/rack, 5 racks/stack	1000	7621.00
C151150R-96	Universal, graduated, natural, 50mm, 96/rack, 10 racks/case	1000	9395.00
C151150RS-96	Universal, graduated, natural, 50mm, 96/rack, 10 racks/case, STERILE	1000	10692.00
C151150Y	Universal, graduated, yellow, 50mm, 1000/pack	1000	4982.00
C151150YR	Universal, graduated, yellow, 50mm, 200/rack, 5 racks/stack	1000	7621.00
C151150YR-96	Universal, graduated, yellow, 50mm, 96/rack, 10 racks/case	1000	9395.00
C151150YRS-96	Universal, graduated, yellow, 50mm, 96/rack, 10 racks/case, STERILE	1000	20692.00
(Graduated at 10uL, 20uL, 50uL, 100uL and 200uL)			
1 - 300uL			
C151151	Universal, graduated, natural, 58mm, 1000/pack	1000	5551.00
C151151R	Universal, graduated, natural, 58mm, 200/rack, 5 racks/stack	1000	7962.00
C151151R-96	Universal, graduated, natural, 58mm, 96/rack, 10 racks/case	1000	11488.00

Filter Tips (Sterile)

Globe Scientific's line of filter pipette tips are designed for accuracy and aerosol protection. They are universal and fit most popular brands of pipettors.



All filter tips are sterile and packaged in lip-top boxes. The heavy duty box is produced from polypropylene and can be autoclaved for repeated use.

- Sterile
- Rnase, DNase and Pyrogen Free
- Universal Fit - designed for use on all popular pipettors
- Pipette tips are not autoclavable

Pipette Tips: Polypropylene (PP) Filter: Polyethylene (PE)

Cat No	Size	Length	Graduations	Unit	Price
C153800	0.1-10uL	31mm	Non-graduated	96 x 10 Rack	26208.00
C153805	0.5-10uL	45mm	Non-graduated	96 x 10 Rack	24570.00
C153810	1-20uL	50mm	10, 20, 50, 100 & 200uL	96 x 10 Rack	24570.00
C153815	1-100uL	50mm	10, 20, 50, 100 & 200uL	96 x 10 Rack	24570.00
C153820	1-200uL	50mm	10, 20, 50, 100 & 200uL	96 x 10 Rack	24570.00
C153825	1-250uL	58mm	10, 20, 50, 100 & 200uL	96 x 10 Rack	24570.00
C153830	1-1000uL	70mm	100, 250, 500, 750 & 1000uL	100 x 10 Rack	27927.00
C153835	1-1000uL	86mm	100, 250, 500, 750 & 1000uL	100 x 10 Rack	27927.00

COMPOSITE LAB LINE PVT LTD

UGSR-1, Ansal Plaza, Vaishali

Delhi, NCR, Ghaziabad - 201010 (India)

Ph. : 0120-4335636 Mail : sales@compositelablin.com

Visit.: www.compositelabline.com

Reg Office :

A-1, OKKP Complex, Near I.T.Crosing

Faizabad Road, Babujanj

Lucknow - 226020 (India)

Ph. : +919580613355, 0522 - 4068342

Mail : sales@compositelablin.com



Authorised Dealer:

