



Professional in

- CLS Clinical / Pathology Product's

Composite Lab Line Pvt Ltd

Wide Mouth Bottles

Wide Mouth, Square Bottles (Rigid), LDPE
Bottle = Low Density Polyethylene (LDPE),
Cap = Polypropylene (PP)

- Graduated
- Dual-cap feature. Each bottle includes a screw cap and a plug insert, providing a completely leak-resistant seal.
- Tagging points for use with optional wire security ties
- Food grade



Cat No	Product	Pack Size	Price
C600608	25mL (0.8oz) wide mouth, square, LDPE	100/case	29484.00
C600609	50mL (1.5 oz) wide mouth, square, LDPE	100/case	33196.00
C600610	100mL (3 oz) wide mouth, square, LDPE	100/case	43734.00
C600611	250mL (8 oz) wide mouth, square, LDPE	50/case	25825.00
C600612	500mL (16 oz) wide mouth, square, LDPE	25/case	17226.00
C600613	1000mL (32 oz) wide mouth, square, LDPE	20/case	22577.00
C600614	2000mL (64 oz) wide mouth, square, LDPE	20/case	41796.00

Wide Mouth Bottles

Wide Mouth, Round Bottles (Extra Flexible), LDPE
Bottle = Low Density Polyethylene (LDPE),
Cap = Polypropylene (PP)

- Graduated
- Dual-cap feature. Each bottle includes a screw cap and a plug insert, providing a completely leak-resistant seal.
- Tagging points for use with optional wire security ties
- Food grade



Cat No	Product	Pack Size	Price
C600408	50mL (1.5 oz) wide mouth, LDPE	100/case	28746.00
C600409	100mL (3 oz) wide mouth, LDPE	100/case	30739.00
C600410	250mL (8 oz) wide mouth, LDPE	50/case	20447.00
C600411	500mL (16 oz) wide mouth, LDPE	20/case	10346.00
C600412	1000mL (32 oz) wide mouth, LDPE	20/case	17007.00
C600413	2000mL (64 oz) wide mouth, LDPE	20/case	36882.00

Wide Mouth Bottles

Wide Mouth, Round Bottles (Flexible),
LDPE

Bottle = Low Density Polyethylene (LDPE),
Cap = Polypropylene (PP)

- Graduated
- Deep threaded neck for excellent leak-resistance
- Food Grade



Cat No	Product	Pack Size	Price
C601608	50mL (1.5 oz) wide mouth, LDPE	100/case	40950.00
C601610	100mL (3 oz) wide mouth, LDPE,	100/case	52361.00
C601612	250mL (8 oz) wide mouth, LDPE	50/case	37237.00
C601614	500mL (16 oz) wide mouth, LDPE	25/case	22113.00
C601616	1000mL (32 oz) wide mouth, LDPE	17/case	29265.00
C601618	2000mL (64 oz) wide mouth, LDPE	15/case	38247.00

Wide Mouth Bottles

Wide Mouth, Round Bottles (Rigid), PP
Bottle and Cap = Polypropylene (PP)

- Completely rigid and rupture-resistant
- Graduated
- Deep threaded neck for excellent leak-resistance
- Food grade
- Autoclavable at 121°C. Do not autoclave with caps on.



Cat No	Product	Pack Size	Price
C601621	50mL (1.5 oz) wide mouth, PP	100/case	46628.00
C601623	100mL (3 oz) wide mouth, PP	100/case	57603.00
C601625	250mL (8 oz) wide mouth, PP	50/case	39039.00
C601627	500mL (16 oz) wide mouth, PP	25/case	25252.00
C601629	1000mL (32 oz) wide mouth, PP	17/case	29647.00
C601631	2000mL (64 oz) wide mouth, PP	15/case	44034.00

Narrow Mouth Bottles

Narrow Mouth, Round Bottles (Rigid), PP
Bottle and Cap = Polypropylene (PP)

- Completely rigid and rupture-resistant
- Graduated
- Tagging points for use with optional wire security ties
- Food grade
- Autoclavable at 121°C. Do not autoclave with caps on



Cat No	Product	Pack Size	Price
C601595	50mL (1.5 oz) narrow mouth, PP	100/case	34206.00
C601597	100mL (3 oz) narrow mouth, PP	100/case	41059.00
C601599	250mL (8 oz) narrow mouth, PP	50/case	3558.00
C601601	500mL (16 oz) narrow mouth, PP	25/case	18181.00
C601603	1000mL (32 oz) narrow mouth, PP	20/case	25771.00

Narrow Mouth Bottles

Narrow Mouth, Round Bottles (Flexible), LDPE
Bottle = Low Density Polyethylene (LDPE),
Cap = Polypropylene (PP)

- Graduated
- Dual-cap feature. Each bottle includes a screw cap and a plug insert, providing a completely leak-resistant seal
- Tagging points for use with optional wire security ties
- Food grade



Cat No	Product	Pack Size	Price
C60031	50mL (1.5oz) narrow mouth, LDPE	100/case	30849.00
C600319	100mL (3 oz) narrow mouth, LDPE	100/case	31449.00
C600323	250mL (8 oz) narrow mouth, LDPE	50/case	22604.00
C600324	500mL (16 oz) narrow mouth, LDPE	25/case	17144.00
C600325	1000mL (32 oz) narrow mouth, LDPE	20/case	20884.00
C600326	2000mL (64 oz) narrow mouth, LDPE	20/case	42369.00

Wash Bottles

Wash Bottle with Integrated Spout, LDPE

Bottle = Low Density Polyethylene (LDPE),
Cap = Polypropylene (PP)

- Easy grip, stable shape
- Dispensing tube with closure cap
- Dual-cap feature. Each wash bottle includes a screw cap and a plug insert, providing a completely leak-resistant seal.
- Tagging points for use with optional wire security ties
- Food grade



Cat No	Product	Pack Size	Price
C601633	250mL (8 oz) wash bottle, LDPE	100/case	8872.00
C601634	500mL (16 oz) wash bottle, LDPE	10/case	10756.00

Centrifuge Tubes

12 mL Flared-Top Urine Centrifuge Tubes

This line of centrifuge tubes has been designed specifically for urinalysis procedures. The tubes offer a flared-opening which is designed to make pouring of easy.

- Dimensions: 20 x 115mm
- Polystyrene (PS): 1500 RCF



Cat No	Product	Pack Size	Price
C112017	12mL polystyrene (PS), graduated, conical bottom, Ames-Type	1500	26567.00
C112015	12mL polystyrene (PS), graduated, round bottom, Kova-Type	2000	30272.00
C112030	12mL polystyrene (PS) with sediment bulb, graduated, Whale-Type	1500	24383.00

12 mL Centrifuge Tubes

- Dimensions: 16 x 100mm
- Graduated at 0.5, 1, 2.5, 5 and 10mL
- Polystyrene (PS): 1500 RCF
- Polypropylene (PP): 3000 RCF



Cat No	Product	Pack Size	Price
C112010	12mL (16 x100mm), conical bottom, graduated, polystyrene (PS)	2000	31354.00
C112010-500	12mL (16 x100mm), conical bottom, graduated, polystyrene (PS)	500	8130.00
C112040	12mL (16 x100mm), conical bottom, graduated, polypropylene (PP)	2000	26667.00
C112040-500	12mL (16 x100mm), conical bottom, graduated, polypropylene (PP)	500	7775.00
C112009	12mL (16 x100mm), conical bottom, polystyrene (PS) w/ separate plug cap	1000	21061.00

10 mL Centrifuge Tubes

- Dimensions: 16 x 100mm
- Printed graduations
- Spin Rate: 3000 RCF
- Tube and Cap: Polypropylene (PP)



Cat No	Product	Pack Size	Price
C6292	10mL (16 x100mm), with separate natural screw cap, 500/bg, 2 bags/cs	1000	28137.00
C6293	10mL (16 x100mm), with attached yellow screw cap, STERILE, 100/bag, 10 bags/case	1000	28878.00

15mL Centrifuge Tubes

- Dimensions: 16.6 x 120mm
- Graduated from 0.5mL to 15mL in 0.5mL increments
- Polystyrene (PS): 1500 RCF
- Polypropylene (PP): 3000 RCF




Cat No	Product	Pack Size	Price
C6260	15mL polystyrene (PS), 250/pack, 4 packs/case	1000	12762.00
C6260-250	15mL polystyrene (PS), 250/pack	250	4586.00
C6265	15mL polystyrene (PS), 100/dispenser box, 10 boxes/case	1000	14414.00
C6261	15mL polypropylene (PP), 250/pack, 4 packs/case	1000	12244.00
C6266	15mL polypropylene (PP), 100/dispenser box, 10 boxes/case	1000	14004.00
C6262	15mL polystyrene (PS), with separate plug cap	1000	17654.00
C6263	Plug cap for 15mL tubes	1000	4868.00

Stainless Steel Racks

These specially shaped "Z" racks are rust proof and will hold up to repeated use. The wide base provides stability and prevents tipping. When not in use, the racks

stack neatly inside one another to conserve benchtop space

- Autoclavable



Cat No	Product	Pack Size	Price
C457100	For 12 & 13mm tubes, 4 x 4 x 1.75", 25-place	Each	3021.00
C457105	For 12 & 13mm tubes, 8 x 4 x 2", 50-place	Each	4067.00
C457200	For 15, 16 & 17mm tubes, 4.5 x 4.25 x 1.5", 25-place	Each	3021.00
C457205	For 15, 16 & 17mm tubes, 8.75 x 4.5 x 2", 50-place	Each	4172.00

15 mL Centrifuge Tubes

- 15 mL Centrifuge Tubes
- Dimensions: 17 x 120mm
 - Cap is made from High Density Polyethylene (HDPE)
 - Sterile tubes are RNase/Dnase and Pyrogen free.
 - White graduations and writing area.

Spin Rate:

Polystyrene (PS) - 1700 RCF (3000 x g)

Acrylic (AC) - 5000 RCF (6,160 x g)

Polypropylene (PP) - 10,000 RCF (12,320 x g)



Cat No	Product	Pack Size	Price
C6268	15mL polystyrene (PS) centrifuge tube, non-sterile, with separate screwcap, bulk	500	22463.00
C6272	15mL polystyrene (PS) centrifuge tube, STERILE, 25/bag, 500/case	500	27389.00
C6274	15mL polystyrene (PS) centrifuge tube, STERILE, 50/rack, 500/case	500	34516.00
C6267	15mL acrylic (AC) centrifuge tube, non-sterile, with separate screwcap, bulk	500	27418.00
C6270	15mL acrylic (AC) centrifuge tube, STERILE, 25/bag, 500/case	500	28924.00
C6273	15mL acrylic (AC) centrifuge tube, STERILE, 50/rack, 500/case	500	26554.00
C6264	15mL polypropylene (PP) centrifuge tube, non-sterile, with separate screwcap, bulk	500	22463.00
C6269	15mL polypropylene (PP) centrifuge tube, STERILE, 25/bag, 500/case	500	25389.00
C6271	15mL polypropylene (PP) centrifuge tube, STERILE, 50/rack, 500/case	500	32305.00

50 mL Self-Standing Polypropylene - Printed Graduations

- Dimensions: 30 x 115 mm
- Cap is made from High Density Polyethylene (HDPE)
- Sterile tubes are RNase/Dnase and Pyrogen free
- White graduations and writing area
- Autoclavable (without screwcap)



Cat No	Product	Pack Size	Price
C255	50mL, non-sterile, bulk	500	25994.00
C256	50mL, STERILE, 25/bag, 500/case	500	31545.00
C255MF	50mL, metal free, neutral cap, non-sterile	500	25994.00

Microcentrifuge Tubes



50 mL Centrifuge Tubes

- 50 mL Centrifuge Tubes
- Dimensions: 29.5 x 114mm
 - Cap is made from High Density Polyethylene (HDPE)
 - Sterile tubes are RNase/Dnase and Pyrogen free
 - White graduations and writing area
 - Polypropylene tubes are autoclavable (without screwcap)

Spin Rate:

Polystyrene (PS) - 1700 RCF (3000 x g)

Acrylic (AC) - 5000 RCF (6,160 x g)

Polypropylene (PP) - 10,000 RCF (12,320 x g)



Cat No	Product	Pack Size	Price
C6240	50mL polypropylene (PP) centrifuge tube, non-sterile, 25/bag, 500/case	500	28319.00
C6241	50mL polypropylene (PP) centrifuge tube, STERILE, 25/bag, 500/case	500	32182.00
C6242	50mL polypropylene (PP) centrifuge tube, non-sterile, 25/rack, 500/case	500	42315.00
C6243	50mL polypropylene (PP) centrifuge tube, STERILE, 25/rack, 500/case	500	49126.00
C6244	50mL polypropylene (PP) centrifuge tube, non-sterile, with separate screwcap, bulk	500	26521.00
C6253	50mL polystyrene (PS) centrifuge tube, non-sterile, with separate screwcap, bulk	500	28146.00
C6254	50mL polystyrene (PS) centrifuge tube, STERILE, 25/bag, 500/case	500	35485.00
C6254R	50mL polystyrene (PS) centrifuge tube, STERILE, 25/rack, 500/case	500	52188.00

Non-Graduated, Attached Cap

Cat No	Product	Pack Size	Max. RPM	Price
C10540	500uL (0.5mL) polyethylene, non-autoclavable	1000	12000	6383.00

Graduated, Attached Cap

C10550	500uL (0.5mL), Eppendorf Type, polypropylene	1000	18000	8800.00
C11550	600uL (0.6mL), polypropylene	1000	30000	9113.00
C11558	1.5mL, polypropylene	1000	26000	3894.00
C11558B	1.5mL, polypropylene, blue	1000	18000	3894.00
C11558P	1.5mL, polypropylene, pink	1000	18000	3894.00
C11558Y	1.5mL, polypropylene, yellow	1000	18000	3949.00
C11553A	1.5mL, polypropylene, AMBER	1000	18000	10769.00
C11568	2mL, polypropylene, round bottom, 1000/pack	1000	18000	4627.00
C11552	2mL, polypropylene, 500/pack, 2 packs/box	1000	18000	10528.00
C11552A	2mL, polypropylene, AMBER, 500/pack, 10 packs/case	5000	18000	50600.00

Screw Cap Microcentrifuge Tubes

Tube = Polypropylene

- Cap = Polypropylene
- O-Ring = Nitrile

Centrifuge Maximum RCF = 18,000 RCF



Conical Bottom


Cat No	Product	Pack Size	Price
C11690	0.5mL tube with o-ring cap, STERILE	1000	42137.00
C11692	0.5mL tube only, non-sterile	1000	9400.00
C11710	1.5mL tube with o-ring cap, STERILE	1000	44089.00
C11712	1.5mL tube only, non-sterile	1000	10751.00
C11730	2.0mL tube with o-ring cap, STERILE	1000	44435.00
C11732	2.0mL tube only, non-sterile	1000	11070.00

Self Standing

Cat No	Product	Pack Size	Price
C11700	0.5mL tube with o-ring cap, STERILE	1000	40053.00
C11702	0.5mL tube only, non-sterile	1000	7152.00
C11720	1.5mL tube with o-ring cap, STERILE	1000	40781.00
C11722	1.5mL tube only, non-sterile	1000	7730.00
C11740	2.0mL tube with o-ring cap, STERILE	1000	44876.00
C11742	2.0mL tube only, non-sterile	1000	11274.00

Cryogenic Vials

CryoClear™ - Internal & External Threads



Internal Thread Cap
OD: 12.50 mm

Cat No	Volume	Dimensions (L x W)	Style	Bottom	Unit	Price
C301	1 mL	42.7 mm x 12.5 mm	Self Standing	Conical	500	33200.00
C302	2 mL	48.0 mm x 12.5 mm	Self Standing	Round	500	33200.00
C303	2 mL	47.0 mm x 12.5 mm	--	Round	500	33200.00
C304	3 mL	75.1 mm x 12.5 mm	--	Round	500	33200.00
C305	4 mL	76.1 mm x 12.5 mm	Self Standing	Round	500	33200.00
C306	5 mL	89.3 mm x 12.5 mm	--	Round	500	33200.00
C308	5 mL	90.4 mm x 12.5 mm	Self Standing	Round	500	33200.00




External Thread Cap
OD: 14.30 mm

Cat No	Volume	Dimensions (L x W)	Style	Bottom	Unit	Price
3001	1 mL	42.7 mm x 12.5 mm	Self Standing	Conical	500	33200.00
3002	2 mL	48.0 mm x 12.5 mm	Self Standing	Round	500	33200.00
3003	2 mL	47.0 mm x 12.5 mm	-	Round	500	33200.00
3004	3 mL	75.1 mm x 12.5 mm	-	Round	500	33200.00
3005	4 mL	76.1 mm x 12.5 mm	Self Standing	Round	500	33200.00
3006	5 mL	89.3 mm x 12.5 mm	-	Round	500	33200.00
3008	5 mL	90.4 mm x 12.5 mm	Self Standing	Round	500	33200.00

Workstation Racks

- 40-Place (10 x 4 format)
- Produced from Polycarbonate (PC)
- Anti-slip rubber feet
- Autoclavable at 121°C, 15 psig (1 bar) for 20 minutes
- Temperature range: -196°C to +121°C
- Dimensions: 215 x 100 x 25mm




Cat No	Product	Pack Size	Price
C349	Workstation	5 Unit	19893.00

Available Color: Blue, Green, Red

BioBOX™ Storage Boxes with Transparent Lids

- Polycarbonate (PC)
- Temperature range: -196°C to +121°C
- Transparent lid for convenient viewing of contents
- The lid is printed with numbers for easy sample identification
- Lid has two angled corners to ensure the correct alignment of the well numbers
- Autoclavable at +121°C, 15 psig (1 bar) for 20 minutes
- Fits in stainless steel storage racks
- A plastic tube picker is supplied with each storage box



Cat No	Product	Pack Size	Price
C339	25-place for 1 & 2mL tubes	5	11465.00
C340	81-place for 1 & 2mL tubes	5	10198.00
C342	81-place for 3, 4 & 5mL external & 5mL internal	5	210297.00
C350	100-place for 1 & 2mL tubes	5	10551.00

Available Color: Blue, Green, Red, Yellow

Carboys

Store, mix and dispense bulk quantities of liquids, solvents & cleaning solutions. The carboys are available in four sizes and feature a standardized industrial packaging design. Not only can they be used in industrial packaging, but they are also useful in everyday laboratory practices. Produced from High Density Polyethylene (HDPE), they feature comfortable handles and secure seal caps.




Cat No	Product	Pack Size	Price
C61780	Carboys, 2 Liter (1/2 Gallon), HDPE	Each	3200.00
C61781	Carboys, 5 Liter (1.3 Gallon), HDPE	Each	5600.00
C61782	Carboys, 10 Liter (2.5 Gallon), HDPE	Each	8800.00
C61783	Carboys, 20 Liter (5 Gallon), HDPE	Each	16800.00
C61660	Carboys, with Spigot, Heavy Duty, 5 Liter (1.33 Gallon), HDPE, Spigot Included	1/Pack	15200.00
C61662	Carboys, with Spigot, Heavy Duty, 10 Liter (2.62 Gallon), HDPE, Spigot Included	1/Pack	20000.00
C61664	Carboys, with Spigot, Heavy Duty, 25 Liter (6.5 Gallon), HDPE, Spigot Included	1/Pack	38000.00
C61666	Carboys, with Spigot, Heavy Duty, 50 Liter (13 Gallon), HDPE, Spigot Included	1/Pack	62000.00

Serological Pipettes

One-Piece Plastic Serological Pipettes

Features:

- Smooth one-piece construction
- Calibration is certified for accuracy
- Aerosol barrier filter plug
- Ergonomic design for use in tight spaces
- Reverse graduations (except for 1mL & 2mL sizes)
- Negative graduations (not on the 25mL size)
- Fits all major brands of pipettors
- Safe alternative to glass pipettes
- Calibrated to deliver (TD). Meets ASTM E934



Technical Specifications

- Material: Polystyrene (PS)
- Plug: Bonded polyolefin fiber
- Inaccuracy: <2% at full volume
- Pyrogen-free, LAL tested (0.5Eu/mL), (AAMI ST72) and Alternatives to Batch Testing (2002/(R)2010), non-hemolytic (ISO 10993-4), non-cytotoxic (ISO 10993-5)
- Sterility: Sterilized by gamma irradiation
- RNase-free and DNase-free
- Sterility assurance level (SAL): 10⁻⁶
- Sizes: 1mL, 2mL, 5mL, 10mL and 25mL

Cat No	Product	Pack Size	Price
C700	1mL, individually wrapped, sterile	500	17390.00
C710	1mL, 25 per pack, sterile	1000	31827.00
C715	1mL, 25 per pack, non-sterile	1000	31617.00
C720	2mL, individually wrapped, sterile	500	18573.00
C730	2mL, 25 per pack, sterile	1000	35435.00
C735	2mL, 25 per pack, non-sterile	1000	31840.00
C740	5mL, individually wrapped, sterile	250	13167.00
C750	5mL, 25 per pack, sterile	375	18746.00
C755	5mL, 25 per pack, non-sterile	375	16493.00
C760	10mL, individually wrapped, sterile	200	11334.00
C770	10mL, 25 per pack, sterile	250	12881.00
C775	10mL, 25 per pack, non-sterile	250	11058.00
C780	25mL, individually wrapped, sterile	100	11739.00
C784	25mL, individually wrapped, non-sterile	100	9442.00
C380	Diamond SeroLogic™ portable pipette controller	1	56875.00

50mL Serological Pipettes

- Material: Polystyrene
- Fiber filter plug
- Non-pyrogenic
- Ascending and descending graduations
- Negative graduations for extra working volume
- Sterilized by gamma irradiation
- Sterility assurance level (SAL): 10⁻⁶
- Individually wrapped in thermoformed paper/plastic peel packaging
- Safe alternative to glass pipettes
- Calibrated to deliver (TD) at +/-2%. Meets ASTM E934



Cat No	Product	Pack Size	Price
C820	50mL serological pipette, sterile, individually wrapped	100	21157.00

Open-End Pipettes

- Material: Polystyrene (PS)
- Fiber filter plug
- Sterilized by gamma irradiation
- Sterility assurance level (SAL): 10⁻⁶
- Non-pyrogenic
- Individually wrapped in thermoformed paper/plastic peel packaging
- Safe alternative to glass pipettes
- Industry standard color coding meets ASTM E1380
- Calibrated to deliver (TD) at +/-2%. Meets ASTM E934



Cat No	Product	Pack Size	Price
C860	1mL open end pipette, sterile, individually wrapped	1000	POR
C864	2mL open end pipette, sterile, individually wrapped	800	POR
C868	5mL open end pipette, sterile, individually wrapped	300	POR
C872	10mL open end pipette, sterile, individually wrapped	200	POR

Bacteriological Pipettes

Designed for use in the dairy field, these sterile pipettes meet the A.P.H.A. standard for the examination of dairy products. They are available in 1.1mL and 2.2mL sizes.

- Material: Polystyrene (PS)
- Fiber filter plug
- Reverse graduations
- Negative graduations
- Sterilized by gamma irradiation
- Sterility assurance level (SAL): 10⁻⁶
- Non-pyrogenic
- Individually wrapped in thermoformed paper/plastic peel packaging
- Safe alternative to glass pipettes
- Calibrated to deliver (TD) at +/-2%. Meets ASTM E934



Cat No	Product	Pack Size	Price
C825	1.1mL bacteriological pipette, individually wrapped	1000	POR
C826	1.1mL bacteriological pipette, 25/pack	1000	POR
C828	2.2mL bacteriological pipette, individually wrapped	800	POR
C829	2.2mL bacteriological pipette, 25/pack	800	POR

Aspirating Pipettes


- Material: Polystyrene (PS)
- Transparent, no graduation marks
- Non-plugged
- Sterilized by gamma irradiation
- Sterility assurance level (SAL): 10⁻⁶
- Non-pyrogenic
- Individually wrapped in thermoformed paper/plastic peel packaging
- Safe alternative to glass pipettes



Cat No	Product	Pack Size	Price
C830	1mL aspirating pipette, sterile, individually wrapped	1000	35000.00
C834	2mL aspirating pipette, sterile, individually wrapped	800	31000.00
C838	5mL aspirating pipette, sterile, individually wrapped	300	16800.00
C842	10mL aspirating pipette, sterile, individually wrapped	200	13200.00
C846	25mL aspirating pipette, sterile, individually wrapped	200	23200.00
C850	50mL aspirating pipette, sterile, individually wrapped	100	19500.00

Large Capacity Rack

This line of plastic test tube racks are shipped flat in order to reduce shelf space and then simply snap together for use. The special design allows all samples to be incubated at the same temperature. These racks can be used in freezers without deforming or going brittle. They will not rot, rust or deform. Alphanumeric markings are molded into the racks to provide easy test tube identification. They are available in four colors for color coding test samples, or simply for aesthetics in the modern laboratory. They are stackable when empty.



Easy to assemble

- Produced from polypropylene (PP)
- Autoclavable at 121°C for 20 minutes
- Non-loading in a water bath

For Microcentrifuge Tubes

100 holes, rack dimensions: 263 x 109 x 45mm

Cat No	Product	Pack Size	Price
C56303	Microcentrifuge Tube Rack	Each	2893.00

For 12 & 13mm Tubes

90 holes, 6 x 15 format, rack dimensions: 246 x 105 x 64mm

C56403	13mm size, 90-place	Each	1342.00
--------	---------------------	------	---------

For 16 & 17mm Tubes

60 holes, 5 x 12 format, rack dimensions: 246 x 105 x 72mm

C56503	16mm size, 60-place	Each	1342.00
--------	---------------------	------	---------

For 20 & 21mm Tubes

40 holes, 4 x 10 format, rack dimensions: 246 x 105 x 72mm

C56603	20mm size, 40-place	Each	1342.00
--------	---------------------	------	---------

For 25mm Tubes

40 holes, 4 x 10 format, rack dimensions: 295 x 125 x 85mm

C56703	25mm size, 40-place, white	Each	3831.00
--------	----------------------------	------	---------

For 30mm Tubes

24 holes, 3 x 8 format, rack dimensions: 300 x 112 x 85mm

C56803	30mm size, 24-place, white	Each	4122.00
--------	----------------------------	------	---------

Color Available: White, Blue, Red, Yellow

10 & 12 Place Racks

These traditional test tube holders are lightweight, compact and unbreakable. They take up less space on the benchtop to provide a larger work area. They are produced from durable polypropylene (PP) and can be autoclaved at 121°C for 20 minutes.



Cat No	Product	Pack Size	Price
C56105	For 12mm & 13mm tubes, 12-place, (2 x 6 configuration), white	Each	3239.00
C56106	For 16 to 18mm tubes, 10-place, (1 x 10 configuration), white	Each	2802.00

10-Place Polystyrene Rack

This rack has been designed for use with the Sedi-Rate™ and Sedigren™ ESR Systems as well as holding a variety of test tubes.



Cat No	Product	Pack Size	Price
C56105	10-place rack	Each	3500.00

Grip Racks

These specially designed test tube racks feature rubber grips that hold test

tubes in place when decanting liquid (see photo A). Item #C56925 & #C56926 feature a unique tube ejector that allows the user to easily eject all tubes into a waste receptacle (see photo B)

- 50-place (5 x 10 configuration)
- Single rack accommodates a variety of tube sizes
- Holds tubes securely while decanting



Grip Rack

For up to 17mm tubes. Dimensions: 254 x 140 x 51mm.

- Low profile rack allows for easy viewing of labels on tubes.

Cat No	Product	Pack Size	Price
C56920	Blue	Each	2052.00
C56921	Orange	Each	2052.00


Grip Rack with Tube Ejector

For up to 16mm tubes. Dimensions: 260 x 140 x 82mm.

Cat No	Product	Pack Size	Price
C56925	Blue	Each	3125.00
C56926	Orange	Each	3125.00

Centrifuge Tube Racks

Designed to hold both 15mL and 50mL centrifuge tubes in one rack! The rack ships flat and is simple to assemble. It holds twenty 50mL and thirty 15mL centrifuge tubes.



Cat No	Product	Pack Size	Price
C56930	For 15mL & 50mL centrifuge tubes, orange Dimensions: 209 x 171 x 63mm	Each	1500.00

COMPOSITE LAB LINE PVT LTD

UGSR-1, Ansal Plaza, Vaishali

Delhi, NCR, Ghaziabad - 201010 (India)

Ph. : 0120-4335636 Mail : sales@compositelablin.com

Visit.: www.compositelabline.com

Reg Office :


A-1, OKKP Complex, Near I.T.Crosing

Faizabad Road, Babujanj

Lucknow - 226020 (India)

Ph. : +919580613355, 0522 - 4068342

Mail : sales@compositelablin.com



Authorised Dealer:

--