



Professional in

- CLS Cell Culture Plasticware

Composite Lab Line Pvt Ltd

Cell Culture Flask with Filter or Vent/Close Caps (Cell Culture Treated)

For more than a century, cells have provided humanity with deep and fascinating insights into the nature. Being a common laboratory technique, cell culture is one of the most critical aspects in biology-related research and industrial. We provide our customers with a wide range of high quality cell culture products, including cell culture flasks, cell culture dishes, microplates, cryovials and other essential accessories. All our products meet the international guidelines including those of ISO 9001 and USP

Manufactured from optically clear virgin polystyrene, Ergonomic design to facilitate easy handling and minimize contamination, Short & wide neck with angled design allow easy access, Plugs & filter caps are available for all flask models, Sterilized, Non – pyrogenic & Non – cytotoxic, Cell growth area ranging from 25 cm² to 175 cm²



Cat No.	Description	Pack	Unit Price	HSN
CLSCF001	Cell Culture Flask 25Cm ² , cap Vent close	100	8362.00	39269099
CLSCF002	Cell Culture Flask 25 Cm ² , Fiter cap,	100	8461.00	39269099
CLSCF003	Cell Culture Flask 75 Cm ² , cap vent close,	100	5678.00	39269099
CLSCF004	Cell Culture Flask 75 Cm ² , Fiter cap,	100	8561.00	39269099
CLSCF005	Cell Culture Flask vent/close cap 175/225Cm ²	100	10932.00	39269099
CLSCF006	Cell Culture Flask 25cm with Filter, sterile	200	19494.00	39269099
CLSCF007	Cell Culture Flask 75cm with Filter, sterile	100	15214.00	39269099
CLSCF008	Cell Culture Flask 175Cm ² with Filter, sterile	30	8974.00	39269099
CLSCF009	Cell Culture Flask 225Cm ² with Filter, sterile	30	20123.00	39269099
CLSCF010	Cell Culture Flask 25 vent close, sterile	200	16711.00	39269099
CLSCF011	Cell Culture Flask 75 vent close, sterile	100	15214.00	39269099
CLSCF012	Cell Culture Flask 175 vent close, sterile	30	8974.00	39269099
CLSCF013	Cell Culture Flask 225 vent close, sterile	30	15587.00	39269099



Cell Culture Dishes (Cell Culture Treated)

CLS cell culture products allow researchers to acquire high resolution microscopic images of cells in 35 mm culture dishes or 6 well plates. We provide two types of confocal products with cover glasses of different thicknesses.


Cat No.	Description	Pack	Unit Price	HSN
CLSCD001	TC Dish 35X10 Delta Sterile	500	19681.00	39269099
CLSCD002	35X10mm Suspension Dish	500	33680.00	39269099
CLSCD003	TC Dish 60X15 Delta Sterile	400	11282.00	39269099
CLSCD004	TC Dish 92X17 Delta Sterile	400	11282.00	39269099
CLSCD005	TC Dish 92X17 Sterile	400	11282.00	39269099



Multidishes (Cell Culture Treated)

Designed to prevent cross contamination, Packed individually, Optically clear and uniform surface for microscopy Sterilized & Treated plates, Crystal – grade polystyrene (PS), Non – pyrogenic, Flat bottoms

Cat No.	Description	Pack	Unit Price	HSN
CLSM001	6 Well Multidish W/Lid	75	10183.00	39269099
CLSM002	6 Well Multidish W/Lid	85	16274.00	39269099
CLSM003	8 Well Dish Ion Sterile Lid	100	68456.00	39269099
CLSM004	12 Well Multidish W/Lid TC Treated	75	11232.00	39269099
CLSM005	24 Well Multidish W/Lid TC Treated	75	12156.00	39269099
CLSM006	24 Well Multidish W/Lid TC Treated	85	18170.00	39269099
CLSM007	48 Well Multidish W/Lid TC Treated	75	20592.00	39269099
CLSM008	48 Well Multidish W/Lid	85	22776.00	39269099
CLSM009	96 Well Multidish W/Lid TC Treated	50	22776.00	39269099
CLSM010	06 Well Multidish PS/Lid Non Treated	100	13451.00	39269099
CLSM011	12 Well Multidish PS/Lid Non Treated	100	16758.00	39269099
CLSM012	24 Well Multidish PS/Lid Non Treated	100	14333.00	39269099
CLSM013	48 Well Multidish PS/Lid Non Treated	75	18462.00	39269099
CLSM014	96 Well Multidish PS/Lid Non Treated	108	22157.00	39269099



Cell Scrapper

Cell Scrapers are specially designed to facilitate collection attached cells on culture dish. Sterilization by irradiation and individual packing eliminate possible chances of contamination during cell culture and scraper handling. Rotatable blades, Sterilized, Non – pyrogenic, Non – cytotoxic


Cat No.	Description	Pack	Unit Price	HSN
CLSCS001	CELL SCRAPER, 23CM, SMALL	250	62773.00	39269099
CLSCS002	CELL SCRAPER, 32CM, LARGE	250	70736.00	39269099



Serological Pipetes

Serological Pipettes are classified by sample volume. Ascending & descending scales facilitate reading both dispensing and remaining volume.

Cat No.	Description	Pack	Unit Price	HSN
CLSSP001	1mL Pipette Individually Wrapped Paper/Plastic	1000	14741.00	39269099
CLSSP002	2mL Pipette Individually Wrapped Paper/Plastic	500	6823.00	39269099
CLSSP003	5mL Pipette Individually Wrapped Paper/Plastic	200	4086.00	39269099
CLSSP004	10mL Pipette Individually Wrapped Paper/Plastic	200	4837.00	39269099
CLSSP005	25mL Pipette Individually Wrapped Paper/Plastic	200	10362.00	39269099
CLSSP006	50mL Pipette Individually Wrapped Paper/Plastic	100	12487.00	39269099
CLSSP007	1mL Pipette Individually Wrapped Plastic	1000	12245.00	39269099
CLSSP008	2mL Pipette Individually Wrapped Plastic	500	7550.00	39269099
CLSSP009	5mL Pipette Individually Wrapped Plastic	200	3653.00	39269099
CLSSP010	10mL Pipette Individually Wrapped Plastic	200	4531.00	39269099
CLSSP011	25mL Pipette Individually Wrappd Plastic	200	9535.00	39269099
CLSSP012	50mL Pipette Individually Wrappd Plastic	100	14551.00	39269099
CLSSP013	1mL Serological Pipette, Bulk Packed	1000	14863.00	39269099
CLSSP014	2mL Serological Pipette, Bulk Packed	500	6614.00	39269099
CLSSP015	5mL Serological Pipette, Bulk Packed	500	11320.00	39269099
CLSSP016	10 mL Serological Pipette, Bulk Packed	500	14389.00	39269099
CLSSP017	25 mL Serological Pipette, Bulk Packed	200	11681.00	39269099
CLSSP018	50mL Serological Pipette, Bulk Packed	100	11878.00	39269099



Cell Strainer

Cell Strainers are ideal for obtaining uniform single cell suspension from various sources. Cell Strainers are made from nylon with 3 different mesh sizes, showing optimal performances in a variety of applications such as stem cell, primary cell preparation.

Cat No.	Description	Pack	Unit Price	HSN
CLSCS50	CELL STRAINER, 40 UM, (Ind)	50	9812.00	39269099
CLSCS70	CELL STRAINER, 70 UM, (Ind)	50	9812.00	39269099
CLSCS100	CELL STRAINER, 100 UM, (Ind)	50	9812.00	39269099



Erlenmeyer Flask

The flasks are individually wrapped in an easy-to-tear, sterile bag, certified to be free of particulate matter, pyrogens, and nucleases (DNase/RNase). These flasks are ideal for all suspension cell cultures including baculovirus cultures, microbial cultures, algae cultures, media preparation, storage, and all related applications. They fit all standard shaking incubator systems and can be placed on any platform.


Cat No.	Description	Pack	Unit Price	HSN
CLSEF125	FLASK, 125ml, 26MM, VC,PC,S,IND	50	9812.00	39269099
CLSEF250	FLASK, 250ml, 31MM, VC,PC,S,IND	50	24024.00	39269099
CLSEF500	FLASK, 500ml, 43MM, VC,PC,S,IND	25	21336.00	39269099
CLSEF1000	FLASK, 1000ml, 43MM, W/VENT,CAO,IND	25	29222.00	39269099



Immunoplate (Elisa Plate)

Microplates for diagnostics and immunological researches. We provide researchers with a variety of microplates dedicated to immunoassays.

Cat No.	Description	Pack	Unit Price	HSN
CLSEP96N	Elisa Plate 96 Well. Flat Bottom, PS, Nonsterile 96 Well	60	7480.00	39269099
CLSEP96S	Elisa Plate 96 Well. Flat Bottom, PS, Sterile 96 Well	60	8794.00	39269099



For more information

+91-120-4968617

For more information

www.compositelabline.com

Vacuum Filtration Units Disposable

The units contain a high quality membrane and combine highest flow-rates and throughput with extremely low protein binding and extractables. They are therefore the product of choice for filtration of tissue culture media, buffers, additives and other aqueous biologicals. They offer improved filtration speed and throughput and are therefore ideally serving the requirement of cost effective filtration.

Low protein binding, Fast flow and high throughput, Easy to use, Engraved graduation.
Ergonomically designed sidewalls and collar can simplify tightening / loosening and adjustments
Designed hose connector can fit multiplicate hose diameters, Non-pyrogenic.
Packaged in easy peel-to-open plastic bag, and receiver bottle cap is individually wrapped.
Each individual unit is lot-numbered for easy identification and tracking.
Sterilized through Gamma Radiation.



Technical Product Information

Material of Const.	Filter Upper Cup Capacity(mL)	Full Unit Overall Height	Filter Diameter (mm)	Working Volume (ml)	Hold-up Volume (mL)	Cap Material	Fitting Outlet (mm)	Maximum Operating Temp. (°C)
GPPS	150	156	50	150	3	PP	45	45
GPPS	250	200	50	250	3	PP	45	45
GPPS	500	245	75	500	3	PP	45	45
GPPS	1000	310	90	1000	3	PP	45	45

Cat Number	Funnel Capacity	Filter Pour Size(µm)	Membrane Material	Membrane Dia. (mm)	Per/case	Unit Price	HSN
CLSPE15075	150 ml	0.22	PES	75	12	7,800.00	84212900
CLSMC15090	150 ml	0.22	MCE	90	12	7,800.00	84212900
CLSPE15050	150 ml	0.45	PES	50	12	7,800.00	84212900
CLSPV25050	250 ml	0.22	PVDF	50	12	8,736.00	84212900
CLSPE25075	250 ml	0.22	PES	75	12	8,736.00	84212900
CLSPV25050	250 ml	0.45	PVDF	50	12	8,736.00	84212900
CLSPE25050	250 ml	0.45	PES	50	12	8,736.00	84212900
CLSPV50050	500 ml	0.22	PVDF	50	12	10,920.00	84212900
CLSPE50075	500 ml	0.22	PES	75	12	10,920.00	84212900
CLSPV50050	500 ml	0.45	PVDF	50	12	10,920.00	84212900
CLSPE50050	500 ml	0.45	PES	50	12	10,920.00	84212900
CLSPE1L75	1000 ml	0.22	PES	75	12	15,600.00	84212900



Bottle Top Filter

Cat No.	Product	Membrane	Pore size (µm)	Sterile	Packing	Price /pack	HSN
CLSPE22150	150ml bottle top filter	PES	0.22	GAMA	24 pcs.	10,649.00	84212900
CLSPE22250	250ml bottle top filter	PES	0.22	GAMA	24 pcs.	13,260.00	84212900
CLSPE22500	500ml bottle top filter	PES	0.22	GAMA	24 pcs.	14,976.00	84212900
CLSPE221L	1000ml bottle top filter	PES	0.22	GAMA	24 pcs.	19,968.00	84212900
CLSPE45150	150ml bottle top filter	PES	0.45	GAMA	24 pcs.	10,649.00	84212900
CLSPE45250	250ml bottle top filter	PES	0.45	GAMA	24 pcs.	13,260.00	84212900
CLSPE45500	500ml bottle top filter	PES	0.45	GAMA	24 pcs.	14,976.00	84212900
CLSPE451L	1000ml bottle top filter	PES	0.45	GAMA	24 pcs.	19,968.00	84212900
CLSMC22150	150ml bottle top filter	MCE	0.22	GAMA	24 pcs.	10,649.00	84212900
CLSMC22250	250ml bottle top filter	MCE	0.22	GAMA	24 pcs.	13,260.00	84212900
CLSMC22500	500ml bottle top filter	MCE	0.22	GAMA	24 pcs.	14,976.00	84212900
CLSMC221L	1000ml bottle top filter	MCE	0.22	GAMA	24 pcs.	19,968.00	84212900
CLSMC45150	150ml bottle top filter	MCE	0.45	GAMA	24 pcs.	13,260.00	84212900
CLSMC45250	250ml bottle top filter	MCE	0.45	GAMA	24 pcs.	13,260.00	84212900
CLSMC45500	500ml bottle top filter	MCE	0.45	GAMA	24 pcs.	14,976.00	84212900
CLSMC451L	1000ml bottle top filter	MCE	0.45	GAMA	24 pcs.	19,968.00	84212900
CLSPV22150	150ml bottle top filter	PVDF	0.22	GAMA	24 pcs.	10,649.00	84212900
CLSPV22250	250ml bottle top filter	PVDF	0.22	GAMA	24 pcs.	13,260.00	84212900
CLSPV22500	500ml bottle top filter	PVDF	0.22	GAMA	24 pcs.	14,976.00	84212900
CLSPV221L	1000ml bottle top filter	PVDF	0.22	GAMA	24 pcs.	19,968.00	84212900
CLSPV45150	150ml bottle top filter	PVDF	0.45	GAMA	24 pcs.	13,260.00	84212900
CLSPV45250	250ml bottle top filter	PVDF	0.45	GAMA	24 pcs.	13,260.00	84212900
CLSPV45500	500ml bottle top filter	PVDF	0.45	GAMA	24 pcs.	14,976.00	84212900
CLSPV451L	1000ml bottle top filter	PVDF	0.45	GAMA	24 pcs.	19,968.00	84212900



COMPOSITE LAB LINE PVT LTD

UGSR-1, Ansal Plaza, Vaishali

Delhi, NCR, Ghaziabad - 201010 (India)

Ph. : 0120-4335636 Mail : sales@compositelablin.com

Visit.: www.compositelabline.com

Reg Office :


A-1, OKKP Complex, Near I.T.Crosing

Faizabad Road, Babujanj

Lucknow - 226020 (India)

Ph. : +919580613355, 0522 - 4068342

Mail : sales@compositelablin.com



Authorised Dealer:

