



Laboratory Essentials

- Professional in
- CLS General Lab Essentials

Composite Lab Line Pvt Ltd

Retort Stand

Non-corrosive, moulded Retort Stands provide excellent replacement for the metal ones. Even the metallic rod has a polypropylene covering making it non- corrosive. The rod measuring 12.5 mm x 750 mm can be screwed into the heavy base giving way to a perfectly stable Retort Stand.
Material : Polystyrene/ Polypropylene

Cat. No.	Size	Pcs/Pack	Price / Pack	HSN Code
CLS-117-01	225 x 150 mm- Side Hold	1	590.00	6903
CLS-117-02	300 x 200 mm-Center Hold	1	975.00	6903
CLS-117-03	300 x 200 mm- Side Hold	1	975.00	6903



Retort Stand

Material : M.S & Blue Powder Coated Rod with L- Key AS PER DIN 12892


Cat. No.	Size	Size (Length x dia)	Pcs/Pack	Price / Pcs	HSN Code
CLS-309-01	150 x 100mm	500mm x 12mm	2	610.00	6903
CLS-309-02	200 x 125mm	600mm x 12mm	2	840.00	6903
CLS-309-03	250 x 150mm	750mm x 12mm	2	1150.00	6903



Retort Stand

Material : Cast Iron & Hammer tone rod with L- Key AS PER DIN 12892

Cat. No.	Size	Size (Length x dia)	Pcs/Pack	Price / Pcs	HSN Code
CLS-310-01	150 x 100 mm	950 x 10 mm	2	510.00	6903
CLS-310-02	175 x 125 mm	600 x 12.5 mm	2	605.00	6903
CLS-310-03	200 x 125mm	600 x 12.5 mm	2	650.00	6903
CLS-310-04	200 x 125mm	600 x 12 mm	2	705.00	6903
CLS-310-05	225 x 150mm	750 x 12 mm	2	1110.00	6903
CLS-310-06	250 x 150mm	600 x 12 mm	2	1310.00	6903
CLS-310-07	250 x 200mm	750 x 12 mm	2	1690.00	6903
CLS-310-08	300 x 175mm	900 x 12 mm	2	1870.00	6903
CLS-310-09	250 x 150mm	750 x 12 mm	2	1410.00	6903



Retort Rings

Material : Mild Steel, Hammer Tone Painted, with Boss Head

Cat. No.	Size	Pcs/Pack	Price / Pack	HSN Code
CLS-313-01	50 mm	4	165.00	6903
CLS-313-02	62.5mm	4	165.00	6903
CLS-313-03	75mm	4	185.00	6903
CLS-313-04	100mm	4	186.00	6903
CLS-313-05	125mm	2	205.00	6903
CLS-313-06	150mm	2	209.00	6903
CLS-313-07	200mm	2	215.00	6903



Retort Ring

Made of steel, chrome plated. With rod from 70 mm to 220 mm length or With boss head to suit rod dia of 14 mm

Cat. No.	L	Ø	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-393-01	160 mm	50mm	Rod	4	222.00	6903
CLS-393-02	200 mm	70mm	Rod	4	236.00	6903
CLS-393-03	220 mm	100mm	Rod	4	270.00	6903
CLS-393-04	220 mm	120mm	Rod	2	302.00	6903
CLS-393-05	220 mm	140mm	Rod	2	280.00	6903
CLS-393-06	70 mm	50mm	Bosshead	4	260.00	6903
CLS-393-07	70 mm	70mm	Bosshead	4	246.00	6903
CLS-393-08	70 mm	100mm	Bosshead	4	306.00	6903
CLS-393-09	70 mm	120mm	Bosshead	2	332.00	6903
CLS-393-10	70 mm	140mm	Bosshead	2	360.00	6903



Retort Clamp

Cross pattern, die-casted, to hold upto 65 mm AS PER DIN 12894

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-324-01	Aluminium, Cork Lined - Cross Pattern,	2	320.00	6903
	to hold upto 50 - 65 mm circles			
CLS-324-04	Zinc, Cork Lined - Cross Pattern,	2	300.00	6903
	to hold upto 50 - 65 mm circles			



Retort Clamp

Material : M.S & Blue Powder Coated Rod with L- Key AS PER DIN 12892

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-326-02	Three Prongs, to hold upto 50 - 65 mm circles	2	344.00	6903



Universal Clamp

Aluminium, Cork Lined, AS PER DIN 12894, Material : Aluminium diecasted, powder coated, chrome plated with MS Steel Rod

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-328-01	to hold upto 50 - 65 mm circles	2	335.00	6903



Speak to our expert
+91-120-4968617

Condenser Clamp - Four Prong

Silicone Sleeved, Material: Aluminium, powder coated with steel rod

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-334-02	to hold upto 125 mm circles	2	830.00	6903



Condenser Clamp - Three Prong

Silicone Sleeved, Material: Aluminium, powder coated with steel rod

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-335-02	to hold upto 125 mm circles	2	835.00	6903



Retort Clamp

Three Prong, EURO DESIGN

Cat. No.	Type	Type	Price / Pack	HSN Code
CLS-390-01	PVC sleeve with rod	2	400.00	6903
CLS-390-02	Silicone sleeve with rod	2	425.00	6903
CLS-390-03	PVC sleeve with boss Head	2	525.00	6903
CLS-390-04	SILICONE sleeve with boss Head	2	560.00	6903



Retort Clamp

Two Prong Standard, EURO DESIGN, Zinc alloy, Chrome plated to hold upto 80 mm

Cat. No.	Type	Pcs/Pack	Price / Pack	HSN Code
CLS-391-01	CORK lined with rod	2	415.00	6903
CLS-391-02	PVC sleeve with rod	2	418.00	6903
CLS-391-03	SILICONE sleeve with rod	2	418.00	6903
CLS-391-04	CORK lined with boss head	2	520.00	6903
CLS-391-05	PVC sleeve with boss head	2	520.00	6903
CLS-391-06	SILICONE sleeve with boss head	2	545.00	6903



BOSS HEAD

Material : Cast Iron, painted, for use retort Stand and Condenser Clamp AS PER DIN 12895

Cat. No.	Size (Dia x Length)	Pcs/Pack	Price / Pack	HSN Code
CLS-337-01	hold upto 25mm heavy, with round screw	4	360.00	6903
CLS-337-02	hold upto 25mm heavy, with thumb screw	4	360.00	6903
CLS-337-03	hold upto 25mm heavy, with rod screw	4	360.00	6903
CLS-338-01	hold upto 20mm heavy with round screw	4	205.00	6903
CLS-338-03	hold upto 20mm heavy with rod screw	4	205.00	6903



BOSS HEAD

Material : Powder coated, AS PER DIN 12895

Cat. No.	Description	Pcs/Pack	Price / Pack	HSN Code
CLS-341-01	Aluminium, Round Screw, hold upto 16 mm AS PER DIN 12895	4	160.00	6903
CLS-341-02	Aluminium, Thumb Screw, hold upto 16 mm AS PER DIN 12895	4	156.00	6903
CLS-341-03	Aluminium, Rod Screw, hold upto 16 mm AS PER DIN 12895	4	158.00	6903
CLS-341-04	Zinc, Round Screw, hold upto 16 mm AS PER DIN 12895	4	136.00	6903
CLS-341-05	Zinc, Thumb Screw, hold upto 16 mm AS PER DIN 12895	4	152.00	6903
CLS-341-06	Zinc, Rod Screw, hold upto 16 mm AS PER DIN 12895	4	158.00	6903


Tripod Stand

Material : Cast iron base, Hammer Tone, Finished Round

Cat. No.	Size	Top	Pcs/Pack	Price / Pcs	HSN Code
CLS-344-01	100 x 150 mm	with round top	1	210.00	6903
CLS-344-02	100 x 175 mm	with round top	1	220.00	6903
CLS-344-03	100 x 200 mm	with round top	1	230.00	6903
CLS-344-04	125 x 150 mm	with round top	1	230.00	6903
CLS-344-05	125 x 175 mm	with round top	1	240.00	6903
CLS-344-06	125 x 200 mm	with round top	1	240.00	6903
CLS-344-07	100 x 150 mm	with triangle top	1	205.00	6903
CLS-344-08	100 x 175 mm	with triangle top	1	210.00	6903
CLS-344-09	100 x 200 mm	with triangle top	1	220.00	6903
CLS-344-10	125 x 150 mm	with triangle top	1	240.00	6903
CLS-344-11	125 x 175 mm	with triangle top	1	235.00	6903
CLS-344-12	125 x 200 mm	with triangle top	1	245.00	6903


Burette Clamp

 Made of Aluminium, die casted, blue painted, black rubber sleeves for firm gripping of buret tes.
 With adjustable clamps to hold upto 30 mm, spring type, thumb screw key for retort rods upto 16 mm dia.

Cat. No.	Type	Pcs/Pack	Price / Pcs	HSN Code
CLS-394-01	Single	12	570.00	6903
CLS-394-02	Double	6	510.00	6903


Burette Clamp, Mild Steel, Single

Material : Mild Steel, Chrome Plated, Fisher Type with Boss Head, Single.

Cat. No.	Type	Max. Temp.	Pcs/Pack	Price / Pcs	HSN Code
CLS-318-01	PVC Sleeve	80°C	12	250.00	6903
CLS-318-02	Silicon Sleeve	280°C	12	255.00	6903



Burette Clamp, Mild Steel, Single

Material : Mild Steel, Black Powder Coated, Fisher Type with Boss Head, Single.

Cat. No.	Type	Max. Temp.	Pcs/Pack	Price / Pcs	HSN Code
CLS-319-01	PVC Sleeve	80°C	12	235.00	6903
CLS-319-02	Silicon Sleeve	280°C	12	230.00	6903



Burette Clamp, Mild Steel, Double

Material : Mild Steel, Chrome Plated, Fisher Type with Boss Head, Double

Cat. No.	Type	Max. Temp.	Pcs/Pack	Price / Pcs	HSN Code
CLS-320-01	PVC Sleeve	80°C	6	380.00	6903
CLS-320-02	Silicon Sleeve	280°C	6	375.00	6903



Burette Clamp, Mild Steel

Material : Mild Steel, Powder Coated, Fisher Type with Boss Head, Double


Cat. No.	Type	Max. Temp.	Pcs/Pack	Price / Pcs	HSN Code
CLS-321-01	PVC Sleeve	80°C	6	365.00	6903
CLS-321-02	Silicon Sleeve	280°C	6	360.00	6903



Wire Gauge

Ceramic Center, Maximum Temperature Upto 800°C

Cat. No.	Type	Pcs/Pack	Price / Pcs	HSN Code
CLS-524-01	100 x 100 mm	10	55.00	6903
CLS-524-02	125 x 125 mm	10	60.00	6903
CLS-524-03	150 x 150 mm	10	64.00	6903
CLS-524-04	200 x 200 mm	10	98.00	6903



Bunsen Burner

With Stop Cock L.P.G.

Cat. No.	Description	Pcs/Pack	Price / Pack	HSN Code
CLS-347-01	Casted Iron painted AS PER DIN 30665	1	505.00	6903
CLS-347-02	Casted Iron painted, easy flame regulator double walled	1	545.00	6903
Brass pipe to give a concentrated higher flame, Eco-Friendly				
AS PER DIN 30665				



Meaker Burner

Material : Casted Iron painted AS PER DIN 30665

Cat. No.	Description	Pcs/Pack	Price / Pack	HSN Code
CLS-348-01	With Stop Cock	1	910.00	6903



Meaker Burner

Material : Casted Iron painted AS PER DIN 30665

Cat. No.	Description	Pcs/Pack	Price / Pack	HSN Code
CLS-349-01	With Stop Cock L.P.G	1	1125.00	6903
CLS-349-02	ECONOMIC With Stop Cock L.P.G, Bunsen	1	1075.00	6903
Burner type with cap, for school purpose				


Laboratory Jack

All sizes are in millimeters

Cat No.	Type	Size (mm)	Pcs/Pack	Price / Pack	HSN Code
CLS-370-09	Stainless Steel	120 x 140	1	4210.00	6903
CLS-370-10	Stainless Steel	200 x 200	1	6810.00	6903
CLS-370-11	Stainless Steel	250 x 250	1	9210.00	6903


Silica Crucible

With Lid Upto 1000°C

Cat. No.	Capacity	Dia (mm)	Height (mm)	Pcs/Pack	Price / Pcs	HSN Code
CLS-518-01	15	41	25	6	690.00	6903
CLS-518-02	25	47	28	6	770.00	6903
CLS-518-03	50	57	37	6	1350.00	6903
CLS-518-04	80	67	45	4	2250.00	6903
CLS-518-05	100	71	50	4	2800.00	6903
CLS-518-06	150	82	51	4	3250.00	6903


Silica Crucible

Without Lid Upto 1000°C

Cat. No.	Capacity	Dia (mm)	Height (mm)	Pcs/Pack	Price / Pcs	HSN Code
CLS-519-01	15	41	25	6	340.00	6903
CLS-519-02	25	47	28	6	370.00	6903
CLS-519-03	50	57	37	6	740.00	6903
CLS-519-04	80	67	45	4	1260.00	6903
CLS-519-05	100	71	50	4	1600.00	6903
CLS-519-06	150	82	51	4	1950.00	6903


Quartz Crucible

Made of Quartz Gel & Powder, Tolerance upto 2000°C

Cat. No.	Size	Type	Pcs/Pack	Price / Pcs	HSN Code
CLS-521-01	15 ML	Without Lid	4	680.00	7017
CLS-521-02	25 ML	Without Lid	4	760.00	7017
CLS-521-03	50 ML	Without Lid	2	1190.00	7017
CLS-521-04	100 ML	Without Lid	2	2290.00	7017
CLS-521-05	15 ML	With Lid	4	1170.00	7017
CLS-521-06	25 ML	With Lid	4	1310.00	7017
CLS-521-07	50 ML	With Lid	2	1870.00	7017
CLS-521-08	100 ML	With Lid	2	3650.00	7017



Crucible Porcelain (Euro Design)

Made of Ceramic With Lid Maximum temperature upto 1150°C

Cat. No.	Capacity	Dia	Height	Pcs/Pack	Price / Pcs	HSN Code
CLS-522-13	2 ml	18 mm	12 mm	6	105.00	6903
CLS-522-14	8 ml	32 mm	20 mm	6	115.00	6903
CLS-522-15	17 ml	41 mm	26 mm	4	140.00	6903
CLS-522-16	30 ml	50 mm	31 mm	4	165.00	6903
CLS-522-17	50 ml	61 mm	37 mm	4	200.00	6903
CLS-522-18	100 ml	76 mm	46 mm	2	305.00	6903
CLS-522-19	150 ml	84 mm	52 mm	2	350.00	6903
CLS-522-20	250 ml	102 mm	60 mm	2	480.00	6903


Crucible Porcelain, Tall form (Euro Design)

Made of Ceramic With Lid Tall Form Maximum temperature upto 1150°C

Cat. No.	Capacity	Dia	Height	Pcs/Pack	Price / Pcs	HSN Code
CLS-523-11	1.5 ml	15 mm	15 mm	6	105.00	6903
CLS-523-12	6 ml	24 mm	20 mm	6	115.00	6903
CLS-523-13	10 ml	31 mm	26 mm	4	120.00	6903
CLS-523-14	15 ml	35 mm	29 mm	4	140.00	6903
CLS-523-15	30 ml	43 mm	37 mm	4	165.00	6903
CLS-523-16	40 ml	50 mm	40 mm	2	185.00	6903
CLS-523-17	50 ml	53 mm	43 mm	2	220.00	6903
CLS-523-18	100 ml	65 mm	54 mm	2	300.00	6903
CLS-523-19	250 ml	88 mm	72 mm	2	450.00	6903


Crucible Tong, SS


Stainless Steel. Grade 304

Cat. No.	Size (mm)	Pcs/Pack	Price / Pcs	HSN Code
CLS-355-01	150	4	140.00	6903
CLS-355-02	200	4	180.00	6903
CLS-355-03	250	4	270.00	6903
CLS-355-04	300	4	330.00	6903
CLS-355-05	375	4	450.00	6903
CLS-355-06	450	4	530.00	6903
CLS-355-07	600	4	1150.00	6903


Crucible Tong, Mild Steel

Made of mild steel, black paint/chrome plated.

Cat. No.	Size (mm)	Type	Pcs/Pack	Price / Pcs	HSN Code
CLS-367-01	150	black painted	4	120.00	6903
CLS-367-02	200	black painted	4	130.00	6903
CLS-367-03	300	black painted	4	260.00	6903
CLS-367-04	150	chrome plated	4	130.00	6903
CLS-367-05	200	chrome plated	4	140.00	6903
CLS-367-06	300	chrome plated	4	245.00	6903



Evaporating Dish, EURO DESIGN

Made of Porcelain, maximum temperature upto 1000°C

Cat No.	Capacity	OD x H	Pcs/Pack	Price / Pcs	HSN Code
CLS-528-07	35 ML	60 x 24 mm	6	110.00	6903
CLS-528-08	70 ML	75 x 30 mm	6	125.00	6903
CLS-528-09	80 ML	80 x 30 mm	6	130.00	6903
CLS-528-10	120 ML	94 x 42 mm	6	140.00	6903
CLS-528-11	150 ML	100 x 42 mm	6	150.00	6903
CLS-528-12	250 ML	115 x 45 mm	4	160.00	6903
CLS-528-13	385 ML	145 x 48 mm	4	390.00	6903
CLS-528-14	525 ML	162 x 51 mm	2	590.00	6903
CLS-528-15	765 ML	185 x 54 mm	2	1220.00	6903
CLS-528-16	1285 ML	215 x 75 mm	2	1540.00	6903
CLS-528-17	2100 ML	265 x 80 mm	1	1750.00	6903
CLS-528-18	3250 ML	305 x 95 mm	1	3590.00	6903


Pestle & Mortar, EURO DESIGN

Made of Porcelain, maximum temperature up to 1000°C

Cat. No.	Size (ID x H)	Pcs/Pack	Price / Pcs	HSN Code
CLS-529-08	60 x 44 mm	6	190.00	6903
CLS-529-09	75 x 48 mm	4	210.00	6903
CLS-529-10	85 x 55 mm	4	250.00	6903
CLS-529-11	100 x 70 mm	4	350.00	6903
CLS-529-12	117 x 75 mm	4	475.00	6903
CLS-529-13	125 x 88 mm	2	430.00	6903
CLS-529-15	150 x 100 mm	1	680.00	6903
CLS-529-16	175 x 115 mm	1	1410.00	6903
CLS-529-17	200 x 140 mm	1	1510.00	6903
CLS-529-18	250 x 175 mm	1	5950.00	6903
CLS-529-19	300 x 205 mm	1	7150.00	6903
CLS-529-20	350 mm ID	1	14750.00	6903
CLS-529-21	400 mm ID	1	21950.00	6903


Rubber Tubing

Material : Natural Rubber, Temperature range -30 to 70°C Red Superior 10 Meter Coil

Cat. No.	Bore	Wall Thickness	Pcs/Pack	Price / Pack	HSN Code
CLS-406-01	4.0 mm	1.0 mm	10 Mtr	240.00	6903
CLS-406-02	5.0 mm	1.5 mm	10 Mtr	320.00	6903
CLS-406-03	6.0 mm	1.5 mm	10 Mtr	380.00	6903
CLS-406-04	7.0 mm	1.5 mm	10 Mtr	460.00	6903
CLS-406-05	8.0 mm	1.5 mm	10 Mtr	530.00	6903
CLS-406-06	9.0 mm	2.0 mm	10 Mtr	650.00	6903
CLS-406-07	10.0 mm	2.0 mm	10 Mtr	710.00	6903
CLS-406-08	12.0 mm	2.5 mm	10 Mtr	920.00	6903


Dissecting Set

Cat. No.	Description	Pcs/Pack	Price / Pack	HSN Code
CLS-536-01	Made of Stainless Steel, Instruments are Dressing Scissor, Cuticle Scissor, Dressing Forceps, Dissecting Forceps, Flapper, Seeker, Arrow Needle, Needle, Blow Pipe, B.P. Handle with Scalpel Blade, Dropper, Chain Hook, Scalpels, Magnifier Lense.	1 set	790.00	6903



Forceps

Stainless Steel straight fine points, 410 SS Grade

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-531-01	100 mm	10	44.00	6903
CLS-531-02	125 mm	10	60.00	6903
CLS-531-03	150 mm	10	72.00	6903
CLS-531-04	200 mm	10	96.00	6903
CLS-531-05	250 mm	10	105.00	6903
CLS-531-06	300 mm	10	130.00	6903



Forceps

Curved Fine Points, 410 SS Grade

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-532-01	100 mm	10	48.00	6903
CLS-532-02	125mm	10	72.00	6903
CLS-532-03	150 mm	10	85.00	6903
CLS-532-04	200 mm	10	95.00	6903
CLS-532-05	250 mm	10	115.00	6903
CLS-532-06	300 mm	10	140.00	6903



Forceps

Stainless Steel, Straight Blunt Points, 410 SS Grade

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-533-01	100 mm	10	44.00	6903
CLS-533-02	125mm	10	60.00	6903
CLS-533-03	150 mm	10	72.00	6903
CLS-533-04	200 mm	10	98.00	6903
CLS-533-05	250 mm	10	106.00	6903
CLS-533-06	300 mm	10	130.00	6903



Forceps

Stainless Steel, Straight Blunt Points, 410 SS Grade

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-534-01	100 mm	10	48.00	6903
CLS-534-02	125mm	10	65.00	6903
CLS-534-03	150 mm	10	80.00	6903
CLS-534-04	200 mm	10	90.00	6903
CLS-534-05	250 mm	10	110.00	6903
CLS-534-06	300 mm	10	130.00	6903



Spatula

Material : Stainless Steel, Grade 202


Cat. No.	Size	Type	Pcs/Pack	Price / Pcs	HSN Code
CLS-374-01	125 mm	one end spoon and one end flat	10	35.00	6903
CLS-374-02	150 mm	one end spoon and one end flat	10	45.00	6903
CLS-374-03	200 mm	one end spoon and one end flat	10	62.00	6903
CLS-374-04	250 mm	one end spoon and one end flat	10	75.00	6903
CLS-374-05	300 mm	one end spoon and one end flat	10	95.00	6903
CLS-374-06	125 mm	both end flat	10	35.00	6903
CLS-374-07	150 mm	both end flat	10	44.00	6903
CLS-374-08	200 mm	both end flat	10	62.00	6903
CLS-374-09	250 mm	both end flat	10	70.00	6903



Rubber Corks

Autoclavable with temperature range -30 to 70°C, Shore Hardness 45

Cat. No.	Type	Size (Bottom Dia x Top Dia x Height)	Pcs/Pack	Price / Pack	HSN Code
CLS-410-01		9 x 6 x 16 mm	10	28.00	6903
CLS-410-02		13 x 9 x 15.5 mm	10	30.00	6903
CLS-410-03		14 x 11 x 18.5 mm	10	36.00	6903
CLS-410-04		16.5 x 13 x 20 mm	10	42.00	6903
CLS-410-05		19 x 14 x 21 mm	10	51.00	6903
CLS-410-06		21 x 15 x 22 mm	10	56.00	6903
CLS-410-07		22.5 x 17.5 x 23 mm	10	70.00	6903
CLS-410-08		24.5 x 20 x 24 mm	10	94.00	6903
CLS-410-09		26 x 20 x 27.5 mm	10	102.00	6903
CLS-410-10		26.5 x 21 x 27.5 mm	10	125.00	6903
CLS-410-11		28 x 21.5 x 31 mm	10	155.00	6903
CLS-410-12		30 x 26 x 33 mm	10	180.00	6903
CLS-410-13		33 x 29 x 32 mm	10	210.00	6903
CLS-410-14		34.5 x 30 x 34.5 mm	10	230.00	6903
CLS-410-15		38.5 x 31 x 37 mm	10	300.00	6903
CLS-410-16		41.5 x 32 x 42.5 mm	10	340.00	6903
CLS-410-17		42 x 38 x 40.5 mm	10	438.00	6903
CLS-410-18		45 x 35 x 36.5 mm	10	440.00	6903
CLS-410-19		49.5 x 40 x 40.5 mm	10	530.00	6903
CLS-410-20		53.5 x 48 x 33.5 mm	10	590.00	6903
CLS-410-21		58 x 44 x 45 mm	10	690.00	6903
CLS-410-22		60 x 51 x 56.5 mm	10	970.00	6903
CLS-410-23	EURO DESIGN	3.5 x 6.5 x 15 mm	10	30.00	6903
CLS-410-24	EURO DESIGN	5 x 9 x 20 mm	10	30.00	6903
CLS-410-25	EURO DESIGN	8 x 12 x 20 mm	10	32.00	6903
CLS-410-26	EURO DESIGN	10.5 x 14.5 x 20 mm	10	36.00	6903
CLS-410-27	EURO DESIGN	12.5 x 16.5 x 20 mm	10	46.00	6903
CLS-410-28	EURO DESIGN	14 x 18 x 20 mm	10	66.00	6903
CLS-410-29	EURO DESIGN	17 x 22 x 25 mm	10	72.00	6903
CLS-410-30	EURO DESIGN	18 x 24 x 30 mm	10	84.00	6903
CLS-410-31	EURO DESIGN	21 x 27 x 30 mm	10	140.00	6903
CLS-410-32	EURO DESIGN	23 x 29 x 30 mm	10	150.00	6903
CLS-410-33	EURO DESIGN	26 x 32 x 30 mm	10	166.00	6903
CLS-410-34	EURO DESIGN	29 x 35 x 30 mm	10	192.00	6903
CLS-410-35	EURO DESIGN	31 x 38 x 35 mm	10	292.00	6903
CLS-410-36	EURO DESIGN	36 x 44 x 40 mm	10	390.00	6903
CLS-410-37	EURO DESIGN	41 x 49 x 40 mm	10	480.00	6903
CLS-410-38	EURO DESIGN	47 x 55 x 40 mm	10	540.00	6903
CLS-410-39	EURO DESIGN	50.5 x 59.5 x 45 mm	10	630.00	6903
CLS-410-40	EURO DESIGN	56 x 65 x 45 mm	10	1120.00	6903
CLS-410-41	EURO DESIGN	60 x 70 x 50 mm	10	1410.00	6903
CLS-410-42	EURO DESIGN	64.5 x 75.5 x 55 mm	10	1610.00	6903
CLS-410-43	EURO DESIGN	71 x 83 x 60 mm	10	1720.00	6903
CLS-410-44	EURO DESIGN	79 x 92 x 65 mm	10	1780.00	6903
CLS-410-45	EURO DESIGN	87 x 100 x 65 mm	10	2030.00	6903
CLS-410-46	EURO DESIGN	94 x 107 x 65 mm	10	2250.00	6903



Silicone Cork

Autoclavable with temperature range -50 to 250°C, Shore Hardness 45-50

Cat. No.	Type	Size (Bottom Dia x Top Dia x Height)	Pcs/Pack	Price / Pack	HSN Code
CLS-405-01	TRANSLUCENT	9 x 11.5 x 20	10	240.00	6903
CLS-405-02	TRANSLUCENT	11 x 14 x 24	10	380.00	6903
CLS-405-03	TRANSLUCENT	13 x 16 x 24	10	450.00	6903
CLS-405-04	TRANSLUCENT	15 x 18 x 24	10	550.00	6903
CLS-405-05	TRANSLUCENT	17 x 20 x 26	10	670.00	6903
CLS-405-06	TRANSLUCENT	18 x 21 x 26	10	740.00	6903
CLS-405-07	TRANSLUCENT	19 x 22.5 x 28	10	850.00	6903
CLS-405-08	TRANSLUCENT	23 x 25.5 x 28	10	1090.00	6903
CLS-405-09	TRANSLUCENT	25 x 28.5 x 28	10	1310.00	6903
CLS-405-10	TRANSLUCENT	27 x 31 x 32	10	1580.00	6903
CLS-405-11	TRANSLUCENT	29 x 33 x 32	10	2020.00	6903
CLS-405-12	TRANSLUCENT	31 x 36 x 35	10	2400.00	6903
CLS-405-13	TRANSLUCENT	33 x 38.5 x 38.5	10	3100.00	6903
CLS-405-14	EURO DESIGN	3.5 x 6.5 x 15	10	110.00	6903
CLS-405-15	EURO DESIGN	5 x 9 x 20	10	140.00	6903
CLS-405-16	EURO DESIGN	8 x 12 x 20	10	170.00	6903
CLS-405-17	EURO DESIGN	10.5 x 14.5 x 20	10	260.00	6903
CLS-405-18	EURO DESIGN	12.5 x 16.5 x 20	10	320.00	6903
CLS-405-19	EURO DESIGN	14 x 18 x 20	10	340.00	6903
CLS-405-20	EURO DESIGN	17 x 22 x 25	10	590.00	6903
CLS-405-21	EURO DESIGN	18 x 24 x 30	10	640.00	6903
CLS-405-22	EURO DESIGN	21 x 27 x 30	10	760.00	6903
CLS-405-23	EURO DESIGN	23 x 29 x 30	10	940.00	6903
CLS-405-24	EURO DESIGN	26 x 32 x 30	10	1140.00	6903
CLS-405-25	EURO DESIGN	29 x 35 x 30	10	1770.00	6903
CLS-405-26	EURO DESIGN	31 x 38 x 35	10	2130.00	6903
CLS-405-27	EURO DESIGN	36 x 44 x 40	10	2840.00	6903
CLS-405-28	EURO DESIGN	41 x 49 x 40	10	3260.00	6903
CLS-405-29	EURO DESIGN	47 x 55 x 40	10	2910.00	6903
CLS-405-30	EURO DESIGN	50.5 x 59.5 x 45	10	4290.00	6903
CLS-405-31	EURO DESIGN	56 x 65 x 45	10	4990.00	6903
CLS-405-32	EURO DESIGN	60 x 70 x 50	10	7810.00	6903
CLS-405-33	EURO DESIGN	64.5 x 75.5 x 55	10	9790.00	6903
CLS-405-34	EURO DESIGN	71 x 83 x 60	10	12510.00	6903
CLS-405-35	EURO DESIGN	79 x 92 x 65	10	16200.00	6903
CLS-405-36	EURO DESIGN	87 x 100 x 65	10	21600.00	6903
CLS-405-37	EURO DESIGN	94 x 107 x 65	10	23390.00	6903



Spatula Knife

Wooden Handle

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-541-01	100 mm	6	80.00	6903
CLS-541-03	150 mm	6	90.00	6903
CLS-541-04	200 mm	6	100.00	6903
CLS-541-05	250 mm	6	140.00	6903
CLS-541-06	300 mm	6	160.00	6903

Double Spatula

Material : Stainless steel, Grade 202, one end Chattaway, one end flat

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-375-01	155 mm	10	50.00	6903
CLS-375-02	210 mm	10	55.00	6903
CLS-375-03	250 mm	10	70.00	6903



Double Spatula

Material : Stainless steel, one end tapered, one end flat, flexible type, Grade 202

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-376-01	155 mm	10	60.00	6903
CLS-376-02	260 mm	10	95.00	6903
CLS-376-03	305 mm	10	100.00	6903




Micro Spatula

Material : Stainless steel, one end spoon and one side spatula, Grade 202

Cat. No.	Size	Pcs/Pack	Price / Pcs	HSN Code
CLS-377-01	155 mm	10	56.00	6903



For more information



+91-120-4968617

For more information



www.compositelabline.com

COMPOSITE LAB LINE PVT LTD

UGSR-1, Ansal Plaza, Vaishali

Delhi, NCR, Ghaziabad - 201010 (India)

Ph. : 0120-4335636 Mail : sales@compositelablin.com

Visit.: www.compositelabline.com

Reg Office :


A-1, OKKP Complex, Near I.T.Crosing

Faizabad Road, Babujanj

Lucknow - 226020 (India)

Ph. : +919580613355, 0522 - 4068342

Mail : sales@compositelablin.com



Authorised Dealer:

