




Chromatography Accessories



- Professional in
- CLS HPLC Septa & Vials

Composite Lab Line Pvt Ltd

Analytical Vials and Septa Selections Guide Along with Technical Information



Capacity ml	3ml	6ml	8ml	10ml	15ml	25ml	30ml	16ml	20ml
Dimension H X D mm	41 x 11.6 mm	41 x 22 mm	61 x 16.6mm	50 x 22mm	48 X 26mm	65 X 26mm	86 X 30mm	71 x 20.6mm	86 x 22.7
Mouth Type	Crimp Neck	Crimp Neck	Screw Neck	Flip Neck	Flip Neck	Flip Neck	Flip Neck	Screw Neck	Screw Neck
Mouth Diameter mm	11 mm	20 mm	15mm	18mm	22mm	22mm	30mm	18mm	20mm
Clear	CLS-CB03C01100C000	CLS-CB06CF2000C000	CLS-CB08C01500S000	CLS-CB10C01800N000	CLS-CB15C02200N000	CLS-CB25C02200N000	CLS-CV30C130PCN000	CLS-CB16C01800S000	CLS-CB20C02000C000
Amber	CLS-CB03A011C000	CLS-CB06AF2000C00	CLS-CB08A01500S000	-	-	-	CLS-CV30A130PCN000	CLS-CB16A01800S000	CLS-CB20A02000C000
Suggested Cap	CLS-BC000011ASC011	CLS-BC000020ASC051	CLS-BC000015PKS011	CLS-BC000018PCN000	CLS-BC000022PCN000	CLS-BC000022PCN000	Supplied with cpas	CLS-BC000018PKS011	CLS-BC000020PWS011
Application	Storage	Storage	Storage	Storage/ Sampling	Storage/ Sampling	Storage/ Sampling	Process sampling (IPOA)	Storage	Storage

Vials made of Type Ist Borosilicate glass, with different volumes in clear and amber. These are non-stock items and may require a minimum order quantity. Please contact for pricing, closure and ordering information.



Speak to our expert
+91-120-4968617

9mm Screw Vials

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CB02C00900S000	9mm 12X32 Clear Screw Vials	10,150.00	70179090
CLS-CB02C10900S000	9mm 12X32 Clear Screw Vials with Write on Patch	12,140.00	70179090
CLS-CB02A00900S000	9mm 12X32 Amber Screw Vials	13,410.00	70179090
CLS-CB02A10900S000	9mm 12X32 Amber Screw Vials with Write on Patch	16,330.00	70179090

9 mm Screw caps with septa

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-BC000009PBS011	9mm Screw Caps with PTFE-SILICONE SEPTA	11,780.00	70179090
CLS-BC000009PBS012	9mm Screw Caps with PRE-SLIT PTFE-SILICONE SEPTA	13,280.00	70179090
CLS-BC000009PBS111	9mm Screw Caps with BONDED PTFE-SILICONE SEPTA	14,340.00	70179090
CLS-BC000009PBS112	9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	15,100.00	70179090
CLS-BC000009PBS113	9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	16,490.00	70179090
CLS-BC000009PBS115	9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	16,490.00	70179090
CLS-BC000009PKS001	9mm Screw Closed Top Screw Caps with PTFE Lining	13,590.00	70179090

9 mm Combipack of Clear Screw Vials with matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02C009PBS011	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with PTFE-SILICONE SEPTA	21,930.00	70179090
CLS-CV02C109PBS011	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with PTFE-SILICONE SEPTA	23,920.00	70179090
CLS-CV02C009PBS012	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with PRE-SLIT PTFE-SILICONE SEPTA	23,430.00	70179090
CLS-CV02C109PBS012	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with PRE-SLIT PTFE-SILICONE SEPTA	25,420.00	70179090
CLS-CV02C009PBS111	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with BONDED PTFE-SILICONE SEPTA	24,490.00	70179090
CLS-CV02C109PBS111	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PTFE-SILICONE SEPTA	26,480.00	70179090
CLS-CV02C009PBS112	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	25,250.00	70179090
CLS-CV02C109PBS112	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	27,240.00	70179090
CLS-CV02C009PBS113	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	26,630.00	70179090
CLS-CV02C109PBS113	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	28,620.00	70179090
CLS-CV02C009PBS115	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	28,620.00	70179090

CLS-CV02C109PBS115	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	28,620.00	70179090
CLS-CV02C009PKS001	Combi pack of 9mm 12X32 Clear Screw Vials and 9mm Screw Closed Top Screw Caps with PTFE Lining	23,740.00	70179090
CLS-CV02C109PKS001	Combi pack of 9mm 12X32 Clear Screw Vials with Write on Patch and 9mm Screw Closed Top Screw Caps with PTFE Lining	25,730.00	70179090

Combi pack of 9 mm Amber Screw Vials with matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02A009PBS011	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with PTFE-SILICONE SEPTA	25,190.00	70179090
CLS-CV02A109PBS011	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with PTFE-SILICONE SEPTA	28,110.00	70179090
CLS-CV02A009PBS012	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with PRE-SLIT PTFE-SILICONE SEPTA	26,690.00	70179090
CLS-CV02A109PBS012	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with PRE-SLIT PTFE-SILICONE SEPTA	29,610.00	70179090
CLS-CV02A009PBS111	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with BONDED PTFE-SILICONE SEPTA	27,750.00	70179090
CLS-CV02A109PBS111	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PTFE-SILICONE SEPTA	30,670.00	70179090
CLS-CV02A009PBS112	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	28,510.00	70179090
CLS-CV02A109PBS112	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	31,430.00	70179090
CLS-CV02A009PBS113	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	29,890.00	70179090
CLS-CV02A109PBS113	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	32,820.00	70179090
CLS-CV02A009PBS115	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	29,890.00	70179090
CLS-CV02A109PBS115	Combi pack of 9mm 12X32 Amber Screw Vials with Write on Patch and 9mm Screw Caps with BONDED PRE-SLIT PTFE-SILICONE SEPTA	32,820.00	70179090
CLS-CV02A009PKS001	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Closed Top Screw Caps with PTFE Lining	27,000.00	70179090
CLS-CV02A109PKS001	Combi pack of 9mm 12X32 Amber Screw Vials and 9mm Screw Caps with BONDED PTFE Lining	29,920.00	70179090

11 mm Snap Vials

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CB02C01100N000	11mm Clear Glass Snap Vials 12 X 32	8,040.00	70179090
CLS-CB02C11100N000	11mm Clear Glass Snap Vials 12 X 32 with Write On Patch	8,520.00	70179090
CLS-CB02A01100N000	11mm Amber Glass Snap Vials 12 X 32	10,310.00	70179090



Speak to our expert
+91-120-4968617

CLS-CB02A11100N000	11mm Amber Glass Snap Vials 12 X 32 with Write On Patch	10,560.00	70179090
--------------------	---	-----------	----------

11 mm Snap Caps with Septa

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-BC000011PCN011	11mm Snap Caps with PTFE-SILICONE SEPTA	11,930.00	70179090
CLS-BC000011PCN012	11mm Snap Caps with PRE-SLIT PTFE-SILICONE SEPTA	13,580.00	70179090

11 mm combipacks of clear snap vials with matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02C011PCN011	Combi pack of 11mm Clear Glass Snap Vials 12 X 32 and 11mm Snap Caps with PTFE-SILICONE SEPTA	19,980.00	70179090
CLS-CV02C011PCN012	Combi pack of 11mm Clear Glass Snap Vials 12 X 32 and 11mm Snap Caps with PRE-SLIT PTFE-SILICONE SEPTA	21,630.00	70179090
CLS-CV02C111PCN011	Combi pack of 11mm Clear Glass Snap Vials 12 X 32 with Write On Patch and 11mm Snap Caps with PTFE-SILICONE SEPTA	20,450.00	70179090
CLS-CV02C111PCN012	Combi pack of 11mm Clear Glass Snap Vials 12 X 32 with Write On Patch and 11mm Snap Caps with PRE-SLIT PTFE-SILICONE SEPTA	22,100.00	70179090

Combi pack of 11 mm Snap Amber Vials with Matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02A011PCN011	Combi pack of 11mm Amber Glass Snap Vials 12 X 32 and 11mm Snap Caps with PTFE-SILICONE SEPTA	22,250.00	70179090
CLS-CV02A011PCN012	Combi pack of 11mm Amber Glass Snap Vials 12 X 32 and 11mm Snap Caps with PRE-SLIT SILICONE SEPTA	2,390.00	70179090
CLS-CV02A111PCN011	Combi pack of 11mm Amber Glass Snap Vials 12 X 32 with Write On Patch and 11mm Snap Caps with PTFE-SILICONE SEPTA	22,490.00	70179090
CLS-CV02A111PCN012	Combi pack of 11mm Amber Glass Snap Vials 12 X 32 Caps with PRE-SLIT PTFE-SILICONE SEPTA	24,140.00	70179090

11 mm Crimp Vials

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CB02C01100C000	11mm Clear Glass Crimp Vials 12 X 32	8,920.00	70179090
CLS-CB02C11100C000	11mm Clear Glass Crimp Vials 12 X 32 with Write On Patch	9,640.00	70179090
CLS-CB02A01100C000	11mm Amber Glass Crimp Vials 12 X 32	10,450.00	70179090
CLS-CB02A11100C000	11mm Amber Glass Crimp Vials 12 X 32 with Write On Patch	12,830.00	70179090

11 mm Crimp Caps with Septa

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-BC000011ASC011	11mm Aluminum Crimp Caps with PTFE-SILICONE SEPTA	10,870.00	70179090
CLS-BC000011ASC012	11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	12,380.00	70179090



Speak to our expert
+91-120-4968617

CLS-BC000011ASC021	11mm Aluminum Crimp Caps with PTFE-ORANGE SILICONE SEPTA	8,450.00	70179090
CLS-BC000011ASC022	11mm Aluminum Crimp Caps with PRE-SLIT PTFE-ORANGE SILICONE SEPTA	9,960.00	70179090

11 mm Combipacks of Clear Crimp Vials with Matching Caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02C011ASC011	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	19,780.00	70179090
CLS-CV02C011ASC012	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	21,300.00	70179090
CLS-CV02C011ASC021	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	17,370.00	70179090
CLS-CV02C011ASC022	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	18,880.00	70179090
CLS-CV02C111ASC011	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PTFE-SILICONE SEPTA	20,500.00	70179090
CLS-CV02C111ASC012	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	22,020.00	70179090
CLS-CV02C111ASC021	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PTFE-ORANGE SILICONE SEPTA	18,090.00	70179090
CLS-CV02C111ASC022	Combi pack of 11mm Clear Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-ORANGE SILICONE SEPTA	19,600.00	70179090

Combipack of Amber Crimp Vials with Matching Caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV02A011ASC011	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PTFE-SILICONE SEPTA	21,320.00	70179090
CLS-CV02A011ASC012	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	22,830.00	70179090
CLS-CV02A011ASC021	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PTFE-ORANGE SILICONE SEPTA	18,900.00	70179090
CLS-CV02A011ASC022	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-ORANGE SILICONE SEPTA	20,420.00	70179090
CLS-CV02A111ASC011	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PTFE-SILICONE SEPTA	23,690.00	70179090
CLS-CV02A111ASC012	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-SILICONE SEPTA	25,210.00	70179090
CLS-CV02A111ASC021	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PTFE-ORANGE SILICONE SEPTA	21,280.00	70179090



Speak to our expert
+91-120-4968617

CLS-CV02A111ASC022	Combi pack of 11mm Amber Glass Crimp Vials 12 X 32 with Write On Patch and 11mm Aluminum Crimp Caps with PRE-SLIT PTFE-ORANGE SILICONE SEPTA	22,790.00	70179090
--------------------	--	-----------	----------

13 mm Screw Vials

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CB04C11300S000	13mm Clear Glass Screw Vials	17,730.00	70179090
CLS-CB04A11300S000	13mm Amber Glass Screw Vials	18,110.00	70179090

13 mm Screw Caps with Septa

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-BC000013PKS011	13mm Screw Caps with PTFE-SILICONE SEPTA	18,230.00	70179090

Combipacks of 13 mm Screw Vials with Matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV04C113PKS011	Combi pack of 13mm Clear Glass Screw Vials and 13mm Screw Caps with PTFE-SILICONE SEPTA	35,960.00	70179090

Combipacks of 13 mm Screw Amber Vials with Matching caps

Product Code	Description	Price Case of 5 x 100	HSN Code
CLS-CV04A113PKS011	Combi pack of 13mm Amber Glass Screw Vials and 13mm Screw Caps with PTFE-SILICONE SEPTA	36,340.00	70179090

EPA Pre-Assembled Vials

Product Code	Description	Price / 72 Pcs	HSN Code
CLS-CV20C124PWS111	20ml Clear Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	8,277.00	70179090
CLS-CV30C124PWS111	30ml Clear Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	9,167.00	70179090
CLS-CV40C124PWS111	40ml Clear Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	11,045.00	70179090
CLS-CV60C124PWS111	60ml Clear Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	15,087.00	70179090
CLS-CV20A124PWS111	20ml Amber Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	9,192.00	70179090
CLS-CV30A124PWS111	30ml Amber Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	10,278.00	70179090
CLS-CV40A124PWS111	40ml Amber Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	12,570.00	70179090
CLS-CV60A124PWS111	60ml Amber Glass EPA Vial with PTFE SILICONE SEPTA Bonded to a 24mm White Screw Cap, 72/pack	17,502.00	70179090



Speak to our expert
+91-120-4968617

20 ml Head Space GC Vials

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CB20F12000C000	20ml Clear Glass Flat Bottom Headspace Crimp Vial	13,760.00	70179090
CLS-CB20R12000C000	20ml Clear Glass Round Bottom Headspace Crimp Vial	13,760.00	70179090
CLS-CB20F12000S000	20ml Clear Glass Flat Bottom Headspace Screw Vial	18,610.00	70179090
CLS-CB20R12000S000	20ml Clear Glass Round Bottom Headspace Screw Vial	20,535.00	70179090

20 mm GC caps with septa

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-BC000020ASC051	Headspace Aluminum Crimp Caps with PTFE-SILICONE SEPTA	14,370.00	70179090
CLS-BC000020ASC041	Headspace Aluminum Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	16,120.00	70179090
CLS-BC000020MSC041	Headspace Magnetic Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	17,880.00	70179090

Headspace Magnetic Screw Cap Closure

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-BC000020MSS061	Headspace Magnetic Screw Caps with PTFE-SILICONE SEPTA	34,785.00	70179090

Combi packs of 20 mm GC Flat Bottom Vials with Matching Caps

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV20F120ASC051	Combi pack of 20ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with PTFE-SILICONE SEPTA	28,130.00	70179090
CLS-CV20F120ASC041	Combi pack of 20ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	29,880.00	70179090
CLS-CV20F120MSC041	Combi pack of 20ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Magnetic Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	31,635.00	70179090

Combi packs of 20 mm GC Round Bottom Vials with Matching Caps

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV20R120ASC051	Combi pack of 20ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with PTFE-SILICONE SEPTA	28,130.00	70179090
CLS-CV20R120ASC041	Combi pack of 20ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	29,880.00	70179090
CLS-CV20R120MSC041	Combi pack of 20ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Magnetic Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	31,635.00	70179090



Speak to our expert
+91-120-4968617

Combipacks of 20 mm GC Flat Bottom Vials with Matching Caps

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV20F120MSS061	Combi pack of 20ml Clear Glass Flat Bottom Headspace Screw Vial and Headspace Magnetic Screw Caps with PTFE-SILICONE SEPTA	50,155.00	70179090

Combipacks of 20 ml Round bottom screw Vials with Matching Caps

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV20R120MSS061	Combi pack of 20ml Clear Glass Round Bottom Headspace Screw Vial and Headspace Magnetic Screw Caps with PTFE-SILICONE SEPTA	55,320.00	70179090

10 ml Headspace Vials

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CB10F12000C000	10ml Clear Glass Flat Bottom Headspace Crimp Vial	7,005.00	70179090
CLS-CB10R12000C000	10ml Clear Glass Round Bottom Headspace Crimp Vial	1,669.00	70179090
CLS-CB10F12000S000	10ml Clear Glass Flat Bottom Headspace Screw Vial	2,330.00	70179090
CLS-CB10R12000S000	10ml Clear Glass Round Bottom Headspace Screw Vial	2,506.00	70179090

10 ml Flat Bottom Crimp Vials with closure

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV10F120ASC051	Combi pack of 10ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with PTFE-SILICONE SEPTA	21,375.00	70179090
CLS-CV10F120ASC041	Combi pack of 10ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	23,125.00	70179090
CLS-CV10F120MSC041	Combi pack of 10ml Clear Glass Flat Bottom Headspace Crimp Vial and Headspace Magnetic Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	26,715.00	70179090

10 ml Round Bottom Crimp Vials with closure

Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV10R120ASC051	Combi pack of 10ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with PTFE-SILICONE SEPTA	24,190.00	70179090
CLS-CV10R120ASC041	Combi pack of 10ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Aluminum Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	24,465.00	70179090
CLS-CV10R120MSC041	Combi pack of 10ml Clear Glass Round Bottom Headspace Crimp Vial and Headspace Magnetic Crimp Caps with Ultra Low Bleed PTFE-SILICONE SEPTA	28,055.00	70179090

10 ml HS screw Vials with closure


Product Code	Description	Price / Case lot of 2 x 100 Pcs	HSN Code
CLS-CV10F120MSS061	Combi pack of 10ml Clear Glass Flat Bottom Headspace Screw Vial and Headspace Magnetic Screw Caps with PTFE-SILICONE SEPTA	43,195.00	70179090

CLS-CV10R120MSS061	Combi pack of 10ml Clear Glass Round Bottom Headspace Screw Vial and Headspace Magnetic Screw Caps with PTFE-SILICONE SEPTA	47,315.00	70179090
--------------------	---	-----------	----------

Technical Information on Closures for Special purpose/ Storage vials

Product Code	Description	Pack Size	Compatible with
BC000015PKS011	15mm Closed top Screw Cap with PTFE-SILICONE septa 100/pk	100/pk X 5	CB08C01500S000
			CB08A01500S000
BC000018PCN000	18mm Clear Flip Top Cap 1000/pk	1000/pk x 5	CB10C01800N000
BC000022PCN000	22mm Clear Flip Top Cap 1000/pk	1000/pk x 5	CB15C02200N000
			CB25C02200N000
BC000018PKS011	18mm Black Closed top Screw Cap with PTFE-SILICONE Septa 100/pk	100/pk X 5	CB16C01800S000
			CB16A01800S000
BC000020PWS011	20mm White Closed Top Screw Cap with PTFE - SILICONE Septa 100/pk	100/pk x 5	CB20C02000C000
			CB20A02000C000

For more information



+91-120-4968617

For more information



www.compositelabline.com

COMPOSITE LAB LINE PVT LTD

UGSR-1, Ansal Plaza, Vaishali

Delhi, NCR, Ghaziabad - 201010 (India)

Ph. : 0120-4335636 Mail : sales@compositelablin.com

Visit.: www.compositelabline.com

Reg Office :


A-1, OKKP Complex, Near I.T.Crosing

Faizabad Road, Babujanj

Lucknow - 226020 (India)

Ph. : +919580613355, 0522 - 4068342

Mail : sales@compositelablin.com



Authorised Dealer: